



## 28R140N Dorso Arexa

<b>DE</b>	Gebrauchsanweisung .....	5
<b>EN</b>	Instructions for use .....	11
<b>FR</b>	Instructions d'utilisation .....	16
<b>IT</b>	Istruzioni per l'uso .....	22
<b>ES</b>	Instrucciones de uso .....	28
<b>PT</b>	Manual de utilização .....	34
<b>NL</b>	Gebruiksaanwijzing .....	40
<b>SV</b>	Bruksanvisning .....	46
<b>DA</b>	Brugsanvisning .....	51
<b>NO</b>	Bruksanvisning .....	57
<b>FI</b>	Käyttöohje .....	63
<b>PL</b>	Instrukcja użytkowania .....	68
<b>HU</b>	Használati utasítás .....	74
<b>CS</b>	Návod k použití .....	80
<b>BG</b>	Инструкция за употреба .....	85
<b>TR</b>	Kullanma talimatı .....	91
<b>RU</b>	Руководство по применению .....	97
<b>JA</b>	取扱説明書 .....	103
<b>ZH</b>	使用说明书 .....	109

**Material** Rostfreier Stahl / stainless steel, PE Schaumstoff / PE foam,  
Frotteebezüge / Terry Cloth Covers  
Kunststoff / plastic: PA, TPU, PP

## INFORMATION

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## INFORMATION

Date of last update: 2020-11-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 28R140N Dorso Arexa hyperextension orthosis.

## 2 Product description

### 2.1 Components

Scope of delivery (##02034)			
Item	Component	Item	Component
1	Sternal support, padded	8	Textile cover
2	Base plate, padded	9	Waistband, two-piece
3	Ratchet closure	10	Screws for waistband
4	Belt guiding element	11	Torx wrench
5	Extension element	12	Clip rivet for pads
6	Symphysis support, padded	13	Waistband pad
7	Back support, padded	–	Sternal support, short (to order if needed, article no.: 29R331=S)

### 2.2 Design

The ventral section consists of a base plate (2) with adaptable extension elements (5), two ratchet closures (3), a sternal support (1) and a symphysis support (6). The supports (1, 6) are continuously height-adjustable. The supports with the tubes are attached to the base plate by special clamping elements. The sternal support is also available in a short version (option).

The flexible, two-piece waistband (9) is attached to the back support (7). The waistband slides into the ratchet closures (3) and engages audibly. The

base plate and waistband as well as the sternal, symphysis and back supports are padded.

### 3 Intended use

#### 3.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

#### 3.2 Indications

- Stable traumatic vertebral body fractures of the lower thoracic spine and upper lumbar spine (also due to osteoporosis) without neurological damage (Th10 to L2).

Indications must be determined by the physician.

#### 3.3 Contraindications

##### 3.3.1 Absolute Contraindications

Not known.

##### 3.3.2 Relative Contraindications

The following conditions require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the support is fitted, tumours, lymphatic congestion – including unclear soft tissue swelling distal to the support, and dysaesthesia in the trunk.

#### 3.4 Effects

The orthosis straightens the spinal column in the sagittal plane and stabilises it. Based on the three-point principle with the back support (7), sternal, and symphysis supports (1, 6), the orthosis causes a lordosis, thereby relieving the anterior vertebral bodies.

### 4 Safety

#### 4.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 General safety instructions

### CAUTION

#### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

### CAUTION

#### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

### NOTICE

#### **Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)**

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

### NOTICE

#### **Improper moulding or application**

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be moulded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

## 5 Handling

### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.

- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

## 5.1 Size Selection

The orthosis is available in a universal size. It can be used for fitting of waist-band circumference up to 150 cm/59 inches.

## 5.2 Fitting and Application

### CAUTION

#### Loosening the ratchet closures

Loss of the support function, injuries due to wear or unintentional release of the ratchet closures

- ▶ At each application, ensure that the ratchet closures audibly engage in the waistband. The waistband must not slip out under load.
- ▶ Practice the application and removal of the product with the patient.
- ▶ In case of damage or noticeable changes in the product, the patient should consult a specialist workshop.

### CAUTION

#### Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

### INFORMATION

Continuous, direct skin contact with the product can lead to skin irritation. Only wear the product over clothing (e.g. a T-shirt or undershirt) in daily use. For limited times (e.g. water therapy), the product can also be worn directly on intact skin.

- > **Requirement:** The patient must be either lying down or standing up.
- 1) Loosen the screws in the base plate (see fig. 2).
  - 2) Taking into account the maximum extension length marked on the tubes, adjust the height of the sternal and symphysis supports to the patient (see fig. 3).
  - 3) Tighten the screws in the base plate.

- 4) If adaptation to the anatomical body contours of the patient is required, shape the tubes between the base plate and the supports manually over a suitable edge (see fig. 4). **Functional risks are associated with, for example, repeated cold forming and the resulting reduction in breaking elongation (strain hardening).**
- 5) If required, remove the same number of extension elements not needed to adjust the circumference from both sides and then reattach the belt guiding elements (see fig. 5). Shorten the base plate pad accordingly (see fig. 6).
- 6) Secure the two-piece waistband to the back support in the optimum position for the patient using countersunk screws (see fig. 7).
- 7) Position the back support dorsally and the base plate ventrally on the patient.
- 8) Insert the waistband symmetrically into the ratchet closures (see fig. 8). Fine adjustments can be made by means of the ratchet closures (see fig. 9).
- 9) Check the proper fit of the orthosis and correct if required.
- 10) Adjust the length of the waistband pad (see fig. 10) and fasten it with clip rivets (see fig. 11).
- 11) Shorten the ends of the waistband protruding from the ratchet closures, leaving some in reserve, and round them off to remove any burrs (see fig. 12).
- 12) To open the orthosis, press one of the buttons on the ratchet closure (see fig. 13).

### 5.2.1 Checks

Before handing over the product to the patient, check the following after applying the orthosis to the patient:

- Pressure, placement, and fit of the three supports (1, 6, 7).
- Tension and proper positioning of the waistband (9) and base plate (2).
- Fit and hyperextension effect in the sitting position.

## 6 Cleaning

The orthosis is not damaged by water, which means that personal hygiene (e.g. showering) and therapy in water that contains chlorine or salt are, in principle, possible when wearing the orthosis. The orthosis must subsequently be cleaned as follows:

- 1) Remove the textile covers from the supports and wash at **30 °C** in the washing machine.
- 2) Clean the orthosis with a damp cloth and neutral soap.

- 3) Rinse the orthosis with water and dry it with a cloth.

## 7 Disposal

Dispose of the product in accordance with national regulations.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 8.2 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.3 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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