2R30, 4R56, 4R56=1, 4R56=2

DE Gebrauchsanweisung	4
EN Instructions for Use	11
FR Instructions d'utilisation	18
IT Istruzioni per l'uso	25
ES Instrucciones de uso	33
PT Manual de utilização	40
NL Gebruiksaanwijzing	47
SV Bruksanvisning	54
DA Brugsanvisning	61
NO Bruksanvisning	69
FI Käyttöohje	75
PL Instrukcja użytkowania	83
HU Használati utasítás	90
CS Návod k použití	97
RO Instrucțiuni de utilizare	104
HR Upute za uporabu	111
SL Navodila za uporabo	118
SK Návod na používanie	125
BG Инструкция за употреба	132
TR Kullanma talimatı	139
EL Οδηγίες χρήσης	146
RU Руководство по применению	154
JA 取扱説明書	162
ZH 使用说明书	169
KO 사용 설명서	175

4R56


4R56=1


4R56=2


2R30

1

4R56, 4R56=1, 4R56=2

4


2R30

5

5 Nm

15 Nm

6

1 Produktbeschreibung

INFORMATION

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INFORMATION

Date of the last update: 2014-02-07

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

1.1 Construction and Function

The 4R56* Tube Clamp Adapters are specially designed for an adjustable connection between the pyramid adapter of the prosthetic hip joint and the pyramid adapter of the prosthetic knee joint or the rotator. The 2R30 Tube is used as a connection between two 4R56* Tube Clamp Adapters.

1.2 Combination Possibilities

For combination possibilities, please see the catalogue 646K2* or consult the manufacturer.

2 Application

2.1 Intended Use

The product is intended solely for lower limb prosthetic fittings.

2.2 Area of Application

Area of application according to the MOBIS® mobility system:

Approved for a body weight of up to **100 kg**.

The 2R30 Tube may be used only as a connecting element between two 4R56* Tube Clamp Adapters.

2.3 Environmental Conditions

Allowable Environmental Conditions
Temperature range for use: -10 °C to +60 °C (14 °F to 140 °F)
Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the user's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

3.2 General Safety Instructions

 CAUTION
Using the product without following the instructions for use Deteriorating health condition and damage to the product because of failure to follow the safety information
<ul style="list-style-type: none">▶ Observe the safety information in these instructions for use.▶ Make the patient aware of all safety instructions marked with "Inform the patient".

 CAUTION
Excessive strain on the product Fall due to breakage of load-bearing components
<ul style="list-style-type: none">▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").▶ Inform the patient.

⚠ CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

⚠ CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ **Inform the patient.**

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthesis components relative to each other and by the development of noises.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

2R30 Tube					
Fig.	Item No.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647H9=1

4R56, 4R56=1, 4R56=2 Tube Clamp Adapter					
Fig.	Item No.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647H9
4	1	▲	4	Set Screw	506G3=M8x12-V
4	–	●	–	Single-component pack	4D4
	Consisting of:				
	2	–	1	Cap Screw	–
	3	–	1	Straight pin	–
	4	–	1	Straight pin	–
4	6	■	1	Plastic ring	4X15

5 Preparation for Use

CAUTION

Incorrect alignment or assembly

Injuries due to fall resulting from damaged prosthesis components

- ▶ Observe the alignment and assembly instructions.

5.1 Adapting and Mounting the Tube

CAUTION

Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ For shortening the tube, use only a tube cutter.

CAUTION

Incorrect mounting of the tube

Fall due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthesis component when mounting.

CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

> **Recommended tools and materials:**

710D4 Torque Wrench, 719R3 Tube Cutter, 718R1 Tube Deburrer, degreasing cleaner (e.g. 634A3 Acetone)

- 1) Use the tube cutter to shorten the tube (see Fig. 5) to fit the user's measurements (see Fig. 2).
- 2) Use the deburring knife to carefully deburr the inside and outside of the cut edge (see Fig. 3).
- 3) Slacken the cap screw by **2 turns** (see Fig. 4, Item 2).
- 4) For definitive mounting, clean the insertion area of the connecting component on the inside and the tube on the outside using a degreasing cleaner.
- 5) While slightly turning the tube, slide it **5 mm** into the connecting component.
- 6) Take the connecting component with the inserted tube and place it on the work desk in a vertical position.
- 7) Slide the tube into the connecting component until the tube rests against the bottom of the clamp fitting.

- 8) Using the torque wrench, tighten the cap screw (for definitive mounting: **5 Nm** - see Fig. 6).

5.2 Mounting in the Modular Prosthesis

CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the installation torques provided.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

Mounting

> **Recommended tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

1) **For definitive mounting:**

Use Loctite® to secure the set screws (see Fig. 4, Item 1).

2) Turn in the set screws.

3) **For trial fitting:**

Use the torque wrench to tighten the set screws (**10 Nm**).

For definitive mounting:

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten (**15 Nm** - see Fig. 6).

4) Use suitable set screws to replace any set screws that extend out too far or screw in too deeply (see selection table).

Selection Table for Set Screws

Reference number	Length (mm)
506G3=M8x12-V	12
506G3=M8x14	14
506G3=M8x16	16

Adjustment

The 4 set screws can be used to make static adjustments at any time during alignment, trial fittings, and even after the prosthesis is finished.

Replacement and Disassembly

When replacing or disassembling the prosthesis component, the previously adjusted position can be maintained if the two deepest screwed-in set screws, which are NOT located opposite to each other, are removed.

6 Maintenance

- The prosthesis components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

8 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93 / 42 / EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH