WalkOn

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You have received the WalkOn, a dynamic lower leg orthosis made of carbon-fibre reinforced plastic. The name WalkOn designates the product range and the various types are the „Trimable“, „Flex“, „Reaction“, „Reaction plus“ and „Reaction junior“. The orthosis type is selected by your orthotist and the doctor who prescribed it

The orthosis is a medical device that provides support when you walk and stand to compensate impaired physical functions. We have compiled a few general questions and answers on the proper use and handling of the WalkOn: Your orthotist will be happy to answer questions about your specific situation and fitting of your orthosis.

What are the differences between the various types of WalkOn? (1)

- ▶ The classic 28U11 WalkOn and 28U23 WalkOn Trimable lift the foot when you swing your leg through to walk. The orthoses are made of a carbon fibre prepreg material. This makes them very lightweight and gives them moderate energy return* when walking. There is a carbon fibre spring at the back that extends from the Achilles tendon to the calf. Both orthoses stabilise the ankle joint.
- ▶ The 28U22 WalkOn Flex also ensures that your foot is lifted when you swing your leg through to walk. Constructed from a glass/carbon-fibre composite material and designed with a spiral extending from the inner ankle to the outer calf, the orthosis is particularly suitable if you want more flexibility in the ankle and knee. It is also recommended if you need little energy return.*
- ▶ The 28U24 WalkOn Reaction not only lifts the foot, it also supports the knee joint. For this reason its design differs considerably from the other orthoses. It is made entirely of a carbon-fibre prepreg material and has high energy return.*
- ▶ The 28U25 WalkOn Reaction plus has a similar functionality to the 28U24 WalkOn Reaction, but also offers enhanced forefoot stiffness, which permits a higher activity level and enables greater support for knee extension.
- ▶ The WalkOn Reaction junior was developed especially for use on children. Its characteristics are similar to those of the WalkOn Reaction plus. The children's version can also be combined with other orthoses.

*Energy return: the energy stored at heel strike and toe-off that is returned at foot roll-over. The aim is an almost natural gait.

Straightforward answers to key questions

The same principle applies to all orthoses: the better the fit, the better the effect. Consequently: if possible, insoles that are attached to the orthosis should be used. Combinations with foot and ankle orthoses are also possible with the WalkOn Reaction junior. (2)

How do I put on the WalkOn?

To be effective, the WalkOn should be worn with a stable shoe. We recommend that you first place the orthosis in the shoe and then insert your foot. The hook and loop closures should fit firmly on your leg without constricting it. To prevent friction and for reasons of hygiene, we recommend that you wear cotton socks. (3)

Do I need a special shoe for the WalkOn?

No, but you should select a shoe with a reinforced heel. It is important that the orthosis fits firmly so that it moves as little as possible. In addition, due to its pre-defined shape, correct function can be achieved only with the correct heel height. The recommended heel height for your shoe is 10 mm (+/- 5 mm), or as determined by your orthotist. Heels in any other shoes that you wear should have the same height.

Are there any movements should I avoid while wearing the WalkOn?

Basically there are no restrictions, but you should avoid extreme movements that flex the orthosis too much. For example, when walking and climbing stairs you should use your whole forefoot, not just the tips of your toes. (4)

Is it possible to take part in sport with the WalkOn?

Yes, in fact it is recommended that you are active. Speak to your orthotist to get their advice on which sports are the most beneficial and how to make the best use of your orthosis. You must avoid any extreme or abrupt movements. However, if you wish to take part in sporting activities that require special movements, please consult your orthotist.

Can I ride a bicycle or drive a car or motorcycle with the WalkOn?

It is not possible to make any specific recommendation. You must be able to safely control any vehicle for which you require a driving licence and insurance, or which is used in traffic. The following principles apply: The ability to operate a vehicle when wearing a WalkOn is determined on a case-by-case basis. Criteria include the type of fitting, clinical assessment and individual abilities.

All users are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.

How should I care for and clean the WalkOn?

The orthosis should be wiped regularly using water or a mild cleaning agent. The pad should be washed at least twice a week. We recommend machine laundering at 40°C (104°F). Please do not use fabric softener. To prolong the life time of the calf band, close the hook and loop closure each time you remove your WalkOn.

Can I get extra pads for the WalkOn?

Replacement pads for all orthoses are of course available and can be ordered by your orthotist. The orthotist should adapt the pad to your needs and make any necessary adjustments.

When should I contact my orthotist?

You should consult your orthotist regularly, at least every 6 months, in order to have the orthosis checked for any damage. The pads and closures should also be replaced regularly.

Please check the condition of your foot, lower leg and orthosis on a daily basis. Immediately report any unusual changes of your skin condition or the orthosis to your orthotist. The same applies if you detect any sharp edges on the orthosis.

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Otto Bock HealthCare GmbH