
17B66=16, 17B66=20 17B66=A-16, 17B66=A-20

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	26
PT	Manual de utilização	32
NL	Gebruiksaanwijzing	38
SV	Bruksanvisning	43
DA	Brugsanvisning	49
NO	Bruksanvisning	54
FI	Käyttöohje	59
PL	Instrukcja użytkowania	65
HU	Használati utasítás	71
CS	Návod k použití	76
RO	Instrucțiuni de utilizare	82
HR	Upute za uporabu	88
SL	Navodila za uporabo	93
SK	Návod na používanie	99
BG	Инструкция за употреба	105
TR	Kullanma talimatı	111
EL	Οδηγίες χρήσης	116
RU	Руководство по применению	123
JA	取扱説明書	130
ZH	使用说明书	136
KO	사용 설명서	140

INFORMATION

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INFORMATION

Last update: 2015-09-22

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.

► Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17B66 ankle joints.

2 Product description

2.1 Available sizes

The bilateral system ankle joint is available in 2 versions supplied in pairs.

	17B66=*		17B66=A*	
Joint bar connection	Inside		Outside	
System width [mm]	16	20	16	20

2.2 Components/design

Item	Component	Quantity	Item	Component	Quantity
1	Set screw	4	5	Bearing ball	4
2	Truss head screw	2	6	Stop pin (small)	8
3	Bearing nut	2	7	Spring	4
4	Upper joint section	2	8	Oval head screw	4

The scope of delivery does NOT include:

Item	Component	Article number
9	System lamination bar	17Y128=16x80, =20x80
10	System foot stirrup	17B113
-	System foot stirrup	17B114
-	System foot stirrup	17B115
-	System shoe stirrup	17B116

3 Application

3.1 Indications for use

The system ankle joints are intended to be used **exclusively** in pairs for the fitting of lower limb orthoses for patients.

The manufacturer recommends that the system ankle joint be processed by means of the lamination resin technique, thermoplastic technique or using carbon prepreg technology.

3.2 Indications

Partial or total paralysis of the leg muscles.

Indications must be determined by the physician.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Incorrect alignment, assembly or adjustment

Injuries caused by falling and damage to orthosis components

- ▶ Observe the alignment, assembly and adjustment instructions.

CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at a temperature above 300 °C (570 °F).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparation for use

5.1 Assembly

Gluing in the system lamination joint bars

After adapting the connecting elements, the lamination joint bars must be glued into the system connection plugs for reasons of stability:

- 1) Degrease the joining surfaces with 634A1 thinner or solvent.
- 2) Mix the 636W28 special adhesive and hardener.
- 3) Apply the mixture inside the system connection plugs.
- 4) Insert the lamination joint bars (9) into the upper joint section (4).
- 5) Fasten with 501T7=7.5x9xM5 oval head screws (8).
- 6) Allow to cure for at least **3-4 hours**. **INFORMATION: Final bonding strength will be reached after 16 hours.**

Mounting the foot stirrup

- 1) Insert the foot component with foot stirrup (10) into the upper joint section (4).
- 2) Lubricate the bearing nut (3) with 633F7 special lubricant and insert it.
- 3) Apply Loctite® 241 to the thread of the truss head screw (2) and screw it into the bearing nut (3).

5.2 Adaptation

The system ankle joint is adapted to the needs of the patient through the individual combination of the bearing ball (5), stop pins (6), and spring (7) with the set screws (1).

The following combinations are possible:

Combination		Application
Installation position anterior	Installation position posterior	
Stop pin		e.g. ICP or spina bifida
Spring	Stop pin	e.g. in case of distinct hyperextension in the knee joint
Spring stop	Stop pin	e.g. with distinct hyperextension in the knee joint, but with additionally adjustable damped dorsal stop
Spring		e.g. in case of weakness of the lower leg muscles
Stop pin	Spring	Dorsiflexion assist, e.g. in case of peroneal paralysis without own muscular knee joint protection
Spring stop	Spring	e.g. in case of pareses of the lower leg muscles
Spring	Spring stop	e.g. in case of weakness of the lower leg muscles with tendency to hyperextension of the knee joint
Stop pin	Spring stop	Dorsiflexion assist, e.g. in case of pareses of the lower leg muscles without muscular knee joint protection
Spring stop		e.g. in case of pareses in the lower leg muscles with weakness of the knee-protecting muscles

- ▶ Install the bearing ball (5), stop pins (6) and spring (7) with the set screws (1) in the combination required. Use Loctite® 241 to secure the set screws (1).

6 Cleaning

After contact with water containing salt, chlorine or soap, or when they get dirty, the system joints must be cleaned promptly.

- 1) Rinse the system joints with clean fresh water.
- 2) Dry the system joints with a cloth or allow to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Performing maintenance

- ▶ Lubricate the system ankle joints. Only use 633F7 special lubricant.

Individual spare parts and service set

Spare part	Single component	Service set	
	Article number	Quantity	Article number
Set screw (1)	506G3=M6x6	4	17B66=S
Truss head screw (2)	501S32=M6x14x10	2	
Bearing nut (3)	17Y93=9x6.25xM6	2	
Bearing ball (5)	509Y1=5.0	4	
Stop pin (6)	506A8=2.5x18	4	
	506A8=5x22	4	
Spring (7)	513D18=4.7x31	4	

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH