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### INFORMATION

Date of last update: 2020-07-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

## 1.1 Construction and Function

The 1C60 Triton, 1C63 Triton Low Profile and 1C64 Triton Heavy Duty prosthetic feet are suitable for fast walking on a variety of surfaces.

The functional properties of the prosthetic feet are achieved through the carbon spring elements. The prosthetic feet permit perceptible plantar flexion at heel strike and a natural rollover movement. Stored energy is released by the spring elements.

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

### Limited combination options for Ottobock components

Combination with mechatronic prosthetic knee joints	
<b>For spring stiffness 4 from size 29</b> <b>For spring stiffness 5 from size 25</b>	Contact Ottobock Customer Service before making a combination with a mechatronic prosthetic knee joint.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

1C60		
<b>Size [cm]</b>	21 to 24	25 to 30
<b>Max. body weight [kg], mobility grade 3</b>	100 (220 lbs)	150 (330 lbs)
<b>Max. body weight [kg], mobility grade 4</b>	100 (220 lbs)	125 (275 lbs)

1C63, 1C64		
<b>Size [cm]</b>	21 to 24	25 to 30
<b>Max. body weight [kg], mobility grade 3 and 4</b>	100 (220 lbs)	150 (330 lbs)

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and the size of the prosthetic foot.

<b>Spring stiffness</b>	
<b>Body weight [kg]</b>	<b>Prosthetic foot size [cm]</b>
	<b>21</b> <b>22</b> <b>23</b> <b>24</b> <b>25</b> <b>26</b> <b>27</b> <b>28</b> <b>29</b> <b>30</b>
<b>Up to 55 (121 lbs)</b>	1 —
<b>56 to 75 (123 lbs to 165 lbs)</b>	2
<b>76 to 100 (167 lbs to 220 lbs)</b>	3
<b>101 to 125 (222 lbs to 275 lbs)</b>	– 4
<b>126 to 150 (277 lbs to 330 lbs)</b>	– 5

## 2.3 Environmental conditions

### 1C60

<b>Allowable environmental conditions</b>
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

<b>Unallowable environmental conditions</b>
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

### 1C63, 1C64

<b>Allowable environmental conditions</b>
<b>Temperature range:</b> -10 °C to +45 °C (14 °F to 113 °F)
<b>Chemicals/liquids:</b> fresh water, salt water, perspiration, urine, soapsuds, chlorine water
<b>Moisture:</b> Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions
<b>Solids:</b> Dust, occasional contact with sand
<b>Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage</b> (see Page 21).

<b>Unallowable environmental conditions</b>
<b>Solids:</b> Highly hygroscopic particles (e.g. talcum), continuous contact with sand
<b>Chemicals/liquids:</b> Acids, continuous use in liquid media

<b>Storage and transport</b>
Temperature range -20 °C to +60 °C (-4 °F to +140 °F), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

## 2.4 Service life

This product has been tested by the manufacturer with 2 million load cycles according to ISO 22675. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

#### **Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

#### **Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 14).

 **CAUTION**

#### **Use under unallowable environmental conditions**

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

**⚠ CAUTION****Exceeding the service life and reuse on another patient**

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

**⚠ CAUTION****Using the prosthetic foot without footshell**

Risk of injury due to lack of surface grip and loss of functionality of the prosthetic foot

- ▶ Always use the prosthetic foot with a suitable footshell.

**⚠ CAUTION****Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**Signs of changes in or loss of functionality during use**

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the carbon spring are indications of loss of functionality.

**4 Scope of delivery**

Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Protective sock (black)	SL=Spectra-Sock-7
1	Heel wedge kit (soft, hard)	2F60*

**Spare parts/accessories (not included in the scope of delivery)**

Designation	Reference number
Footshell	2C6*

Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Connection cap	2C19*, 2C20*
Protective sock (white)	SL=Spectra-Sock

## 5 Preparing the product for use

### CAUTION

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### CAUTION

#### **Use of talcum**

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

### NOTICE

#### **Grinding the prosthetic foot or footshell**

Premature wear resulting from damage to the product

- ▶ Do not grind the prosthetic foot or footshell.

## 5.1 Alignment

### INFORMATION

There is a plastic adapter cover on the proximal connection of the prosthetic foot. This protects the connecting section from scratches during the alignment and trial fitting.

- ▶ Remove the adapter cover before the patient leaves the workshop/fitting area.

### 5.1.1 Applying/removing the footshell

#### INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
- ▶ Always use the prosthetic foot with the footshell.

- ▶ Apply or remove the footshell as described in the footshell instructions for use.

## 5.1.2 Bench Alignment

### TT bench alignment

Bench alignment process		
<b>Required materials:</b> 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)		
Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:		
<b>Sagittal plane</b>		
①	Heel height: <b>Effective heel height</b> (shoe heel height – sole thickness of fore-foot) + <b>5 mm</b>	
②	Exterior foot rotation: <b>approx. 5°</b>	
③	a-p positioning, middle of the prosthetic foot to the alignment reference line:	
	<b>Foot size [cm]:</b>	<b>Anterior placement, middle of foot to alignment reference line[mm]:</b>
	21 to 25	25
	26 to 28	30
	29 to 30	35
④	Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed.	
⑤	Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: <b>Individual residual limb flexion + 5°</b>	
<b>Frontal plane</b>		
⑥	Alignment reference line of prosthetic foot: <b>between big toe and second toe</b> Alignment reference line of prosthetic socket: <b>along the lateral patella edge</b>	
⑦	Observe the abduction or adduction position.	

### TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

#### 5.1.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219\***, TT modular leg prostheses: **646F336\***) may be requested from Ottobock.

### **5.1.4 Dynamic Trial Fitting**

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

#### **5.1.4.1 Optimising the heel characteristics**

The heel characteristics are optimised by using heel wedges. In case the heel strike or heel contact is too soft during the mid-stance phase, the heel can be stiffened by means of a heel wedge.

Two heel wedges are available (transparent=soft, black=hard). Ottobock recommends starting with the transparent heel wedge.

- 1) Align the recesses in the heel wedge, in the proximal and posterior directions.
- 2) Slide the heel wedge into position between the attachment spring and heel spring until it engages (see fig. 2).
- 3) To remove the heel wedge, push it out to the side.

### **5.2 Optional: Installing the foam cover**

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam cover should be stretched as little as possible in order to increase its service life. There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

- > **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive
- 1) Measure the length of the foam cover on the prosthesis and add the length allowance.  
**TT prostheses:** Distal allowance for movement of the prosthetic foot.  
**TF prostheses:** Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.
  - 2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.
  - 3) Pull the foam cover over the prosthesis.
  - 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.
  - 5) Install the prosthetic foot on the prosthesis.
  - 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
  - 7) Disassemble the prosthetic foot and remove the connecting element.
  - 8) Clean the connecting element using a degreasing cleaner.
  - 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
  - 10) Let the glue dry (approx. **10 minutes**).
  - 11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

## 6 Cleaning

### 1C60

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

### 1C63, 1C64

- > **Recommended cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)
- 1) Clean the product with clear fresh water and a pH neutral soap.
  - 2) Rinse the soap away with clear fresh water.
  - 3) Dry the product with a soft cloth.
  - 4) Allow to air dry in order to remove residual moisture.

## 7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 10 Technical data

1C60 Triton											
Sizes [cm]		21	22	23	24	25	26	27	28	29	30
With narrow foot-shell	Heel height [mm]	15 ± 5							–		
	System height [mm]	116	117	118	120	124	125	136	–		
	Weight [g]	490	500	565	590	650	660	720	–		
With normal foot-shell	Heel height [mm]	–			10 ± 5						
	System height [mm]	–			126	129	131	140	141	154	156

1C60 Triton											
Sizes [cm]		21	22	23	24	25	26	27	28	29	30
With normal foot-shell	Weight [g]	–			600	660	680	730	750	810	840
	Max. body weight [kg], mobility grade 3	100 (220 lbs)				150 (330 lbs)					
Max. body weight [kg], mobility grade 4		100 (220 lbs)				125 (275 lbs)					
Mobility grade		3 and 4									

1C63 Triton Low Profile											
Sizes [cm]		21	22	23	24	25	26	27	28	29	30
With narrow foot-shell	Heel height [mm]	15 ± 5							–		
	System height [mm]	35			36	40		45	–		
	Weight [g]	435	450	520	550	600	615	665	–		
With normal foot-shell	Heel height [mm]	–			10 ± 5						
	System height [mm]	–			42	43	45	49		52	55
	Weight [g]	–			560	610	635	675	700	820	855
Max. body weight [kg], mobility grade 3 and 4		100 (220 lbs)				150 (330 lbs)					
Mobility grade		3 and 4									

1C64 Triton Heavy Duty											
Sizes [cm]		21	22	23	24	25	26	27	28	29	30
With narrow foot-shell	Heel height [mm]	15 ± 5							–		
	System height [mm]	116	117	118	120	124	125	136	–		
	Weight [g]	565	575	640	665	725	735	795	–		
With normal foot-shell	Heel height [mm]	–			10 ± 5						
	System height [mm]	–			126	129	131	140	141	154	156
	Weight [g]	–			675	735	755	805	825	885	915
Max. body weight [kg], mobility grade 3 and 4		100 (220 lbs)				150 (330 lbs)					
Mobility grade		3 and 4									

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---

---

**The Triton family 1C60, 1C63 and 1C64 is covered by the following patents:**

European Patent: EP 2 509 540 in DE, FR, GB, IT, SE, TR  
Japan: Patent No. 5 781 531  
Russia: RU 2 550 003  
China: ZL201080053507.3  
Patents pending in Brazil, China, USA, and as European Patent Application.

**The 1C60 Triton and the 1C64 Triton Heavy Duty are covered by following registered design and design patent:**

European Design: No. 001810425  
USA: Patent US D653,759 S  
IR-Design: DM/076079 in TR

Ottobock SE & Co. KGaA  
Max-Näder-Straße 15 · 37115 Duderstadt · Germany  
T +49 5527 848-0 · F +49 5527 848-3360  
healthcare@ottobock.de · www.ottobock.com

Otto Bock HealthCare LP  
3820 West Great Lakes Drive  
Salt Lake City, UT 84120 · USA  
T +1 800 328 4058 · F +1 800 655 4963