## 8E33=\*, 8E34=\*

<b>DE</b> Gebrauchsanweisung (Benutzer) .....	3
<b>EN</b> Instructions for use (user) .....	17
<b>FR</b> Instructions d'utilisation (Utilisateur) .....	30
<b>IT</b> Istruzioni per l'uso (Utilizzatore) .....	45
<b>ES</b> Instrucciones de uso (Usuario) .....	58
<b>PT</b> Manual de utilização (Usuário) .....	73
<b>NL</b> Gebruiksaanwijzing (Gebruiker) .....	87
<b>SV</b> Bruksanvisning (Användare) .....	101
<b>JA</b> 取扱説明書(ユーザー用) .....	114

## INFORMATION

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### INFORMATION

Last update: 2016-11-07

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Receive instruction in the correct and safe use of the product by qualified personnel.
- ▶ Please contact the qualified personnel if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The 8E33=7-1, 8E34=7-1, 8E33=9-1, 8E34=9-1 System Electric Greifers are referred to as the product/prosthesis/terminal device in the following.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Design

The product consists of the following components:

1. **Fingertips**  
The position can be rotated by  $\pm 60^\circ$  to adapt to various tasks. For other types of fingertips, see the section "Accessories" (see Page 23)
2. **Safety lever**  
To unlock the prehensile arm drive (see Page 27)
3. **Dial**  
For manual opening/closing
4. **On/Off switch**  
To turn the electronics and flexion joint on/off  
Position of the rocker switch towards the terminal device: electronics switched on and flexion brake released  
Position of the rocker switch towards the prosthetic socket: electronics switched off and flexion joint blocked
5. **Quick-disconnect wrist unit** (only with 8E33\*)  
Mechanical connection to the lamination ring

### 2.2 Function

The product is a myoelectrically controlled terminal device of the MyoBock prosthesis system.

Various control programmes make optimal adaptation to the individual needs and abilities possible.

### Only with 8E33=9\* and 8E34=9\*

The terminal device is based on the Ottobock DMC system (DMC=Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally control the gripping speed and gripping force in accordance with the muscle signal.

## 3 Application

### 3.1 Indications for use

The System Electric Greifers are intended **exclusively** for exoprosthetic upper limb fittings.

### 3.2 Use/Field of Application

The System Electric Greifers can be used for unilateral or bilateral amputations, and in case of dysmelia.

### 3.3 Environmental conditions

The prosthesis system was developed for everyday activities and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the prosthesis system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see Page 29).

### 3.4 Qualification

Fitting a patient with the product may only be carried out by orthopaedic technicians who have been authorised by Ottobock after completion of a corresponding training course.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

 <b>CAUTION</b>
<p><b>The heading describes the source and/or the type of hazard</b></p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none"><li>&gt; E.g.: Consequence 1 in case of failure to observe the hazard</li><li>&gt; E.g.: Consequence 2 in case of failure to observe the hazard</li><li>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</li></ul>

### 4.3 General safety instructions

#### **WARNING**

##### **Operating a vehicle with the prosthesis**

Accidents due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).

#### **WARNING**

##### **Operating machines with the prosthesis**

Injury due to unexpected actions of the prosthesis.

- ▶ The prosthesis should not be used for the operation of industrial machines or motor-driven equipment.

#### **WARNING**

##### **Operating the prosthetic system near active implanted systems**

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.

- ▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

#### **WARNING**

##### **Use of damaged power supply unit, adapter plug or battery charger**

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

#### **WARNING**

##### **Charging the prosthesis without taking it off**

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

#### **CAUTION**

##### **Independent manipulation of the product**

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### **CAUTION**

##### **Signs of wear to system components**

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of user safety and in order to maintain the operational reliability and the validity of the warranty, the specified service intervals must be observed.

#### 4.4 Information on Proximity to Certain Areas

##### **⚠ CAUTION**

##### **Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

##### **⚠ CAUTION**

##### **Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850 / GSM 900: 0.99 m
  - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points ...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

##### **⚠ CAUTION**

##### **Wearing in extreme ambient temperatures**

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see Page 29).

#### 4.5 Information on Use

##### **⚠ CAUTION**

##### **Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)**

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

##### **⚠ CAUTION**

##### **Risk of pinching between the prehensile arms/fingertips**

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the prehensile arms/fingertips when using the product.

### **⚠ CAUTION**

#### **Overloading due to unusual activities**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional. If necessary, he or she will forward the product to an authorised Ottobock Service Centre.

### **⚠ CAUTION**

#### **Mechanical stress on the product**

Injury due to faulty control or malfunction of the prosthetic system.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

### **⚠ CAUTION**

#### **Penetration of dirt and humidity into the prosthetic components**

Injury due to faulty control or malfunction of the prosthesis system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the prosthesis system or prosthetic components (e.g. the terminal device).

### **NOTICE**

#### **Improper care of the prosthesis system**

Damage to the prosthesis system due to the use of incorrect cleaning agents.

- ▶ Clean the prosthesis system with a damp cloth and mild soap only (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Use only the following products for cleaning/disinfecting the inner socket:  
**Cleaning:** 453H10=1 Ottobock DermaClean  
**Disinfection:** 453H16 Ottobock DermaDesinfect

### **⚠ CAUTION**

#### **Changing prosthetic components when switched on**

Injury due to faulty control or malfunction of the prosthesis system.

- ▶ Before changing prosthetic components (such as the terminal device), remove the rechargeable battery from the battery receptacle or switch off the prosthesis system.

## **4.6 Information on the Power Supply/Battery Charging**

### **⚠ CAUTION**

#### **Charging the prosthesis system with soiled or damaged contacts**

Injury due to unexpected behaviour of the prosthesis system caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.

- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

#### NOTICE

##### **Charging the prosthesis system with an incorrect power supply/battery charger**

Damage to the prosthesis system due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for the prosthesis system by Ottobock (see instructions for use and catalogues).

#### NOTICE

##### **Contact of the charging plug with magnetic data carriers**

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

#### INFORMATION

The prosthesis system cannot be used as long as the charging plug is connected to the charging receptacle. The prosthesis system is deactivated for the duration of charging.

### 4.7 Notice for using the torch

#### ⚠ WARNING

##### **Use of the incorrect battery type and/or failure to note the polarity**

Explosion due to escaping gases and their self-ignition caused by short circuit and excessive temperature.

- ▶ Only use type CR2032 batteries.
- ▶ When inserting the two batteries, note the polarity information on the cover of the battery compartment.

#### ⚠ CAUTION

##### **Looking directly into the beam of light from the torch**

Eye damage due to looking directly into the beam of light.

- ▶ Do not use the torch for examinations, for instance of the eyes (pupil test).
- ▶ Do not aim the beam of light directly at the eyes of people or animals.
- ▶ Do not look directly into the beam of light.
- ▶ If the beam of light is aimed at the eyes, be sure to close your eyes and immediately move the head out of the beam of light.

### 4.8 Notes on using an elbow joint (e.g. 12K44, 12K50, 12K100N, 12K110N)

#### ⚠ CAUTION

##### **Risk of pinching in elbow joint flexion region**

Injury due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the elbow joint.

#### ⚠ CAUTION

##### **Manual unlocking of elbow lock under load**

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.

- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

### **⚠ CAUTION**

#### **Penetration of dirt and moisture in the prosthetic arm**

Injury due to unexpected behaviour of the prosthesis system or malfunction.

- ▶ Do not let foreign particles or liquids get into the prosthetic arm.
- ▶ Do not expose the prosthetic arm, especially the elbow joint, to dripping or splashing water.
- ▶ Wear the prosthetic arm and especially the elbow joint under suitably resistant clothing in the rain.

### **NOTICE**

#### **Coating, gluing or painting the prosthesis**

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

## **5 Scope of Delivery and Accessories**

### **5.1 Scope of delivery**

- 2 pc. CR2032 lithium battery
- 1 pc. 8E33=9-1 DMC VariPlus System Electric Greifer or  
1 pc. 8E34=9-1 DMC VariPlus System Electric Greifer or  
1 pc. 8E33=7-1 Digital Twin System Electric Greifer or  
1 pc. 8E34=7-1 Digital Twin System Electric Greifer
- 1 pc. 647G1275 instructions for use (user)

### **5.2 Accessories**

#### **Fingertips**

- 9S138 tip set
- 9S145 rubber gripping pad set
- 9S149 fingertip blank set
- 9S234=PAA tip pair
- 9S278=PAA tip pair

## **6 Charging the battery**

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product every day when used on a daily basis.
- When used daily, the complete charging unit (power supply – battery charger) may remain plugged into the wall socket.
- The battery should be charged for at least 3 hours prior to initial use.

### **6.1 757L14 battery charger and 757B15 X-ChangePack**

### **NOTICE**

#### **Deep discharge of the battery**

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The red LED lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the red LED flashes ☀.
- 5) Remove the rechargeable battery after the charging process is complete.

### 6.1.1 Display of the current charge level

The current charge level can only be checked on the battery charger with connected power supply unit and inserted rechargeable battery.

LED on pulse charger	Event
●	Charge level 0% to 79%
☀	<b>Long flash:</b> charge level 80% to 99% <b>Short flash:</b> charge level 100%

### 6.1.2 Troubleshooting

LED on pulse charger	Date	Error	Solution
☀	Immediately after the battery is inserted in the charging slot	Defective battery	The battery must be inspected by an authorised Ottobock Service Centre.
●	After 4 or 8 hours of charging	Defective battery	The battery must be inspected by an authorised Ottobock Service Centre.
○	Immediately after the battery is inserted in the charging slot	Defective charging slot/pulse charger or power supply	The pulse charger and power supply must be inspected by an authorised Ottobock Service Centre.

## 6.2 757L20 battery charger and 757B2\* EnergyPack


### NOTICE

#### Deep discharge of the battery

Irreparable damage to the product and resulting loss of functionality in case of irregular and incomplete charging of the battery.

- ▶ In case of long-term storage, fully charge the battery every 4 to 6 months.

- 1) Plug the charging station with power supply unit into the wall socket.  
→ The green LED on the power supply unit and the green LED in the middle of the LED row light up.
- 2) Insert one or two rechargeable batteries into the charging station.
- 3) The green LED for the respective charging socket lights up ● and the charging process starts.
- 4) When the rechargeable battery is fully charged, the green LED for the respective charging socket flashes ☀.
- 5) Remove the rechargeable battery after the charging process is complete.

### 6.2.1 Display of the current charge level

When inserting/engaging the rechargeable battery in the prosthesis, a charge level indicator is activated on the battery for a few seconds.

LED indicator	Event
●	Charge level above 50%
● and ●	Charge level below 50%
●	Charge level below 5%

### 6.2.2 Troubleshooting

Charging socket	LED for the charging socket	Event
Empty	● or ● or ●	<b>The battery charger is defective</b> The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
Empty or rechargeable battery inserted	☀	
Rechargeable battery inserted	●	<b>The rechargeable battery is charging</b>
	☀	<b>The rechargeable battery is fully charged</b>
	●	<b>The temperature of the rechargeable battery is too high</b> Remove the rechargeable battery and let it cool down
	●	<b>The inserted rechargeable battery is defective</b> The rechargeable battery must be inspected by an authorised Ottobock Service Centre.

## 6.3 757L35 battery charger and 757B35=\* rechargeable battery

### 6.3.1 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle of the prosthesis system.  
→ The correct connection of the battery charger to the prosthesis system is indicated by feedback (Status signals).
- 2) The charging process starts.  
→ The prosthesis system is switched off automatically.
- 3) Disconnect the prosthesis system after the charging process is complete.

### 6.3.2 Safety Shutoff

The purpose of the prosthesis system safety shutoff is to protect the rechargeable battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit

After a short circuit the charging plug must be plugged into the charging receptacle and then unplugged in order for the electronics to be activated.

### 6.3.3 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthesis system switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level.


Charging receptacle	Status
●	Battery fully charged
●	Battery 50% charged
●	Battery empty

### 6.3.4 Beep signals

Beep signal	Additional display	Status
1 x long	-	<ul style="list-style-type: none"> <li>• Turn off prosthesis on charging receptacle</li> <li>• Charging begins (charging plug connected to charging receptacle)</li> <li>• Charging ends (charging plug disconnected from charging receptacle)</li> </ul>
2 x short	LED on charging receptacle lights briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis system shuts off automatically

## 7 Use

### 7.1 Switching the terminal device on and off

#### Switching on


Move the rocker switch to position 1 until it engages. The electronics of the terminal device are switched on and the flexion brake is released.

#### Switching off

Move the rocker switch to position 2 until it engages. The electronics of the terminal device are switched off and the flexion joint is blocked.

### 7.2 Torch

(for terminal devices with built-in torch)

#### 1. **2x powerful LEDs**

to light the working environment

#### 2. **On/Off switch**


Pressing this button turns the torch on or off. If this button is pressed immediately after switching on, the brightness can be changed in 3 levels by pressing it again.

After about 15 minutes, the torch is switched off automatically.

#### 3. **Battery compartment**

### Inserting the batteries (2x CR2032)

To insert the batteries, complete the following steps:

- 1) Turn the cover on the battery compartment counter-clockwise and remove it.
- 2) Insert 2 CR2032 button cells with the + pole facing the cover.
- 3) Install the cover. In doing so, make sure the rubber seal is not damaged.
- 4) Close the cover by turning it clockwise.

### 7.3 Opening the terminal device in an emergency

This safety function allows the prehensile arms to be opened regardless of the control signals.

- Swing the unlock lever away from the dial (arrow) and open the prehensile arms manually.

## 7.4 Changing the 8E33=7\*, 8E33=9\* terminal device

### Separating the terminal device from the socket

- 1) Rotate the terminal device once around its own axis until a slight resistance is felt (about 360°).
- 2) Overcome this resistance and pull the terminal device off the socket.

### Attaching the terminal device to the socket

- 1) Insert the quick-disconnect wrist unit into the lamination ring and press it in firmly.
- 2) Turn the terminal device slightly to the left or right.
- 3) Verify proper attachment by pulling on the terminal device.

## 8 Cleaning

### NOTICE

#### Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.  
Ensure that no liquid penetrates into the system component(s).
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

## 9 Maintenance

Since all moving mechanical parts are subject to wear, regular service inspections are recommended. These service inspections include an inspection of the product and, when necessary, the replacement of wear parts.

## 10 Legal information

### 10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### 10.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

The product meets the requirements under the RoHS Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

## 11 Technical data

Environmental conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Storage and transport without packaging	-20 °C/-4 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Operation	-10 °C/+14 °F to +60 °C/+140 °F max. 93 % relative humidity, non-condensing

Reference number	8E33=7-1	8E34=7-1	8E33=9-1	8E34=9-1
Operating voltage [V]	6/7.2			
Opening width [mm]	95			
Gripping force [N]	160		0 - 160	
Weight [g]	540	520	540	520
Service life [years]	5			
Fingertip torque value [Ncm]	40 - 50			
Proportional speed [mm/sec]	180		8 - 200	

Torch	
Illuminants [number, type]	2, LED
Total maximum luminosity of the two LEDs [mcd]	16400-24000
Battery [number, type]	2x CR2032
Maximum operating time with new batteries at highest brightness and at room temperature [h]	8
Brightness levels [number]	3
Automatic switch-off	after 15 minutes

Rechargeable battery reference number	757B15	757B20	757B21	757B35 =0	757B35 =1	757B35 =3
Temperature range for charging [°C]	0°C to +40°C	0°C to +50°C		0°C to +60°C		
Capacity [mAh]	550	900	800	300	600	1150
Rechargeable battery service life [years]	2					
Product behaviour during the charging process	The product is non-functional					
Operating time of the product with fully charged battery [grip cycles]	approx. 1000	approx. 2500 – 3000	approx. 2000 – 2500	approx. 1000	approx. 1000 – 2000	approx. 4000
Charging times (rechargeable battery fully discharged) [hours]	approx. 10	approx. 3.5	approx. 3	approx. 2.5		
Nominal voltage [V]	6	approx. 7.2		approx. 7.4		
Battery technology	Ni-Mh	Li-Ion		Li-Po		Li-Ion
Approved battery chargers	757L14	757L20		757L35		

<b>Power supply</b>	
Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V ==

## 12 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN YYYY WW NNN

Serial number  
YYYY – year of manufacture  
WW – week of manufacture

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Otto Bock Healthcare Products GmbH