





## 50P20 Manu Arexa 50P21 Manu Arexa Pollex

<b>DE</b> Gebrauchsanweisung .....	5	<b>FI</b> Käyttöohje.....	35
<b>EN</b> Instructions for use .....	8	<b>PL</b> Instrukcja użytkowania.....	38
<b>FR</b> Instructions d'utilisation.....	11	<b>HU</b> Használati utasítás .....	41
<b>IT</b> Istruzioni per l'uso .....	14	<b>CS</b> Návod k použití.....	45
<b>ES</b> Instrucciones de uso .....	17	<b>SK</b> Návod na používanie .....	48
<b>PT</b> Manual de utilização.....	20	<b>ET</b> Kasutusjuhis .....	51
<b>NL</b> Gebruiksaanwijzing.....	23	<b>RU</b> Руководство по применению .....	54
<b>SV</b> Bruksanvisning .....	27	<b>JA</b> 取扱説明書.....	57
<b>DA</b> Brugsanvisning .....	29	<b>ZH</b> 使用说明书.....	60
<b>NO</b> Bruksanvisning .....	32		

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### INFORMATION

Last update: 2019-11-28

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information on the fitting and application of the 50P20 Manu Arexa and 50P21 Manu Arexa Pollex wrist braces.

## 2 Intended use

### 2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the upper limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 2.2 Indications

- Sprains of the wrist joint or the thumb saddle joint
- Rheumatoid arthritis
- Wrist drop due to paresis or paralysis of the radial nerve

### 50P20 Manu Arexa additional indications

- Tendinopathy of the wrist joint

- Ligament injury, ligament strain in the wrist
- Wrist instability
- Impingement syndrome
- Irritation in the wrist area
- Carpal tunnel syndrome

### 50P21 Manu Arexa Pollex additional indications

- Fractures of the distal forearm or carpal bone
- Ligament rupture (partial rupture)
- Partial arthrodeseis
- Preoperative/postoperative/posttraumatic conditions that require immobilization of the wrist
- Irritation in the wrist and thumb areas
- Night positioning

The indication must be determined by the physician.

## 2.3 Contraindications

### 2.3.1 Absolute Contraindications

None known.

### 2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

### 2.4 Mechanism of Action

The brace stabilizes the wrist.

### 3 Safety

#### 3.1 Explanation of warning symbols

**⚠ CAUTION**

Warning regarding possible risks of accident or injury.

**NOTICE**

Warning regarding possible technical damage.

#### 3.2 General safety instructions

**⚠ CAUTION****Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

**⚠ CAUTION****Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

**NOTICE****Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

**NOTICE****Improper molding or application**

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

**NOTICE****Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

### 4 Handling

**INFORMATION**

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

#### 4.1 Selecting the Size

- 1) Measure the wrist circumference.
- 2) Determine the brace size (see size chart).

## 4.2 Fitting and Application

### ⚠ CAUTION

#### Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

The 50P20 and 50P21 braces are pre-configured for the right hand upon delivery (see fig. 1, see fig. 2). The metal splints must be modified for use on the left hand (see fig. 3). Two metal splints – a palmar spoon-shaped splint and a dorsal splint – are located in the outer pockets of the brace. Shape the splints by hand to customize them for the user.

- 1) Open the brace.
- 2) Remove the palmar splint (A) from the outer pocket of the brace (see fig. 6, see fig. 4).
- 3) Mold the splint so it fits the anatomical shape of the patient's palm.
- 4) Insert the splint into the pocket with the outer spoon-shaped curve facing the palm.
- 5) **50P20: The dorsal splint (B1) can optionally be used on the ulnar side (B2).** Remove the splint from the outer pocket, shape it by hand, and place it in the dorsal or ulnar pocket (see fig. 4).
- 6) **50P21:** Remove the dorsal splint (B1) from the outer pocket, shape it by hand, and place it in the pocket.
- 7) Pull the brace over the wrist, guiding the thumb through the thumb opening (see fig. 5). Tighten the lacing and fasten the hook-and-loop closure.
- 8) Check the fit of the brace.

- 9) **Optional:** Trim the hook-and-loop closure with scissors to better customize the fit of the brace (see fig. 7).  
→ The brace fits snugly without pinching.

## 4.3 Removal

- 1) Open the hook-and-loop closure of the brace.
- 2) Take off the brace.

## 5 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- 1) Remove all metal splints.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the brace in warm water at **40 °C (104 °F)** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Reinsert the metal splints.

## 6 Disposal

Dispose of the product in accordance with national regulations.

## 7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disre-

garding the information in this document, particularly due to improper use or unauthorized modification of the product.

## 7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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