
# ottobock.

CE

## 13E200

DE	Gebrauchsanweisung (Fachpersonal) .....	3
EN	Instructions for use (qualified personnel) .....	14
FR	Instructions d'utilisation (Personnel spécialisé) .....	25
IT	Istruzioni per l'uso (Personale tecnico specializzato) .....	36
ES	Instrucciones de uso (Personal técnico especializado) ..	47
PT	Manual de utilização (Pessoal técnico) .....	59
NL	Gebruiksaanwijzing (Vakmensen) .....	70
SV	Bruksanvisning (Fackpersonal) .....	82
DA	Brugsanvisning (Faguddannet personale) .....	92
HU	Használati utasítás (szakszemélyzet) .....	103
CS	Návod k použití (Odborný personál) .....	114
RO	Instrucțiuni de utilizare (Personal de specialitate) .....	125
TR	Kullanma talimatı (Uzman personel) .....	136
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό) .....	147
RU	Руководство по применению (Квалифицированный персонал) ..	158
JA	取扱説明書 (有資格担当者) .....	170

## INFORMATION

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### INFORMATION

Date of last update: 2021-04-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "13E200 electrode" is referred to simply as the product/electrode below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Function

The product samples muscle action potential produced by muscle activity on the surface of the skin. Unwanted interference signals are filtered out by the product.

## 3 Intended use

### 3.1 Indications for use

The product is to be used **exclusively** for the control of myoelectric upper limb prosthetic components.

### 3.2 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a patient with the product may only be carried out by an O&P professional.

## 4 Safety

### 4.1 Explanation of warning symbols

 <b>WARNING</b>	Warning regarding possible serious risks of accident or injury.
 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
<b>NOTICE</b>	Warning regarding possible technical damage.

## 4.2 Structure of the safety instructions

### WARNING

#### **The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

## 4.3 General safety instructions

### WARNING

#### **Non-observance of safety notices**

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

### WARNING

#### **Fitting of children**

Injury due to swallowing small parts.

- ▶ Never leave children unattended during the fitting.

### CAUTION

#### **Insufficient skin contact of the electrodes**

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.

- ▶ In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary.
- ▶ If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

**⚠ CAUTION**

**Incorrect electrode settings due to muscle fatigue**

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

**⚠ CAUTION**

**Damage to the electrode cable due to kinking or abrasion**

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ Replace the damaged electrode cable immediately.

**⚠ CAUTION**

**Turning the adjustment regulator too far**

Injury due to faulty control or malfunction of the product.

- ▶ Do not turn the adjustment regulator further than the noticeable stop.

**⚠ CAUTION**

**Slipping electrode**

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ Inform the patient that uncontrolled movement of the prosthesis component may occur if the electrode slips.

**⚠ CAUTION**

**Penetration of liquids**

Injury due to faulty control or malfunction of the product.

- ▶ Make sure that no liquids can penetrate into the product.
- ▶ Do not use the product for bathing prostheses.

**⚠ CAUTION**

**Establishing a connection between the skin and metal components of the prosthesis with the use of carbon fibre**

Injury due to unexpected product behaviour as a result of interference conducted by carbon fibres.

- ▶ During fabrication, make sure that no connection is formed between the skin and metal prosthesis components by carbon fibres.

**NOTICE**

**Improper care of the housing**

Damage to the casing through the use of acetone, white spirit or similar solvents.

- ▶ Only clean the housing with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

**NOTICE**

**Thermal overloading due to external heat supply**

Destruction of the product.

- ▶ Do not expose the electrode to hot air or other sources of heat.

## 4.4 Patient information

**⚠ CAUTION**

**General use of electrodes**

Skin irritation may occur under the electrodes.

- ▶ After removing the prosthetic components, check the skin beneath them for pressure sores and irritation.
- ▶ Observe safety information in the instructions for use.

**⚠ CAUTION**

**Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850/GSM 900: 0.99 m
  - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points,...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

**NOTICE**

**Improper product care**

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Do not bring the surface of the electrodes into contact with a special cleaner for Ottobock cosmetic gloves, e.g 640F12.

**⚠ CAUTION**

**Damaged electrode cable**

Injury due to unexpected product behaviour as a result of uncontrolled activation of the prosthesis component.

- ▶ Power down the system and contact your orthopaedic technician.

**⚠ CAUTION**

**Use of a damaged product**

Injury due to loss of product functionality.

- ▶ Prior to use, conduct an external visual inspection to verify that all parts of the product are undamaged.
- ▶ In case of damage, have the product repaired promptly.

**⚠ CAUTION**

**Penetration of product with dirt and humidity**

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

**⚠ CAUTION**

**Remaining in areas outside the allowable temperature range**

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 23).

**⚠ CAUTION**

**Mechanical stress on the product**

Injury due to faulty control or malfunction of the product.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**⚠ CAUTION**

**Independent manipulation of the product**

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.

- The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

## **5 Scope of Delivery and Accessories**

### **5.1 Scope of delivery**

- 1 pc. 13E200=50/13E200=60 electrode
- 1 pc. 13E80 adjustment tool
- 1 pc. Instructions for use (qualified personnel)

### **5.2 Accessories**

- 1 pc. 13E191 casting template for inner socket
- 1 pc. 13E192 lamination template
- 2 pc. 507S15 lamination plate with bore
- 1 pc. 13E201 electrode accessories for vacuum-formed inner sockets
- 2 pc. 503F3 socket screw with Allen head

#### **The following products are also required:**

- 757M11=X-Change MyoBoy
- 13E129=\* electrode cable

## **6 Preparing the product for use**

### **6.1 Connecting the product**

- 1) Shorten the electrode cable as needed.
- 2) Push the electrode cable with the grey side up into the plug connection all the way to the stop (see fig. 2).
- 3) Grease the connector receptacle of the product and the pin contacts with 633F11 silicone grease.
- 4) Apply 633F11 silicone grease to the electrode cable receptacle.
- 5) Insert the plug connection with connected electrode cable as far as possible into the product (see fig. 3).
- 6) Remove any silicone grease that seeps out.

## Electrode cable replacement

- 1) Pull off the electrode cable at right angles to the product (see fig. 4).
- 2) Cut off at least 5 mm of the electrode cable and push it back into the plug connection.

### 6.2 Inserting the product

- ▶ Insert the product into the inner socket from the outside.  
For slipping the electrode into the mounting channel, slightly press the suspension arms together. The suspension arms engage if slight pressure is applied.

### 6.3 Adjusting the product

- 1) Use the adjustment tool to set the product sensitivity (see fig. 5).
- 2) After a pause of 10 to 15 minutes, readjust the sensitivity.

#### INFORMATION

The patient must be able to hold the maximal value of the control unit for at least 2 seconds. If this is not observed, the product setting is overly sensitive which makes it more difficult for the patient to control the prosthesis.

#### INFORMATION

The scale values of the 13E200 electrode do not match the scale values of the 13E125 electrode. Therefore, do not adjust the electrodes by comparing positions of the two electrodes.

## 7 Legal information

### 7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 7.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

## 7.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

## 8 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing

General information	13E200=50	13E200=60
Frequency	50 Hz	60 Hz
Product service life	5 years	
Power supply	All batteries compatible with the Ottobock Myosystem	
Frequency bandwidth	90 - 450 Hz	
Sensitivity range	2,000 - 100,000 x	
Protection rating	IP67	

## 9 Appendices

### 9.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

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