
## 3R55

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	10
<b>FR</b> Instructions d'utilisation .....	15
<b>IT</b> Istruzioni per l'uso .....	20
<b>ES</b> Instrucciones de uso .....	25
<b>PT</b> Manual de utilização .....	30
<b>NL</b> Gebruiksaanwijzing .....	36
<b>SV</b> Bruksanvisning .....	41
<b>TR</b> Kullanma talimatı .....	45
<b>EL</b> Οδηγίες χρήσης .....	50
<b>RU</b> Руководство по применению .....	56
<b>ZH</b> 使用说明书 .....	61

● 4D19

Einzelteile-Pack  
 Single-Component Pack  
 Kit de pièces de rechange  
 Imballo dei singoli componenti  
 Kit componentes  
 Onderdelenpakket  
 Service-Set  
 Conjunto de peças de reposto  
 Σετ μεμονωμένων εξαρτημάτων  
 Комплект деталей  
 维修组件包



▲




Mindestmenge  
 Minimum order quantity  
 Quantité minimum  
 Quantità minima  
 Cantidad mínima  
 Minimikvantitet  
 Minimum aantel  
 Quantidade minima  
 Ελάχιστη ποσότητα  
 Минимальное количество  
 最低起订量

■

Einzelteile  
 Single components  
 Pièces à l'unité  
 Singoli componenti  
 Componentes  
 Separata delar  
 Onderdelen  
 Componentes  
 Μεμονωμένα εξαρτήματα  
 Отдельная деталь  
 维修件

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Last update: 2020-04-15

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document in a safe place.

## 1 Scope of delivery (Fig. 1)

### 1.1 4D19 Single-component pack for 3R55 (Fig. 1, ●)

Consisting of: 2 knee stops, 1 damper protection, 2 small attachment nipples, 1 large attachment nipple, 4 slotted bushings, 4 Belleville spring washers, 2 lock rings, 2 lock nuts.

### 1.2 Minimum order quantity (Fig. 1, ▲)

- |  |                             |
|--|-----------------------------|
| (6) 513T4=18.8x10.2 Belleville Spring Washer | (10) 501T15=M6x40 Cap Screw |
| (7) 4B86 Knee Axis Bushing                   | (11) 507S16=Ø10x1 Lock Ring |
| (8) 4B82 Slotted Bushing                     | (12) 502S19=M6 Lock Nut     |
| (9) 507U16=Ø6.4 Ball Washer                  |                             |

### 1.2 Single components (Fig. 1, ■)

- |   |  |
|---|--|
| (1) 4G107 Lower Joint Section             | (13) 4Z43=H Knee Stop                    |
| (2) 4G157 Upper Joint Section             | (14) 709S18 Socket Wrench                |
| (3) 4G129 Posterior Linkage Bar           | (15) 4V120 Hydraulic Swing Phase Control |
| (4) 4G130 Posterior Linkage Bar, complete | (16) 4G147 Damper Protection             |
| (5) 4G104 Anterior Linkage Bar            | (17) 4G148=9 Attachment Nipple           |
|   | (18) 4G148=11 Attachment Nipple          |
|   | (19) 4X16 Adjustment Aid                 |


## 2 Description

### 2.1 Application purpose

The 3R55 Modular Polycentric Knee Joint is to be used **exclusively** for the prosthetic fitting of lower limb amputees.

### 2.2 Field of application

Field of application according to the **Ottobock Mobility System MOBIS**:

Recommended for amputees with **mobility grade 3 and 4** (non-restricted outdoor walkers and non-restricted outdoor walkers with especially rigorous demands).

**Approved for a patient weight of up to 125 kg / 275 lbs.**

By means of the Hydraulic Swing Phase Control the gait pattern has been optimised, especially with active patients.

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

**Attention!**

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

**Please inform your patients.**

**2.3 Function**

Upper Joint Section (2) with adjustment pyramid and Lower Joint Section (1) of the four-axes joint are connected by Posterior (3 + 4) and Anterior (5) Linkage Bars. Due to its polycentricity the joint provides a rotational and translational movement. The centre of rotation (instantaneous centre of rotation) changes its position dependent on the degree of flexion. In the extended position the centre of rotation is placed at the height of the femoral condyles. It moves downward and forward with increased flexion.

The small integrated Hydraulic Swing Phase Control (15) provides resistance to motion, thus enabling the lower leg to move harmonically and naturally. Other than mechanical gait controls, the hydraulic system automatically accommodates to the different walking speeds. The amount of resistance to motion (damping resistance) can be regulated by adjusting the size of the orifices through which the fluid flows. The movements of flexion and extension may be separately adjusted (refer to section 5).

**2.4 Service life**

**Attention!**

**Reuse on another patient.**

Risk of injury due to loss of functionality as well as damage to the product.

**Only use the product for a single patient.**

As a basic principle, all Ottobock knee joints are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

**3 Technical data**

Article number	3R55
Proximal connection	Adjustment pyramid
Distal connection	Adjustment pyramid
Knee flexion angle	110°
Weight	720 g
System height	90 mm
Proximal system height up to alignment reference point	9 mm
Distal system height up to alignment reference point	81 mm
Max. user weight	125 kg / 275 lbs
Mobility Grade	3 + 4

## 4 Handling

### 4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis.

The advantages of the 3R55 Knee Joint can only be made optimal use of in case of a correct alignment.

**The optimal residual limb position must be anticipated when positioning the socket connector.**

Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

**To align the prosthesis please proceed in two steps:**

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimisation.

#### 4.1.1 Bench alignment with alignment tool (the following steps refer to Fig. 5)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run through the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee-ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter.
- ⑤ Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

#### **Attention!**

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. If required, use the 4R118 Adapter Plate for optimal positioning of the prosthetic knee joint.

- ⑦ Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

#### 4.1.2 Static alignment optimization using L.A.S.A.R. Posture 743L100

(the following steps refer to Fig. 6)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **35 mm anterior to the upper front knee axis** (see Fig. 6).
- ③ After step 2, dynamic optimisation can take place between parallel bars.

## 4.2 Combination possibilities

### Information

When fabricating a prosthesis for the lower limbs, all of the prosthetic components used have to meet the requirements in regard to the patient's body weight and activity level.

## 4.3 Adjustments and final assembly

### 4.3.1 Adjusting the bearing

The axis pins of the linkage bars are mounted in bushings and are hingedly connected to the Upper and Lower Joint Sections. This type of bearing and assembly allows increase and decrease of the friction as well as readjustment of tolerance variation due to wear. To this end fix Lock Nut (12) using enclosed 709S18 Socket Wrench (14) and tighten Cap Screw (10) with 5 mm Allen wrench (Fig. 2).

**Adjustment to eliminate clearance must take place after 2 to 4 weeks walking. Every time maintenance work is carried out, the adjustment should be controlled and the 4Z43=H Knee Stops (13) replaced. For lubrication of the Knee Axis Bushings (7) and Slotted Bushings (8) use 633G6 Ottobock Special Lubricant.**

### 4.3.2 Adjusting resistance to motion

Resistance to flexion and extension are separately adjusted by moving the pins within the horizontally arranged slots at the front and back side of the piston rod in the Hydraulic Swing Phase Control (15). To this end the enclosed Adjustment Aid (19) should be used.

Moving to + = resistance increases

Moving to - = resistance decreases

### 4.3.3 Adjusting the flexion resistance (Fig. 3)

Moving the posterior pin (at the flexion side) to the **right** = resistance **decreases** = easy running flexion.

Moving to the **left** = resistance **increases** = tight running flexion.

### 4.3.4 Adjusting the extension resistance (Fig. 4)

Flex the joint, so that the anterior pin becomes accessible.

Moving the pin to the **right** = resistance **decreases** = easy running extension.

Moving the pin to the **left** = resistance **increases** = tight running extension.

### Attention!

The elastic cover of the Hydraulic Unit must be protected from mechanical damage. The Hydraulic Unit should not be removed, as its position within the joint has been accurately adjusted. In case of malfunctions, please return the complete joint.

## 4.4 Cosmetic foam cover

Use foam cover 3R6, 3S106 or 3R56 for joint 3R55 .

### Attention!

To eliminate noise in the cosmetic foam cover, use 519L5 Silicone Spray. **Do not use talcum powder!** Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. **Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.**

### **Useful information**

To optimise sliding and for eliminating noises, please apply 519L5 Silicone Spray directly on the friction surfaces of the cosmetic foam covers.

#### **4.5 Maintenance information**

Ottobock recommends the readjustment of the knee joint's settings after the patient has spent a period of time getting used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the wear and functionality of the knee joint at least once a year and make follow-up adjustments if need be. Pay special attention to resistance to movement and the development of unusual sounds. Complete flexion and extension must be guaranteed.

We recommend carrying out regular safety checks once a year.

#### **Attention! Please inform your patients.**

Due to different environmental conditions and conditions of application, the function of the knee joint can be impaired. To avoid the risk of accident and possible injury to the patient, the knee joint must no longer be used after a noticeable loss of function has occurred. This noticeable loss of function can be for example stiffness, non-attainment of the extension stop, decreasing swing phase control or stance phase stability, abnormal noises, etc.

#### **Measure**

Have the knee joint examined at a specialist workshop.

## **5 Disposal**

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## **6 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **6.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **6.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website

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