
3R60, 3R60=ST, 3R60=KD, 3R60=HD

DE Gebrauchsanweisung	6
EN Instructions for use	16
FR Instructions d'utilisation	25
IT Istruzioni per l'uso	36
ES Instrucciones de uso	46
PT Manual de utilização	57
NL Gebruiksaanwijzing	66
SV Bruksanvisning	77
DA Brugsanvisning	86
NO Bruksanvisning	95
FI Käyttöohje	104
PL Instrukcja użytkownika	113
HU Használati utasítás	123
CS Návod k použití	132
RO Instrucțiuni de utilizare	141
HR Upute za uporabu	152
TR Kullanma talimatı	162
EL Οδηγίες χρήσης	171
RU Руководство по применению	182
JA 取扱説明書	193
ZH 使用说明书	202








DE Einzelteile
 EN Single components
 FR Pièces à l'unité
 IT Singoli componenti
 ES Componentes
 PT Componentes
 NL Onderdelen
 SV Separata delar
 DA Komponenter
 NO Enkle komponenter
 PL Podzespoły
 HU Alkatrészek
 CS Jednotlivé díly
 RO Componente individuale
 HR Pojedinačni dijelovi
 TR Müferit parçalar
 EL Μεμονωμένα εξαρτήματα
 RU Отдельные детали
 JA 1個から発注いただける部品
 ZH 单个部件




▲ Mindestmenge
 Minimum order quantity
 Quantité minimum
 Quantità minima
 Cantidad mínima
 Quantidade mínima
 Minimikvantitet
 Minimum oantal
 Mindeste mængde
 Minste ordremenge
 Ilość minimalna
 Minimum mennyiség
 Minimální množství pro objednání
 Cantitate minimă
 Minimalna količina
 Minimum miktar
 Ελάχιστη ποσότητα
 Минимальное количество
 複数入パックで発注いただく部品
 最低起订量

1

INFORMATION

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INFORMATION

Last update: 2020-08-31

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document in a safe place.

1 Single components (Fig. 1) ■

- (1) 710H10=2×3 Adjustment Wrench
- (2) 4Z59=5.5×6 Extension Stop ▲
- (3) 4G444=3R60-NB1 Signature Plate
- (4) 4G375=N Protective Knee Cap
- (5) 4G70 Lamination Anchor with
 - 4 pcs of 501T1=M5×16 Truss Head Screw (a)
 - 4 pcs of 502R1=M5×16 Round Nut with two holes (b) ▲
 - 2 pcs of 506G3=M8×12-V Set Screw (c) ▲
 - 1 pc. of 506G3=M8×10 Set Screw (d) ▲


2 Description

2.1 Indications for use

The 3R60 Modular **EBS** Knee Joint is to be **exclusively** used for the prosthetic fitting of amputations of the lower limb. It is a polycentric joint offering an **E**rgonomically **B**alanced **S**tride and hydraulic swing phase control.

2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

Recommended for amputees with **Mobility Grades 2 and 3** (restricted outdoor walkers, non-restricted outdoor walkers).

Approved for a patient weight of up to **125 kg/275 lbs.**

Caution!

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

2.3 Lifetime

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. The maximum lifetime is 5 years.

2.4 Environmental conditions

Caution!

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table „Unallowable environmental conditions“ in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Non-condensing)

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

2.5 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, environmental conditions and field of application.

2.6 Construction and function

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

The amputee will experience considerably greater comfort and stability during stance and swing phase with the proven joint.

The unique characteristic of the 3R60 is its **Ergonomically Balanced Stride (EBS)** which has proven highly useful during 10 years of practical application. In conjunction with its polycentric joint kinematics, it provides notable stance phase stability. **At heel strike a cushioned stance flexion of up to 15° takes place without initiating normal swing flexion. This stance phase**

flexion means shock-absorbing comfort at heel strike as well as additional stability by shifting the centre of rotation posteriorly and proximally. **All in all, walking with the prosthesis becomes more comfortable and much more physiological.**

Two axial forks connect the upper and lower joint section to one another and create a kinematic chain. A pivoting mount (Fig. 2, A) joins the lower joint section to the rear axial fork, which is also connected to the shock absorbing elements of the EBS unit (Fig. 2, C).

At heel strike, the proximal joint sections swing posteriorly around the lower axes (Fig. 4, D). The elastomer bumper (Fig. 2, B) of the EBS unit is compressed as the pivoting mount rotates (Fig. 3a/3b, Fig. 4, E). This pivoting mount serves as a visual control for measuring the efficiency of the EBS function. The resistance of the EBS unit is continuously adjustable and can be set for the weight and activity level of the patient.

The multiaxial action of the joint is composed of rotary and gliding movements. The centre of rotation (instantaneous centre of rotation) varies with the flexion position (Fig. 4, F).

An individually adjustable hydraulic unit located between the axial forks controls swing phase. It prevents excessive heel rise and dampens extension before initial impact. The hydraulic unit's resistances are independently adjustable.

For optimal function, in addition to proper alignment and adjustments, the patient must be taught to use the stance flexion feature. The comfort provided by the knee is most noticeable when walking long distances or on uneven surfaces. Initially, the patient may be unaccustomed to the cushioned knee flexion provided by the EBS unit at heel strike. With proper instruction and experience, the amputee will learn to take considerable advantage of this more natural gait pattern.

2.7 Different connection systems

The EBS Knee Joint is available in 4 different versions. They differ from each other solely by the type of proximal connection used (see front page):

Article number	Version
3R60	Knee Joint with adjustment pyramid connector
3R60=ST	Version with threaded connector for long residual limbs
3R60=KD	Version with lamination anchor for knee disarticulation amputees
3R60=HD	Version with adjustment pyramid connector anteriorly angled by 10° for hip disarticulation amputees

Attention!

For the prosthetic fitting of hip disarticulation amputees the specially adapted 3R60=HD joint version must be used.

3 Technical data

Article number	3R60	3R60=KD	3R60=ST	3R60=HD
Proximal connection	Adjustment pyramid	Lamination anchor	Threaded connector	Adjustment pyramid anteriorly angled by 10°
Distal connection	Adjustment pyramid			
Knee flexion angle	150°	150°	150°	150°
Weight	845 g	940 g	845 g	880 g
System height	171 mm	193 mm	189 mm	174 mm
Proximal system height up to alignment reference point	-2 mm	20 mm	16 mm	1 mm
Distal system height up to alignment reference point	173 mm			
Max. user weight	125 kg/ 275 lbs			
Mobility Grade	2, 3			

4 Handling

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position of the axes. In extension, the instantaneous centre of rotation is located above the adjustment pyramid and behind the alignment reference line, thus providing knee stability in the stance phase (Fig. 4). The advantages of the EBS Knee Joint can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly; PROS.A. Assembly 743A200.
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimization.

4.1.1 Bench Alignment with Alignment Tool (the following steps refer to Fig. 12)

- ❶ Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ❷ Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ❸ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run through the upper anterior axis (alignment reference point)**. At

that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.

- 4 Connect the foot to the modular knee joint using a tube adapter.
- 5 Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- 6 Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

Attention

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear. If required, use the 4R118 Adapter Plate for optimal positioning of the prosthetic knee joint.

- 7 Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

4.1.2 Static alignment optimization using 743L100 L.A.S.A.R. posture (the following steps refer to Fig. 13)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- 1 To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- 2 Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **10 mm anterior to the lower front knee axis**.
- 3 After step 2, dynamic optimization can take place between parallel bars (see Section 4.3.2).

4.2 Adjustments and final assembly

4.2.1 3R60=KD socket fabrication

4.2.1.1 Lamination prior to trial fitting

- Take a piece of 623T3 Perlon Stockinette twice the length of the plaster model and pull half of it over the isolated socket. Twist the other half of the stockinette on the distal end and reflect it over the socket.
- In order to absorb strong forces in the area of the modular knee joint, the socket must be gradually reinforced with 616G13 Fibreglass Stockinette. To do this, pull over the first layer to 2/3 of the socket length, tie it off and then reflect it up to half of the socket length.
- Put two layers of 616G12 Carbon-Fibre Cloth on the distal area so as to provide an underlying base with a projecting brim of 3 cm around the 4G70 Adapter, which will be mounted later in the procedure. Pull on two layers of 623T3 Perlon Stockinette.
- The socket lamination is performed in two stages. In the first stage the socket is laminated up to 2/3 of its length with 617H19 Orthocryl Lamination Resin. In the second stage the proximal part of the socket is laminated with 617H17 Orthocryl Flexible Resin.
- Once the first casting has hardened, pull on 2 more layers of 623T3 Perlon Stockinette before laminating the proximal part of the socket with 617H17 Orthocryl Flexible Resin.

Attention!

The lamination anchor acts as a flexion stop in the area of the EBS pivoting mount (Fig. 10). To provide this stop function, and in order to avoid damage to the knee joint, an even contact surface in the area of the EBS pivoting mount must be present after final lamination.

4.2.1.2 Applying the lamination anchor

- The correctly fitted lamination anchor must be prepared before it is glued to the socket.
- To do this, fill the coupling opening with 636K8 Plastaband. Clamp the socket and joint into the alignment apparatus. Remove any soft, flexible or porous materials from underneath the lamination anchor's supporting surface.

Attention!

Mix the putty exclusively from 617H21 Orthocryl Sealing Resin and 639A1 Talcum Powder. Then glue on the lamination anchor.

- Secure with 627B2 Adhesive Tape for trial fitting. Test the stop function (see Section 4.3.1.1). If necessary, form the required contact surface with putty. It is also possible to glue a stop bumper of Pedilin onto the outer laminate.

4.2.1.3 Socket completion

Once the trial fitting has been completed, attach the lamination anchor to the socket with truss head screws and round nuts with two holes, and then laminate.

Proceed with further reinforcement as follows:

- Two layers of 616G12 Carbon Fibre Cloth must be sandwiched between two layers of 623T3 Perlon Stockinette. To do this, pull one layer of 623T3 Perlon Stockinette over the whole socket and tie it on the top circularly, apply the two layers of 616G12 Carbon Fibre Cloth over the arms of the 4G70 Adapter, and then reflect the remaining half of the 623T3 Perlon Stockinette over the carbon fibre layers.
- Now reinforce once again with 616G13 Fibreglass Stockinette (as described in Section 4.3.1.1).
- Finally, pull on two layers of 623T3 Perlon Stockinette. Proceed with lamination as performed in the first stage of lamination.

Attention!

Deviating from the processing instructions and the recommended materials for the lamination of the anchor can lead to a loosening or breaking of the adapter.

Attention!

Upon completion of the assembly tighten the adjustment screws with 710D1 Torque Wrench. Torque: **15 Nm (11 ft.lbf/133 in.lbf)**. Once you have attained the final setting for the prosthesis, use 636K13 Loctite to secure the set screws.

4.2.2 Adjusting the knee joint during trial walking

The user will notice that the EBS joint is different from other prosthetic knees because of the Ergonomically Balanced Stride function. With EBS, the user can walk with a natural gait, allowing the knee to slightly flex during stance phase without risk of falling (Fig. 3a + 3b). He or she must learn to avoid the typical extension of the residual limb and encourage the cushioned flexion through compression of the elastomer bumper.

As with other polycentric joints, swing phase flexion occurs when the weight is on the ball of the foot. Alignment, foot design, condition of the residual limb, and the amputee's activity level all

influence knee function. The pendulum effect from the weight of the foot and shin length are additional variables. Initial dynamic alignment should begin from the base setting and the function of the knee should be explained to the patient.

Before changing the factory settings, the following guidelines must be observed:

- For changing the degree of stance flexion, first correct the position (see Section 4.3.2.1), then re-adjust the EBS unit (see Section 4.3.2.2).
- For swing phase, adjust flexion resistance first and then adjust extension resistance (see Section 4.3.2.3).

The knee is factory adjusted to a basic position for initial dynamic alignment

4.2.2.1 Stance phase stability by positioning the joint (Fig. 5 + 6)

Contrary to single-axis knee joints, polycentric knee joints are stable at heel strike. In the EBS Knee Joint, the compression of the EBS unit adds additional stance phase stability. At the factory, the EBS unit is adjusted to provide low resistance.

Initiation of swing flexion depends on the position of the joint and thus of the location of the instantaneous centre of rotation. Tipping the joint in the sagittal plane (changing the angle using the adjustment screws) can be used to determine the location of the centre of rotation and thus to vary the degree of stability.

Stance phase stability too great (flexion can only be initiated with a great deal of effort)	=	instantaneous centre of rotation (ICR) too far posterior	➔	Tip the joint forward using the adjustment screws.	= Fig. 5
Insufficient stance phase stability (knee joint buckles)	=	instantaneous centre of rotation (ICR) too far anterior	➔	Tip the joint backward using the adjustment screws.	= Fig. 6

Any adjustments performed above the joint must be corrected by corresponding counter-adjustments below the joint and on the foot, in order to maintain the position of the foot, as described in Section 4.1.2. Afterwards, repeat static alignment optimization.

Attention! The =KD and =ST connection systems require another way of proceeding during manufacture of the prosthesis. Stance phase stability and initiation of flexion depend on the position of the lamination anchor. Later readjustment of the prosthesis in the sagittal and frontal planes is not possible!

4.2.2.2 Adjusting the ergonomically balanced stride feature (footfall stability)

The stance flexion resistance may be adjusted by compressing the elastomer bumper with the adjustment nut. To do this, insert 710H10=2x3 Adjustment Wrench provided with the knee into the hole (Fig. 7):

Rotation to the left (in - direction)	=	decreased resistance	=	more knee flexion	=	strong raising of the pivoting mount (Fig. 3b)	=	increased stance stability
Rotation to the right (in + direction)	=	increased resistance	=	less knee flexion	=	slight raising of the pivoting mount (Fig. 3a)	=	decreased stance stability

**Attention!**

The EBS unit has an overturn protection. When increasing the resistance (rotation to the right in + direction), do not forcibly turn beyond the point where it becomes more difficult to turn!

If the maximally adjustable resistance of the EBS unit is insufficient (knee flexion too easy), possibly tip the joint forward (see Section 4.3.2.1), but be aware that knee stability is reduced by this measure.

Notice:

The round spring of the EBS unit must be adjusted so that the pivoting mount does not permanently move into end position.

4.2.2.3 Adjusting the swing phase

The hydraulic swing phase control ensures a smooth gait. Movement resistances prevent the prosthetic lower leg from swinging too far when flexing and guarantee dampened extension. Along with the patient's habits, the foot weight and the length of the lower leg, which have a pendular weight, also have an influence on his or her gait.

Before changing the factory settings, the following guidelines must be observed:

- Upon delivery, the valve screws are set to the lowest resistance for extension damping. For flexion damping, a medium setting has been chosen.
- The flexion and extension resistances can be independently adjusted using the 710H10=2x3 Adjustment Wrench provided.

Adjusting flexion (F) on the left joint side (from the posterior direction of view, see Fig. 8)

Rotation of the valve screw of the hydraulic unit to the right (+)	=	greater resistance	=	less heel rise
Rotation of the valve screw of the hydraulic unit to the left (-)	=	less resistance	=	more heel rise

Adjusting extension (E) on the right joint side (from the posterior direction of view, see Fig. 9)

Rotation of the valve screw of the hydraulic unit to the right (+)	=	greater resistance	=	more extension damping
Rotation of the valve screw of the hydraulic unit to the left (-)	=	less resistance	=	less extension damping

**Attention!**

Proceed **with caution and in small increments** when adjusting the resistances, since there is the risk of falling!

When making extension damping adjustments, make sure that the leg still completely extends with each step. When adjusting the swing phase control, the influence of a cosmetic foam cover must be compensated for.

**Attention!**

In extreme cases of dynamic permanent loading, the hydraulic damping systems can heat up. When using the knee joint do not reach into the joint mechanism since there is the risk of injury.

4.3 Cosmetic cover

Combination possibilities if foam covers are used optionally	
3R60	3S107
3R60=KD	3S107, 6R6
3R60=ST	3S107
3R60=HD	3S27

Attention!

Do not use talcum powder to eliminate noise in the cosmetic foam cover. Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. Using this medical product after application of talcum powder will render all claims null and void.

Note:

To reduce friction and to eliminate noise, please apply 519L5 Silicone Spray directly on the contact surfaces of the cosmetic foam cover.

4.4 Maintenance Instructions

Attention - Please inform your patients!

The knee joint function may be affected by both environmental and service conditions. To reduce the risk for the patient, the knee joint should not be used when noticeable functional changes occur. Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

Measure after noticeable functional changes have occurred:

Have the prosthesis examined at a specialist workshop.

Attention!

The knee joint must not be lubricated or greased, since this would result in a risk of bearing damage and loss of function.

Attention!

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the knee joint in for service.

Ottobock recommends readjusting the knee joint's settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the knee joint for wear and proper function at least once a year and make adjustments if necessary. Special attention should be paid to the flexion and extension resistances, bearings, and abnormal noises. Full flexion and extension must be guaranteed.

This component has been tested according to ISO 10328 standard for three million load cycles. Depending on the amputee's activity this corresponds to a duration of use of three to five years.

We recommend carrying out regular yearly safety checks.

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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