




17K43, 17K45, 17K46, 17K47

DE	Gebrauchsanweisung	5
EN	Instructions for use	11
FR	Instructions d'utilisation	17
IT	Istruzioni per l'uso	24
ES	Instrucciones de uso	30
PT	Manual de utilização	37
NL	Gebruiksaanwijzing	44
SV	Bruksanvisning	50
DA	Brugsanvisning	56
NO	Bruksanvisning	62
PL	Instrukcja użytkowania	68
HU	Használati utasítás	75
CS	Návod k použití	82
TR	Kullanma talimatı	88
EL	Οδηγίες χρήσης	94
RU	Руководство по применению	101


5

6

INFORMATION

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INFORMATION

Date of last update: 2021-05-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17K43, 17K45, 17K46 and 17K47 orthotic knee joints.

2 Product description

2.1 Available sizes

The 17K43, 17K45, 17K46 and 17K47 orthotic joints made of stainless steel are available in one size.

2.2 Components/design

Scope of delivery for 17K43, 17K47 (see fig. 1, see fig. 4)			
Item	Quantity [pcs.]	Designation	Reference number
1	1	Joint bar, upper section	

Scope of delivery for 17K43, 17K47 (see fig. 1, see fig. 4)			
Item	Quantity [pcs.]	Designation	Reference number
2	2	Bushing	17Y17=6x8x2
3	1	Joint centre piece	17Y127
4	3	Oval head screw	501S42=M3x6
5	2	Joint screw	501S32=M4x12x8
6	1	Joint bar, lower section	
7	1	Extension stop, 6 °	17Y90=6
8	1	Spacer bushing	17Y17=3.1x6x2.1
9	2	Bearing nut	17Y93=6x5xM4

Scope of delivery for 17K45, 17K46 (see fig. 2, see fig. 3)			
Item	Quantity [pcs.]	Designation	Reference number
1	1	Joint bar, upper section	
2	1	Flexion stop, 100 °	17Y91=100
3	2	Bushing	17Y17=6x8x2
4	1	Joint centre piece	17Y92
5	4	Oval head screw	501S42=M3x6
6	2	Joint screw	501S32=M4x12x8
7	1	Joint bar, lower section	
8	1	Extension stop, 6 °	17Y90=6
9	2	Bearing nut	17Y93=6x5xM4

Scope of delivery for 17K43, 17K47 (see fig. 5)			
Item	Quantity [pcs.]	Designation	Reference number
1	3	Extension stop	17Y90=6 (pre-assembled) 17Y90=10 17Y90=20 17Y90=30

Scope of delivery for 17K45, 17K46, (see fig. 5)			
Item	Quantity [pcs.]	Designation	Reference number
1	3	Extension stop	17Y90=6 (pre-assembled) 17Y90=10

Scope of delivery for 17K45, 17K46, (see fig. 5)			
Item	Quantity [pcs.]	Designation	Reference number
			17Y90=20 17Y90=30
2	3	Flexion stop	17Y91=0 17Y91=60 17Y91=80 17Y91=100 (pre-assembled)

3 Intended use

3.1 Indications for use

The orthotic knee joints are to be used **exclusively** in pairs for the orthotic treatment of the lower limbs in knee orthoses or KAFOs.

3.2 Indications

- Orthopaedic diseases of the lower limbs
- Indications must be determined by the physician.

3.3 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Safety instructions

 CAUTION
Excessive strain on load-bearing components
Injuries due to changes in or loss of functionality
▶ Only use the product for the defined area of application.

- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).

- Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

CAUTION

Repeated bending at the same position

Risk of injury due to breakage of the bar, change in or loss of functionality

- Avoid repeated bending at the same position.

INFORMATION

Increased wear and tear

The orthotic knee joints are intended only for use in combination with free-motion or foot lifting orthotic ankle joints. Wear occurs more rapidly in combinations involving orthotic ankle joints with a dorsal stop.

INFORMATION

Parallel alignment of the orthotic joints

Use the 743R6 alignment fixture to position the orthotic joints on the plaster positive.

Shaping the bars

- > **Required tools:** Vice, 2x 711S* bending irons, grinding machine, screwdriver, Loctite 241
- 1) Loosen the joint screw and disassemble the joint.
 - 2) Use the 711S* bending irons for contouring the bar.
 - 3) Clamp a bending iron in the vice and position the bar.
 - 4) Adapt the bar with the second bending iron.
 - 5) Contour the bar until the desired result is obtained.
 - 6) Observe the minimum bending radius of **10 mm** (see fig. 6).
 - 7) Remove grooves and burrs by grinding.
 - 8) For final assembly, secure the joint screw with Loctite 241.

Replacing the stop

The 6 ° extension stop and the 100 ° flexion stop are pre-installed on the product on delivery.

- 1) Loosen the screw and remove it from the joint centre piece.
- 2) Disassemble the stop.
- 3) Insert the new stop into the joint centre piece.

- 4) Install the screw in the joint centre piece.

Improving corrosion protection

Polish or sinter the bars for improved corrosion protection. Ottobock recommends the 618T40* sintering powder for surface coating.

- ▶ Do not heat the material for more than **5 minutes** at max. **150 °C**.

Using protective devices

- ▶ Use joint protectors to protect the product against soiling.

6 Cleaning

After contact with water containing salt, chlorine or soap, or if they get dirty, the orthotic joints must be promptly cleaned.

- 1) Rinse the orthotic joints with clean fresh water.
- 2) Dry the orthotic joints with a cloth or allow them to air dry. Avoid exposure to direct heat (e.g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

Only use 633F7 special lubricant.

Replacing the bushing

- 1) Loosen the joint screws and remove the bearing nuts.
- 2) Loosen and remove the oval head screws.
- 3) Disassemble the joint centre piece.
- 4) Remove the bushing from the bar and replace it with a new bushing.
- 5) For final assembly, secure the screw connection with Loctite 241.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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