

3K9

DE Gebrauchsanweisung	4
EN Instructions for use	6
FR Instructions d'utilisation	8
IT Istruzioni per l'uso	10
ES Instrucciones de uso	12
SV Bruksanvisning	14
NL Gebruiksaanwijzing	16
PT Manual de utilização	18
HU Használati utasítás.....	20
EL Οδηγίες χρήσης.....	22

INFORMATION

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INFORMATION

Last update: 2017-03-31

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1 Application

English

⚠ CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- Only use the product for a single patient.

The 3K9 Plastic Single Axis Knee-Shin Component is to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

2 Parts

- (1) 4K11 Knee Part of Pedilen
- (2) 4K10 Knee Joint Section of Plastic with Manual Lock
- (3) 4K12 Shin Part of Pedilen, calf circumference: 38 cm
- (4) 4F17=N Lock Handle
- (5) 13Z2=2.3 Lock Washer
- (6) 513D25 Compression Spring
- (7) 4F12 Guide Piece for Lock

3 Description

The 3K9 Knee-Shin Component has been designed for the fitting of geriatric amputees, and it has proved useful for this group of patients. Attend to careful selection of patients.

The components used in the known lamination and Pedilen technique enable fabrication of light-weight and thin-walled prostheses with sufficient stability and best possible knee security.

Special attention must be paid to exact processing.

4 Processing Instructions

The 4K10 Knee Joint is made of thermoplastic material offering shape stability. The material must be treated according to its requirements; for example, petroleum benzine only should be used for cleaning.

When processing please follow the items below:

4.1 Before Trial Walking

Bond the knee joint and the Pedilen parts with 636K7 Orthocryl Putty (fig. 1).

Hold the pasted seams in position by means of cramps (bore holes: 4 mm) and secure with linen adhesive tape (fig. 2).

4.2 After Trial Walking

Finishing must comprise reinforcement lamination of the upper and lower joint section.

The following reinforcement materials are recommended:

616G13 Fiberglass Woven Stockinette

623T9 Nylglass Stockinette, and
616G12 Carbon Fiber Webbing.

Detach the joint, loosen the pasted seams thermoplastically (too high temperatures must be avoided), remove putty (fig. 3).

Before reinforcement takes place, carefully seal the openings of the joint parts (fig. 4).

Apply sealing resin to the fillings of plastilin and plaster (fig. 5).

Separate the Pedilen parts with PVA sheeting. Fasten the joint parts with Orthocryl putty (fig. 6).

Tie off the fiberglass woven stockinette at the upper joint part with roving; invert the bag and tie off under tension (fig. 7).

Pull over a seamed layer of nylglass stockinette (fig. 8).

Reinforce as well the lower part of the leg up to the joint part using fiberglass woven stockinette (fig. 9).

Between the two layers of nylglass stockinette circumferentially bond a carbon fiber strip – about 50 mm wide.

After setting of the lamination resin, make sure the laminate is even, if required, smooth and laminate again with a layer of nylglass stockinette (fig. 10).

NOTICE

The structural strength of the component created by the laminate offers sufficient stability for normal weight bearing especially within the scope of geriatric fitting. The prosthesis should not be used with patients requiring a prosthesis for heavy weight bearing.

4.3 Instructions for overriding the Locking mechanism

In order to secure the knee joint at any time, the locking mechanism is provided with an automatic catch.

The guide groove of the 4F17=N Lock Handle is designed to prevent the patient from over-riding the automatic locking mechanism.

For patients who are able to walk safely without knee joint lock, the prosthetist, however, can permit an override by fully opening the groove in the lock handle.

When the knee is used in this way, a clicking noise may be heard at terminal impact. It is caused by the 4F12 Guide Piece. The noise can be eliminated by removing the 4F12 Guide Piece, as follows: Separate the knee from the shin. Remove the 13Z2 Lock Washer with a screw driver (fig. 11).

Pull out the 513D25 Compression Spring (fig. 12).

Remove the 4F12 Guide Piece with a screw driver (fig. 13).

Remove the tab in the groove of the 4F17=N Lock Handle (fig. 14).

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH