






2R45=S, 2R45=34

DE	Gebrauchsanweisung	4
EN	Instructions for Use.....	10
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização.....	34
NL	Gebruiksaanwijzing.....	40
SV	Bruksanvisning	46
DA	Brugsanvisning.....	52
HU	Használati utasítás	58
CS	Návod k použití.....	64
EL	Οδηγίες χρήσης	69

INFORMATION

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INFORMATION

Last update: 2020-12-14

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

1 Description

The Length Adjustable Tube Adapter (2R45=S/2R45=34) facilitates the correct static alignment of lower limb prostheses during the trial fitting. Reproducible adjustments make it possible to optimize gait during the trial walking. The length adjustable tube adapter is mounted between the prosthetic foot and the knee joint. Tube lengths are adapted for use between the prosthetic foot and the knee joint. The total length of the tube adapter can be adjusted with a tube clamp. The outer and inner rotation of the foot is also adjustable. **2R45=34:** The graduated scales ensure reproducibility of adjustment (Fig. 1).

2 Application and Description

The Length Adjustable Tube Adapter (2R45=S/2R45=34) is to be used **exclusively** for the initial fitting or interim prostheses, testing purposes and trial fitting of lower limb prostheses. It is used as an adjusting element for alignment optimization and must be disassembled before completion of the definitive prosthesis.

Field of application according to the **Ottobock MOBIS® Mobility System:**

2R45=S:

Approved for patient weight **up to 100 kgs (220 lbs)**.

2R45=34:

Approved for patient weight **up to 125 kgs (275 lbs)**.

3 Safety

3.1 Explanation of warning symbols

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

1. Use the product according to the specified area of application (Area of application).

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

1. Only combine the product with prosthetic components that are approved for that purpose.
2. Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

1. Do not expose the product to unallowable environmental conditions.
2. If the product has been exposed to unallowable environmental conditions, check it for damage.
3. If damage is apparent or in case of doubt, do not continue using the product.
4. Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

1. Ensure that the approved service life is not exceeded.

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

1. Clean the threads before every installation.
2. Apply the specified torque values.
3. Follow the instructions regarding the length of the screws and about how to secure the screws.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

1. Use caution when working with the product.
2. If the product is damaged, check it for proper function and readiness for use.
3. In case of changes in or loss of functionality, do not continue using the product (see „Signs of changes in or loss of functionality during use“ in this section).
4. Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Technical Data

2R45=S:

Weight of the tube adapter: 200 g

Max. body weight: 100 kg

2R45=34:

Weight of the tube adapter with the 2R56=230 Tube with Scale:
427 g (15.06 oz)

Weight of the tube adapter with the 2R56=300 Tube with Scale:
466 g (16.44 oz)

Total length of the tube adapter with 2R56=230: min. 250 mm/ max.
310 mm, Ø 34 mm (min. 9.84 in/max. 12.2 in, Ø 1.36 in)

Total length of the tube adapter with 2R56=300: min. 320 mm/ max.
380 mm, Ø 34 mm (min. 12.60 in/max. 14.96 in, Ø 1.36 in)

Max. body weight: 125 kgs (275 lbs.)

5 Assembly / Trial Fitting / Disassembly

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

1. Observe the alignment and assembly instructions.

Caution!

Examine the parts for any damage or deformation before re-using them. Damaged tubes and tubes which have been in use for more than 5 years must be replaced. The year of manufacture code is indicated on the tube (Fig. 3).

The Length Adjustable Tube Adapter (2R45=S/2R45=34) has modular connections and is ready for installation (distal = pyramid receiver with set screw, proximal = 2R45=S: 30mm tube/2R45=34: 34-mm tube). The instructions for static alignment according to instructions for use for the specific components should be taken into consideration. It is advisable to set the central vertical scaling line centrally in the tube adapter slot during the static alignment.

Alignment is optimized during trial fitting by making adjustments to the adapter in length and rotation position. **2R45=34:** The graduated scales allow for a precise adjustment (Fig. 1).

Fine changes in the alignment can be performed at any time with 506G3 Adjustment Set Screws. When changing a module or in case of disassembly, the previously determined position can be maintained, if two adjacent set screws (specifically those most deeply screwed in) are removed.

Caution!

2R45=34: Pulling out the tube above the STOP mark and shortening the tube is prohibited (see Fig. 4).

Tighten the socket head cap screw of the tube clamp with a torque of **2R45=S:** 13 Nm/**2R45=34:** 10 Nm and the 506G3 Adjustment Set Screws with a torque of 15 Nm (see Fig. 2).

Once you have finished the trial fitting, mark the height and perpendicular lines. For replacement of the length adjustable tube adapter by modu-

lar components according to MOBIS® the 743A16 Transfer Apparatus is recommended.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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