
3R33, 3R17

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	14
IT Istruzioni per l'uso	19
ES Instrucciones de uso	24
PT Manual de utilização	28
NL Gebruiksaanwijzing	34
SV Bruksanvisning	39
PL Instrukcja użytkowania	44
HU Használati utasítás	49
RO Instrucțiuni de utilizare	54
HR Upute za uporabu	59
TR Kullanma talimatı	64
EL Οδηγίες χρήσης	68
JA 取扱説明書	73
ZH 使用说明书	78

● 4D10 Einzelteile-Pack
 ○ 4D11 Single-Component Pack
 Kit de pièces de rechange
 Imballo dei singoli componenti
 Kit componentes
 Conjunto de peças de reposto
 Onderdelenpakket
 Service-Set
 Alkatrész-csomag-hez
 Zestaw komponentów
 Σετ μεμονωμένων εξαρτημάτων
 组件包

▲ Mindestmenge
 Minimum order quantity
 Quantité minimum
 Quantità minima
 Cantidad mínima
 Quantidade minima
 Minimum aantal
 Minimikvantitet
 Minimalna ilość
 Minimum mennyiség
 Ελάχιστη ποσότητα
 最低定货量

■ Einzelteile
 Single components
 Pièces à l'unité
 Singoli componenti
 Componentes
 Componentes
 Onderdelen
 Separata delar
 Zestaw wymienny
 Alkatrészek
 Μεμονωμένα εξαρτήματα
 零件

INFORMATION

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INFORMATION

Last update: 2020-04-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document in a safe place.

1 Components (fig. 1)

(for 3R33 in brackets)

- | | |
|---|---|
| (1) 4G32=2 (4G99=2) Upper Joint Section | (16) 4G35 Lock Cam |
| (2) 4G33 (4G98) Lower Joint Section | (17) 4G36 Lock Lever |
| (3) 4G34 Extension Stop | (18) 4A46 Axis Pin |
| (4) 4H45 Rubber Bumper | (19) 4A47 Axis Pin |
| (5) 506G3=M8×20 Set Screw | (20) 4A51 Axis Pin |
| (6) 2H31=14×12×15 Joint Axis Bushing | (21) 4V104 Spring Guide Housing |
| (7) 4Y10 Thrust Bearing | (22) 513D18=4.7×31 Compression Spring |
| (8) 4A42=46 Axis | (23) 501S22=8×M3.5 Lock Screw |
| (9) 4A43 Nut for Axis | (24) 21A18=2 Perlon Cord |
| (10) 4Z30 Safety Plate | (25) 21A12 Threaded Fitting, short |
| (11) 501S27=M4×5 Oval Head Screw | (26) 21A8 Cable Clamp |
| (12) 513D35=60×7.6×1.2 Compression Spring | (27) 516S3 Button Head Screw, nickel plated |
| (13) 4V102 Spring Guide Housing | (28) 4F17=N Lock Handle, complete |
| (14) 4V101 Spring Guide Pin | (29) 4G66 Plastic Knee Cap |
| (15) 504R1=2.5×20 Guide Pin | (30) 501Z2=M4×8 Cap Screw |

1.1 Single-component packs

4D10 Single-component pack for 3R33 and 3R17 (fig.1, ●)

Consisting of: 1 plastic knee cap, 2 bushings, 2 spring bushings, 1 guide pin, 1 compression spring, 1 set screw, 1 safety plate, 1 oval head screw, 1 rubber bumper, 2 thrust bearings, 2 compression springs, 2 spring guide pins, 2 lock screws, 1 spring guide housing.

4D11 Single-component pack for 3R33 and 3R17 (fig.1, ○)


Consisting of: 1 lock unit, 1 short threaded fitting, 1 cable clamp and perlon cord.

2 Description

2.1 Application and recommendations

The monocentric Modular Knee Joints 3R33 (titanium) and 3R17 (stainless steel) with adjustable manual lock are to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.


Field of application according to the **Ottobock Mobility System MOBIS**

3R33: Recommendation for **mobility degree 1**

(indoor walkers)

Approved for a patient with a weight of up to 125 kg

3R17: Recommendation for **mobility degree 1**

(indoor walkers)

Approved for a patient with a weight of up to 150 kg

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

Due to the small space requirements above the axis it is suitable for long above knee amputation levels. Further, the 3R17 and 3R33 Knee Joints can be interchanged with our 3R20, 3R36, 3R22 and 3R34 Knee Joints.

Attention!

Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

Please inform your patients.

2.2 Environmental conditions

⚠ CAUTION

Environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Non-condensing)

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

2.3 Service life

⚠ CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

This prosthetic component has been tested by the manufacturer for 1 million load cycles. For patients with mobility grade 1, this corresponds to a service life of a maximum of 5 years.

2.4 Construction

The Upper Joint Section (1) is connected to the Lower Joint Section (2) by the Axis (8). Medio-lateral movement is prevented by Ball Thrust Bearings (7). The axis is a hardened, precision ground bolt and is mounted in special bushings. The Nut for Axis (9) is serrated and secured with a Safety Plate (10).

2.5 Extension assist

The Extension Assist consists of two Spring Guide Pins (14), two Spring Guide Housings (13) and two Compression Springs (12). The primary function of this unit is to aid the knee in reaching full extension to help insure that the lock engages.

2.6 Lock mechanism

When the knee is in full extension, the Compression Spring (22) forces the Lock Cam (16) into position against the distal segment of the Upper Joint Section (1). This blocks flexion of the knee.

The adjustable Extension Stop (3) controls hyperextension of the knee.

When the Lock Handle (28) is pulled upward, the Lock Lever (17) rotates the Lock Cam (16) downward, unlocking the knee and permitting it to flex.

On reaching full extension once again, the Compression Spring (22) forces the Lock Cam (16) back into position locking the knee automatically.

If the patient wishes to walk with a free knee, the Lock Handle (28) is designed so that it will hold the lock mechanism disengaged and override the automatic lock feature.


In place of the Lock Handle (28) the 4F34 Lock Control Unit may be used, which is adjustable left or right and for push or pull.

3 Alignment recommendations and assembly

The three-dimensional arrangement of the prosthetic socket and modular elements affects the prosthesis' static and dynamic functions. The position of the residual limb must be taken into account for the connection to the socket. Plumb lines in the frontal and sagittal planes from the hip joint's center of rotation marked during plaster cast taking and trial fitting of the test socket will facilitate correct positioning of lamination anchor or socket adapter.

The prosthetic foot and knee joint are aligned according to the recommended values. Recommended assembly position: The reference alignment line runs through the anterior, upper axis (i. e. **10-15 mm posterior placement of the joint**, see illustration). Function of the joint is influenced by the position of the axes. Thus when providing basic alignment, the two lower axes should create a horizontal line from the lateral aspect.

We recommend using the Ottobock device for alignment L.A.S.A.R. Assembly 743L200 (or the configuration appliance 743A100). The L.A.S.A.R. posture 743L100 has succeeded in optimizing the alignment of the prosthesis for the patient.

Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

Unallowable combinations (7E10)

Unallowable combinations	
Prosthetic hip joints	7E10*

4 Adjusting the extension stop (fig. 2)

The Extension Stop (3) can be adjusted with the Set Screw (5) to control the amount of hyperextension movement in the knee when the lock is engaged.

Flex the Knee Joint until the Set Screw (5) becomes accessible (4 mm Allen wrench).

Turning clockwise = reduces hyperextension movement

Turning counter-clockwise = increases hyperextension movement

Important:

Reducing the hyperextension movement makes it more difficult for the lock to engage.

5 Exchanging the compressing spring (12) (fig. 3)

5.1 Removal

Press the Spring Guide Pin (14) downward using the flat blade of a screw driver, until it disengages from the Axis Pin (20).

Secure the spring with your thumb to prevent it will jump out!

5.2 Installation

Place the notch of the Spring Guide Pin (14) onto the Axis Pin (20). Push the head of the Spring Guide Housing (13) up and into the receptacle in the Lower Joint Section (2).

6 Compensating for mediolateral movement (fig. 4+5)

If mediolateral side play develops, it may be reduced by tightening the Axis Nut (9). Remove the Safety Plate (10) and tighten the Axis Nut (9). Overtightening the nut may restrict normal movement of the Knee Joint. Be sure to install the Safety Plate (10) after the adjustment has been made.

7 Foam cosmetic cover

The 3R24 or 3S124 Foam Cover is used with the 3R17/3R33 Knee Joints.

Attention!

To eliminate noise in the cosmetic foam cover, use 519L5 Silicone Spray. **Do not use talcum powder!** Talcum powder reduces the lubrication of the mechanical parts, which may lead to a malfunction and thus increase the risk of failure. **Using this medical product after application of talcum powder will render all claims against Otto Bock HealthCare null and void.**

Useful information

To optimize sliding and for eliminating noises, please apply 519L5 Silicone Spray directly on the friction surfaces of the foam cosmetic covers.

**Attention! Please inform your patients.**

Due to different environmental conditions and conditions of application, the function of the knee joint can be impaired. To avoid the risk of accident and possible injury to the patient, the knee joint must no longer be used after a noticeable loss of function has occurred. This noticeable loss of function can be for example stiffness, non-attainment of the extension stop, decreasing swing phase or stance phase stability, abnormal noises, etc.

Measure

Have the knee joint examined by a specialist workshop where service parts or, if necessary, the entire knee joint will be exchanged.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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