
## 10R1, 10R2, 10R3, 10R4

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	7
<b>FR</b>	Instructions d'utilisation .....	12
<b>IT</b>	Istruzioni per l'uso .....	16
<b>ES</b>	Instrucciones de uso .....	21
<b>PT</b>	Manual de utilização.....	26
<b>NL</b>	Gebruiksaanwijzing.....	31
<b>SV</b>	Bruksanvisning.....	36
<b>DA</b>	Brugsanvisning.....	40
<b>NO</b>	Bruksanvisning.....	44
<b>PL</b>	Instrukcja użytkowania.....	49
<b>HU</b>	Használati utasítás .....	54
<b>CS</b>	Návod k použití.....	58
<b>RO</b>	Instrucțiuni de utilizare.....	63
<b>TR</b>	Kullanma talimatı.....	68
<b>EL</b>	Οδηγίες χρήσης .....	72
<b>RU</b>	Руководство по применению .....	77
<b>JA</b>	取扱説明書.....	82
<b>ZH</b>	使用说明书.....	86

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## Expalantion of symbols

English

<b>⚠ WARNING</b>	Warnings regarding possible risks of severe accident or injury.
------------------	---

<b>⚠ CAUTION</b>	Warnings regarding possible risks of accident or injury.
------------------	--

<b>NOTICE</b>	Warnings regarding possible technical damage.
---------------	---

<b>INFORMATION</b>	Additional information on the fitting/use.
--------------------	--

### **INFORMATION**

Last update: 2017-02-13

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

The scope of delivery is shown on the front page.

## 1 Single components

For available single components refer to the catalogue.

## 2 Description

### 2.1 Indications for use

The 10R1 / 10R2=\* / 10R3=\* Ottobock Adapters and the 10R4=\* Knurled Plate are to be used **solely** for the exoprosthetic fitting of the upper extremities in combination with Ottobock arm components.

### 2.2 Field of application

The 10R1 / 10R2=\* / 10R3=\* Ottobock Adapters and the 10R4=\* Knurled Plate are used for attachment of an Ottobock hook, Ottobock system hand or passive prosthetic hand to the modular arm prosthesis.

### 2.3 Safety instructions

Please forward the following safety instructions to your patients:

#### **WARNING**

**Risk of accident when driving a motor vehicle.** The ability of prosthesis users to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

## CAUTION

**Risk of injury as a result of excessive strain.** The Ottobock prosthetic components described in this manual have been developed for everyday life activities and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.). Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service.

## NOTICE

**Risk of corrosion.** Please avoid exposing prosthetic components to surroundings that corrode metal parts, for example, freshwater, saltwater, acids and other liquids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

## NOTICE

**Damage due to incorrect environmental conditions.** Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthesis and its components. Failure to follow these instructions can lead to malfunctioning and damage of the prosthesis.

## NOTICE

**Damage caused by improper cleaning agents.** Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

### 2.4 Function

The adapters are available in the following versions:

- 10R1– with threaded stud M12x1.5

- 10R2=M12x1.5/ 10R3=M12x1.5 – with metric inner thread
- 10R2=1/2"-20/ 10R3=1/2"-20 – with inch inner thread
- 10R4=M12x1.5 – with threaded stud M12x1.5 on both sides
- 10R4=1/2"-20 – with threaded stud M12x1.5 on one side; with threaded stud 1/2"-20 on the other side

The 10R1 Adapter is used for attaching an Ottobock prosthesis hand to the Ottobock modular arm prosthesis.

The 10R2=\* Adapter is used for attaching an Ottobock system hand or an Ottobock hook to the Ottobock modular arm prosthesis.

The 10R3=\* Adapter is used for attaching an Ottobock system hand or an Ottobock hook to the Ottobock modular arm prosthesis.

This adapter allows for rotation and flexion of the terminal device. The restraint of the rotation is adjustable by screwing the screw A (Fig. 1, item A). The restraint of the flexion is adjustable by tightening the screw B (Fig. 2, item B).

The flexion and extension angle amounts to 30°. The rotation and flexion movements are not locked.

The 10R4 Knurled Plate is used for attaching the passive prosthesis hand to the 10R2 or 10R3 Adapters.

### 3 Technical Data

Article number	10R1=*	10R2=*	10R3=*	10R4=*
Weight (g)	25	20	35	15
Total height (mm)	43	30	44	38

## 4 Handling

### 4.1 Assembly

Connect distal tube end to the Modular Adapter (10R2=\* / 10R3=\*).

#### INFORMATION

The correct torque of the cap screw is 3.5 Nm (Fig. 3). Use 710D6 Torque Wrench to tighten the cap screw.

#### **4.1.1 Connecting to a passive prosthesis hand:**

- Connect the Light Metal Tube (13R3/ 13R4) or Modular Arm Component (12R6=\* / 12R7=\*) to the 10R1 Adapter and screw the prosthesis hand onto the adapter.
- Connect the 10R4=\* Knurled Plate to the Adapter (10R2=\* / 10R3=\*).
- Connect the Light Metal Tube (13R3/ 13R4) or Modular Arm Component (12R6=\* / 12R7=\*) to the Adapter (10R2=\* / 10R3=\*) and screw the 10R4=\* Knurled Plate into the prosthesis hand.

#### **4.1.2 Connecting to an Ottobock system hand or Ottobock hook:**

- Connect the Light Metal Tube (13R3/ 13R4) or Modular Arm Component (12R6=\* / 12R7=\*) to the Adapter (10R2=\* / 10R3=\*) and screw into the Ottobock system hand or Ottobock hook.

## **5 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **5.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **5.2 CE conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH