














757B35=1, 757B35=3

DE Gebrauchsanweisung (Fachpersonal)	4
EN Instructions for use (qualified personnel)	13
FR Instructions d'utilisation (Personnel spécialisé)	22
IT Istruzioni per l'uso (Personale tecnico specializzato)	31
ES Instrucciones de uso (Personal técnico especializado)	40
PT Manual de utilização (Pessoal técnico)	50
NL Gebruiksaanwijzing (Vakmensen)	59
SV Bruksanvisning (Fackpersonal)	69
DA Brugsanvisning (Faguddannet personale)	77
NO Bruksanvisning (Fagpersonell)	86
TR Kullanma talimatı (Uzman personel)	95
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	103
RU Руководство по применению (Квалифицированный персонал)	113

INFORMATION

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INFORMATION

Last update: 2015-09-22

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

In the following, the 757B35=* MyoEnergy Integral is referred to simply as the product.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Instruct the patient in the proper use and care of the product. The product may not be transferred to the patient without prior instruction.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Construction and Function

The product serves to provide power to the prosthesis system.

The product consists of several components. These components are inseparable and ready for mounting.

2.1.1 Components

The product consists of the following components:

Charging receptacle

The charging receptacle (see fig. 1, item 1) has the following functions:

- Contact for charging the battery
- LED to indicate the current battery charge level
- LED for providing feedback on operating states
- Button to switch the prosthesis system on and off
- Beeper for providing feedback on operating states

Rechargeable battery

The rechargeable battery (see fig. 1, item 2) consists of 2 cells with integrated electronics. This protects the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range.

Communication cable

The communication cable (see fig. 1, item 3) is used to exchange data and connects the respective prosthesis component to the battery.

Supply cable

The supply cable (see fig. 1, item 4) establishes the connection between the battery and the respective prosthesis component.

3 Application

3.1 Indications for use

The product to be used **exclusively** to supply energy to myoelectric upper limb prosthesis components.

3.2 Area of Application

The product has been developed for use in everyday activities. It must not be used for activities where it would be exposed to unallowable environmental conditions or strains (e.g. extreme sports).

3.3 Conditions of use

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see Page 21).

3.4 Qualification of the prosthetist

The fitting of a patient with the product may only be carried out by an orthopaedic technician.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 CAUTION
<p>The heading describes the source and/or the type of hazard</p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none">> E.g.: Consequence 1 in case of failure to observe the hazard> E.g.: Consequence 2 in case of failure to observe the hazard <p>▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.</p>

4.3 General safety instructions

 CAUTION
<p>Changing prosthetic components when switched on</p> <p>Injury due to faulty control or malfunction of the prosthetic system.</p> <p>▶ Before changing prosthetic components (e.g. gripping components) switch the product off.</p>

⚠ CAUTION

Penetration of dirt and humidity into the product

Injury due to unexpected behaviour of the prosthetic system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ Do not use the product for bathing prostheses.

NOTICE

Improper care of the housing

Damage to the casing through the use of acetone, white spirit or similar solvents.

- ▶ Only clean the housing with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

INFORMATION

Prior to initial operation of the product, charge the supplied battery.

4.4 Patient information

⚠ WARNING

Operating the product near active implanted systems

Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the product.

- ▶ When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.
- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ CAUTION

Penetration of dirt and humidity into the product

Injury due to unexpected product behaviour of the prosthesis system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the prosthetic system.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see Page 21).

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

CAUTION

Independent manipulation of the product

Injury due to malfunction and resulting unexpected prosthesis actions.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

CAUTION

Charging the product with soiled charging receptacle contacts

Injury due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contact surface coating with pointed or sharp objects.

CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850/GSM 900: 0.99 m
 - Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points,...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

5 Scope of Delivery and Accessories

Scope of Delivery

- 1 pc. 757B35=1/757B35=3 MyoEnergy Integral
- 1 pc. 647G533 instructions for use (qualified personnel)

Additional scope of delivery for socket fabrication

- 2 pc. lamination dummy
- 1 pc. lamination dummy – charging receptacle
- 1 pc. drilling template for charging receptacle
- 2 pc. battery installation aid
- 2 pc. Hook and loop closure strap (only included in the scope of delivery for the 757B35=3 MyoEnergy Integral)
- 4 pc. countersunk head screws

Required accessories for operation

- 1 pc. battery charger 757L35 MyoCharge Integral
- 757VAK1 battery adapter cable (only for combination with VASI system)

6 Preparation for use

6.1 Required Tools

- Screwdriver
- Adhesive tape
- 726W9=* conical drill
- Hand drill ø 4 mm

6.2 Recommended materials

- 623T1* nylon tube sock
- 99B71=* PVC sheeting
- 627B5=* Neoprene cell rubber tape
- 616T3=* ThermoLyn Trolene
- 200 617H12=* Pedilen rigid foam
- 617P21=* Pedilen hardener
- 99B81=* PVA bag
- 623T* Perlon stockinette
- 616G12 carbon fibre cloth
- 616G12 carbon fibre cloth

6.3 Laminating instructions

The product can be positioned in the upper arm or forearm socket. Depending on how much space is available, the battery is secured in the outer socket or between the outer and inner socket.

NOTICE

Use of the product with a damaged flat cable

Malfunction due to incorrect or interrupted data transfer.

- ▶ Avoid bending or pressing the flat cable multiple times in the same position.

INFORMATION

The length of the flexible flat cable between the battery and the charging receptacle can be adjusted by carefully bending or folding it. Only light pressure must be applied when pinching and sliding down the folded edge. If the flexible flat cable is bent in the opposite direction at a spot where it has already been folded, serious damage can occur.

6.3.1 Laminating instructions for the upper arm socket

- 1) Position the battery as far proximal as possible on the alignment aid for the ErgoArm elbow component.
- 2) Use the cap screws to screw the battery to the alignment aid from the inside (see fig. 2).
- 3) Connect the 2-pin supply cable to the electronics housing of the elbow joint.

INFORMATION: The 3-pin communication cable is not plugged in.

6.3.2 Laminating instructions for the forearm socket

- 1) Fabricate the inner socket using the familiar technique.
- 2) Cut nylon tube sock to twice the length, pull one half over the inner socket and tie it off.
- 3) Establish the position of the battery dummies (on the inner socket for patients with a long residual limb).
- 4) Pull the rest of the nylon tube sock down so it is free of wrinkles and tie off distally.
- 5) Warm the PVC sheeting and pull it over the residual limb.
- 6) Tie off the PVC sheeting and evacuate.
- 7) Glue on neoprene cell rubber tape in 2 layers according to the course of the outer socket edge.

- 8) Wrap and fasten Trolen around the inner socket.
- 9) Fabricate the foam positive using rigid foam and hardener in the familiar technique.
- 10) Sand the foam positive (until the battery dummies are visible, Schaum beschliffen).
- 11) Position the charging receptacle dummy on the foam positive.
- 12) Pull a PVA bag over the foam positive.
- 13) Pull on several layers of Perlon stockinette and tie off at the lamination ring.
INFORMATION: To reinforce the socket, carbon fibre cloth can be glued between the first and second layers of Perlon stockinette. The carbon fibre cloth is not allowed to touch or overlap the lamination ring.
- 14) Laminate the outer socket using the familiar technique and allow it to harden.
- 15) Remove the outer socket from the inner socket.
- 16) Remove the plaster from the inner socket.
- 17) Sand the contours on the outer socket.
- 18) Align the bore holes of the charging receptacle template with the knobs on the outer socket and transfer the template contours (see fig. 4).
- 19) Burn the knobs free from the inside (see fig. 5).
- 20) After burning them free, use a conical drill to bore the nubs out to the marked contour (see fig. 6).
- 21) Mill the contour and check it with the charging receptacle template (see fig. 7).
- 22) Insert the charging receptacle through the cutout from the inside of the socket (see fig. 8). The charging receptacle must fit flush to the outside of the socket without pressure being applied (see fig. 9).
INFORMATION: If the wall thickness of the outer socket is more than 2 mm, the inside of the socket must be sanded.
- 23) Apply bonding agent to the small rim of the charging receptacle (see fig. 10).
- 24) Insert the charging receptacle into the outer socket and allow the bonding agent to harden.
- 25) For patients with long residual limbs, lay the battery in the laminated recess. For patients with short residual limbs, position the batteries in the available space of the outer socket.

6.4 Fastening the MyoEnergy Integral in the outer socket

Assembly of the separate battery units

- 1) Determine the battery positions using the alignment aid.
INFORMATION: Maintain sufficient distance between the installation aids/batteries (see fig. 11).
- 2) Mark the hole positions on the socket and drill them with the drill bit.
- 3) Pull the double-sided adhesive tape off the installation aids and apply it to the battery (see fig. 13).
- 4) Secure the installation aid and battery with the hook and loop closure strap (757B35=3) or adhesive tape (757B35=1).
INFORMATION: Lay the hook and loop closure strap into the groove on the installation aid (see fig. 14). The overlapping ends of the hook and loop closure strap have to be up against the battery.
- 5) Use the countersunk head screws to screw the batteries to the socket.

Installation of the assembled battery units

- 1) Fasten the batteries and installation aid with the double-sided adhesive tape (see fig. 16).
- 2) Wrap the installation aid and batteries with the hook and loop closure strap (757B35=3) or adhesive tape (757B35=1).
INFORMATION: Lay the hook and loop closure strap into the groove on the installation aid.

7 Handling

7.1 Charging the battery

- 1) Connect the charging plug to the charging receptacle of the product.
 - The correct connection of the battery charger to the product is indicated by feedback (Status signals).
 - The prosthesis system is switched off automatically.
 - The charging process starts.
- 2) Disconnect the product after the charging process is complete.

7.2 Display of the current charge level

The battery charge level can be queried at any time.

- 1) With the prosthetic system switched on, press the charging receptacle button (Taste betätigen) and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

7.3 Safety shutoff

The purpose of the safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature during the charging process
- Short circuit
- Overvoltage or undervoltage

INFORMATION

After a safety shutdown, the battery charger must be connected to reactivate the electronics.

8 Cleaning and Care

NOTICE

Improper care of the housing

Damage to the casing through the use of acetone, white spirit or similar solvents.

- ▶ Only clean the housing with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

9 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

10.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11 Appendices

11.1 Symbols Used

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

LOT:PPPP YYYY WW

Lot number

Legal manufacturer

11.2 Operating states/error signals

11.2.1 Status signals

Charging receptacle	Status
	Battery fully charged
	Battery 50% charged
	Battery empty

11.2.2 Beep signals

The following beep signals are set by default:

Beep signal	Additional display	Status
1 x long	-	<ul style="list-style-type: none">• Turn off prosthesis on charging receptacle• Charging begins (charging plug connected to charging receptacle)• Charging ends (charging plug disconnected from charging receptacle)

Beep signal	Additional display	Status
2 x short	LED on charging receptacle lights briefly	Switch on prosthesis
3 x short	-	Battery voltage too low, prosthesis system shuts off automatically

The following beep signals exist only in conjunction with the 7in1 Controller (9E420=*):

Beep signal	Status
1 x long	Disconnect the MyoLino Link: Prosthesis was previously switched off
2 x short	Connect the MyoLino Link
2 x short	Disconnect the MyoLino Link: Prosthesis was previously switched on
Continuous pulses (as long as the button is being pressed)	Parental access switch

11.3 Technical data

Ambient conditions	
Storage and transport in original packaging	-20 °C/-4 °F to +50 °C/+122 °F
Storage and transport without packaging	-20 °C/-4 °F to +50 °C/+122 °F Max. 93% relative humidity, non-condensing
Operation	0 °C/+32 °F to +60 °C/+140 °F Max. 93% relative humidity, non-condensing
Charging the battery	0 °C/+32 °F to +45 °C/+113 °F

Battery of the product	757B35=1	757B35=3
Battery type	Li-Po	Li-Ion
Output voltage	approx. 7.4 V	approx. 7.4 V
Charging voltage	approx. 8.2 V	approx. 8.2 V
Capacity	600 mAh	1150 mAh
Dimensions	52 x 25 x 9 mm	52 x 36 x 9 mm
Weight	45 g	70 g
Operating time of prosthesis	approx. 1500 – 2000 grip cycles with MyoBock Adult System approx. 4000 grip cycles with MyoBock Hand System for children	approx. 4000 grip cycles with MyoBock Adult System
Service life of the MyoEnergy Integral	2 years	2 years
Battery service life	2 years	2 years
Charging time until battery is fully charged	2.5 h	2.5 h

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Otto Bock Healthcare Products GmbH