
ottobock.

2WR95, 2WR95=1, 4WR95=3

| | | |
|-----------|----------------------------------|-----|
| DE | Gebrauchsanweisung | 4 |
| EN | Instructions for Use | 12 |
| FR | Instructions d'utilisation | 20 |
| IT | Istruzioni per l'uso | 29 |
| ES | Instrucciones de uso | 37 |
| PT | Manual de utilização | 45 |
| NL | Gebruiksaanwijzing | 53 |
| SV | Bruksanvisning | 62 |
| DA | Brugsanvisning | 70 |
| NO | Bruksanvisning | 77 |
| PL | Instrukcja użytkowania | 85 |
| HU | Használati utasítás | 93 |
| CS | Návod k použití | 101 |
| RO | Instrucțiuni de utilizare | 109 |
| HR | Upute za uporabu | 118 |
| TR | Kullanma talimatı | 126 |
| EL | Οδηγίες χρήσης | 134 |
| RU | Руководство по применению | 142 |
| JA | 取扱説明書 | 151 |
| ZH | 使用说明书 | 159 |

| | ● 4D27 | ▲ | ■ |
|-----------|---------------------------------|----------------------------------|------------------------|
| DE | Einzelteile-Pack | Mindestmenge | Einzelteile |
| EN | Single-Component Pack | Minimum order quantity | Single components |
| FR | Kit de pièces de rechange | Quantité minimum | Pièces à l'unité |
| IT | Imballo dei singoli componenti | Quantità minima | Singoli componenti |
| ES | Kit componentes | Cantidad mínima | Componentes |
| PT | Conjunto de peças de reposto | Quantidade mínima | Componentes |
| NL | Onderdelenpakket | Minimikvantitet | Onderdelen |
| SV | Service-Set | Minimum aantal | Separata delar |
| DA | Komponentpakke | Mindste mængde | Komponenter |
| PL | Zestaw naprawczy | Ilość minimalna | Podzespoły |
| HU | Alkatrész-csomag | Minimum mennyiség | Alkatrészek |
| CS | Sada jednotlivých dílů | Minimální množství pro objednání | Jednotlivé díly |
| RO | Set componente individuale | Cantitate minimă | Componente individuale |
| HR | Pakiranje pojedinačnih dijelova | Minimalna količina | Pojedinačni dijelovi |
| TR | Münferit parça ambalajı | Minimum miktar | Münferit parçalar |
| EL | Σετ μεμονωμένων εξαρτημάτων | Ελάχιστη ποσότητα | Μεμονωμένα εξαρτήματα |
| RU | Комплект деталей | Минимальное количество | Отдельные детали |
| JA | セットで発注いただく部品 | 複数入バックで発注いただく部品 | 1個から発注いただける部品 |
| ZH | 维修组件 | 最低起订量 | 单个部件 |

INFORMATION

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INFORMATION

Date of the last update: 2012-08-13

- Please read this document carefully.
- Follow the safety instructions.

1.1 Function

The 2WR95=1 and 2WR95 tube adapters are used in combination with the 4WR95=3 tube clamp adapter as adjustable, variable-length connection elements for modular, waterproof Ottobock walking devices.

The tube clamp adapter must always be installed proximally. The clamping slot must always be positioned anteriorly.

1.2 Combination possibilities

The adapters 2WR95, 2WR95=1 and 4WR95=3 may only be used in combination with Ottobock components for modular waterproof walking devices.

2 Application

2.1 Intended Use

The adapters are intended **exclusively** for exoprosthetic fittings of the lower extremities.

2.2 Area of Application

The adapters 2WR95, 2WR95=1 and 4WR95=3 are intended **exclusively** for use in a waterproof walking device used as a secondary orthosis.

They may only be used in combination with Ottobock components for waterproof walking devices.

The adapters are water resistant.

2.3 Conditions of Use

The adapters 2WR95, 2WR95=1 and 4WR95=3 are approved for a body weight up to 150 kg (330 lbs).

3 Safety

3.1 Explanation of the ANSI symbols

| | |
|---|--|
|  CAUTION | Warnings regarding possible risks of accident or injury. |
|  NOTICE | Warnings regarding possible technical damage. |
|  INFORMATION | Additional information on the fitting / use. |

3.2 General Safety Instructions

| |
|---|
|  CAUTION |
| Non-observance of safety instructions |
| Non-observance of the following safety instructions can lead to malfunction of the 2WR95, 2WR95=1, and 4WR95=3 adapters and result in a risk of injury to the patient and damage to the products. |

3.3 Fitting safety instructions

CAUTION

Non-observance of the torque values

Loss of functionality due to structural failure may result if the tightening torques required for the screw connections are not observed. This can cause the patient to fall.

Observing the specified tightening torques is essential for safe operation.

CAUTION

Risk of falls due to improper use

For reasons of structural integrity, the sections cut off from tube adapters may not be used as components in final prostheses.

3.4 Safety instructions for patients / users

Please make your patients aware of the following safety instructions:

CAUTION

Risk of falls due to improper use

The waterproof walking device is not intended to replace the initial prosthesis and should only be used in wet areas. It may not be used for extreme activities such as extreme sports (water skiing, kite surfing, etc.).

CAUTION

Risk of falls due to improper use

The walking speed must be adjusted accordingly in wet areas and in areas where a higher risk of falls can be expected.

CAUTION

Risk of burns in case of heat exposure

The Ottobock components for waterproof walking devices are not suitable for use in saunas or for activities during which the prosthesis is constantly exposed to high heat above 50 °C (122 °F).

NOTICE

Risk of loss of function due to improper storage/transport

The waterproof walking device may only be stored or transported at temperatures between -10°C (14°F) and +50°C (122°F).

NOTICE

Risk of loss of function due to improper care

The prosthesis must be cleaned with clean fresh water after any contact with saltwater or water containing chlorine/soap or contact with abrasive substances (e.g. sand). To clean the prosthesis, it should be placed in a sink (tub) with clean, fresh water and moved around to rinse the corners, edges, and the tube.

The prosthesis must also be dried immediately after every use.

Contact with abrasive substance (e.g. sand) can lead to premature wear.

4 Scope of Delivery

The scope of delivery is shown on the front page.

5 Preparation for use

INFORMATION

According to the Ottobock alignment recommendations, bench alignment of the waterproof walking device with the 1WR95 Aqua Foot on 0 mm heel height is usually only possible with the aid of the angled tube adapter. Use the 2WR95=1 tube adapter which is angled at the top. If necessary, the 2WR95 straight tube adapter can be used in the thigh area for TF patients.

5.1 Alignment Recommendations

INFORMATION

The following instructions apply to the tube adapters 2WR95 and 2WR95=1 and the tube clamp adapter 4WR95=3.

5.2 Shortening the tube adapter

⚠ CAUTION

Risk of material damage due to deformation

The tube may not be clamped in a vise while working on it.

INFORMATION

The cut must be perpendicular to the longitudinal axis of the tube (see Fig. 2).

The tube adapter must be shortened according to the patient's measurements. Special care must be taken to ensure that the tube is not deformed. A deformed tube cannot be fit easily into the tube clamp adapter.

We recommend using the 719R2 Tube Cutter or a similar tool.

- 1) Use the tube cutter to shorten the tube adapter to fit the patient's measurements.
- 2) Deburr the inside and outside of the cut edge using the 718R1 tube deburring tool (see Fig. 3) or a similar tool.
- 3) Loosen the tube clamp adapter cylinder head screw by approx. 2 turns.
- 4) Slide the tube 5 mm into the tube clamp adapter while turning it slightly.
- 5) Slide the tube adapter into the tube clamp adapter until the tube rests against the bottom of the collar.

5.3 Adjustment

NOTICE

Risk of corrosion

The cylinder head screw and the set screws are protected against corrosion and may only be replaced with set screws from the 4D27 single-component pack.

The four set screws can be used to make static adjustments at any time during alignment, trial fittings, and even after the prosthesis is finished. The previously adjusted position can be maintained if a module is replaced or during disassembly, if only the two set screws which are located beside one another and turned in the furthest are removed.

Use the 710D4 torque wrench to tighten the set screws.

- For trial fittings, tighten the adapter set screws clockwise with a torque of 10 Nm.

5.4 Final assembly

CAUTION

Risk of falls due to changes in the prosthesis geometry

The adapter may turn while the prosthesis is in use if there is any contamination on the contact surfaces.

Prior to final assembly (i.e. prior to inserting the tube adapter into the tube clamp adapter), the inside of the tube clamp adapter and the outside part of the tube adapter which slides into the tube clamp adapter must be cleaned with 634A3 Acetone!

CAUTION

Risk of falling due to loosening of the set screws while the prosthesis is in use

The applied strains can cause the set screws to come loose.

- Secure the set screws with 636K13 Loctite® and observe the prescribed tightening torques.

- 1) Use 636K13 Loctite® to install the set screws.
- 2) Pre-tighten the set screws in the clockwise direction using the 710D4 torque wrench.
Tightening torque: **10 Nm**.
- 3) During final assembly, tighten the set screws in the clockwise direction using the 710D4 torque wrench.
Tightening torque **15 Nm** (see Fig. 1 ●)

Set screws which seem too long or too short for the adjusted position should be replaced by suitable corresponding ones from the 4D27 single-component pack.

Use the 710D4 torque wrench to tighten the cylinder head screw!

Use the 710D4 torque wrench to tighten the cylinder head screw, both for trial fitting and for final assembly.

Tightening torque **10 Nm** (see Fig. 1 ●)

- 4) The end of the M5X30 cylinder head screw must be covered with the included shrink tubing. To do so, cut a suitable piece off at a right angle, place it over the end of the screw, and use a hot air gun or a small flame to shrink it until it is sufficiently tight.

6 Maintenance/Repair

INFORMATION

The fundamental design of all Ottobock adapters allows for a possible period of use of 3 to 5 years, depending on the patient's activity level.

In general, we recommend regular annual safety inspections.

The safety inspections must be performed by a prosthetist.

6.1 Spare parts

Single components

none

Minimum order quantity

none

Single-component pack

4D27 contains:

4 M8x10 corrosion protected set screws

8 M8x12 corrosion protected set screws

8 M8x14 corrosion protected set screws

4 M8x16 corrosion protected set screws

2 M5x22 corrosion protected cylinder head screws

2 rounded washers

1 M5x30 corrosion protected cylinder head screw

1 plastic ring

3 cm shrink tubing

7 Disposal

7.1 Disposal Information

INFORMATION

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

8 Legal information

8.1 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

8.2 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH