
2WR95, 2WR95=1, 4WR95=3

DE	Gebrauchsanweisung	4
EN	Instructions for use	9
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	23
PT	Manual de utilização	28
NL	Gebruiksaanwijzing	33
SV	Bruksanvisning	38
DA	Brugsanvisning	42
NO	Bruksanvisning	47
PL	Instrukcja użytkowania	51
HU	Használati utasítás	56
CS	Návod k použití	61
RO	Instrucțiuni de utilizare	66
HR	Upute za uporabu	71
TR	Kullanma talimatı	75
EL	Οδηγίες χρήσης	80
RU	Руководство по применению	85
JA	取扱説明書	90
ZH	使用说明书	94

	●	▲	■
DE	Einzelteile-Pack	Mindestmenge	Einzelteile
EN	Single-Component Pack	Minimum order quantity	Single components
FR	Kit de pièces de rechange	Quantité minimum	Pièces à l'unité
IT	Imballo dei singoli componenti	Quantità minima	Singoli componenti
ES	Kit componentes	Cantidad mínima	Componentes
PT	Conjunto de peças de reposto	Quantidade mínima	Componentes
NL	Onderdelenpakket	Minimikvantitet	Onderdelen
SV	Service-Set	Minimum aantal	Separata delar
DA	Komponentpakke	Mindste mængde	Komponenter
PL	Zestaw naprawczy	Ilość minimalna	Podzespoły
HU	Alkatrész-csomag	Minimum mennyiség	Alkatrészek
CS	Sada jednotlivých dílů	Minimální množství pro objednání	Jednotlivé díly
RO	Set componente individuale	Cantitate minimă	Componente individuale
HR	Pakiranje pojedinačnih dijelova	Minimalna količina	Pojedinačni dijelovi
TR	Münferit parça ambalajı	Minimum miktar	Münferit parçalar
EL	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα	Μεμονωμένα εξαρτήματα
RU	Комплект деталей	Минимальное количество	Отдельные детали
JA	セットで発注いただく部品	複数入パックで発注いただく部品	1個から発注いただける部品
ZH	维修组件	最低起订量	单个部件

INFORMATION

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INFORMATION

Date of the last update: 2014-09-05

- Please read this document carefully.
- Follow the safety instructions.

INFORMATION

Please read these instructions for use carefully.

Pay special attention to the safety instructions!

1.1 Function

The 2WR95=1 and 2WR95 tube adapters are used in combination with the 4WR95=3 tube clamp adapter as adjustable, variable-length connection elements for modular, waterproof Ottobock walking devices.

The tube clamp adapter must always be installed proximally. The clamping slot must always be positioned anteriorly.

1.2 Combination possibilities

The adapters 2WR95, 2WR95=1 and 4WR95=3 may only be used in combination with Ottobock components for modular waterproof walking devices.

2 Application**2.1 Intended use**

The adapters are intended **exclusively** for exoprosthetic fittings of the lower extremities.

2.2 Area of application

The adapters 2WR95, 2WR95=1 and 4WR95=3 are intended **exclusively** for use in a waterproof walking device used as a secondary prosthesis.

They may only be used in combination with Ottobock components for waterproof walking devices. The adapters are water resistant.

2.3 Conditions of use

The adapters 2WR95, 2WR95=1 and 4WR95=3 are approved for a body weight up to 150 kg (330 lbs).

3 Safety**3.1 Explanation of the ANSI symbols**

 **CAUTION** Warnings regarding possible risks of accident or injury.

 **NOTICE** Warnings regarding possible technical damage.

 **INFORMATION** Additional information on the fitting / use.

3.2 General safety instructions** CAUTION****Non-observance of safety instructions**

Non-observance of the following safety instructions can lead to malfunction of the 2WR95, 2WR95=1, and 4WR95=3 adapters and result in a risk of injury to the patient and damage to the products.

3.3 Fitting safety instructions

CAUTION

Non-observance of the torque values

Loss of functionality due to structural failure may result if the tightening torques required for the screw connections are not observed. This can cause the patient to fall.

Observing the specified tightening torques is essential for safe operation.

CAUTION

Risk of falls due to improper use

For reasons of structural integrity, the sections cut off from tube adapters may not be used as components in final prostheses.

3.4 Safety instructions for patients / users

Please make your patients aware of the following safety instructions:

CAUTION

Risk of falls due to improper use

The waterproof walking device is not intended to replace the initial prosthesis and should only be used in wet areas. It may not be used for extreme activities such as extreme sports (water skiing, kite surfing, etc.).

CAUTION

Risk of falls due to improper use

The walking speed must be adjusted accordingly in wet areas and in areas where a higher risk of falls can be expected.

CAUTION

Risk of burns in case of heat exposure

The Ottobock components for waterproof walking devices are not suitable for use in saunas or for activities during which the prosthesis is constantly exposed to high heat above 60 °C (122 °F).

NOTICE

Risk of loss of function due to improper storage/transport

The waterproof walking device may only be stored or transported at temperatures between -10 °C (14°F) and +60 °C (122°F).

NOTICE

Risk of loss of function due to improper care

The prosthesis must be cleaned with clean fresh water after any contact with saltwater or water containing chlorine/soap or contact with abrasive substances (e.g. sand). To clean the prosthesis, it should be placed in a sink (tub) with clean, fresh water and moved around to rinse the corners, edges, and the tube.

The prosthesis must also be dried immediately after every use.

Contact with abrasive substance (e.g. sand) can lead to premature wear.

4 Scope of delivery

The scope of delivery is shown on the front page.

5 Preparation for use

INFORMATION

According to the Ottobock alignment recommendations, bench alignment of the waterproof walking device with the 1WR95 Aqua Foot on 0 mm heel height is usually only possible with the aid of the angled tube adapter. Use the 2WR95=1 tube adapter which is angled at the top. If necessary, the 2WR95 straight tube adapter can be used in the thigh area for TF patients.

5.1 Alignment recommendations

INFORMATION

The following instructions apply to the tube adapters 2WR95 and 2WR95=1 and the tube clamp adapter 4WR95=3.

5.2 Shortening the tube adapter

CAUTION

Risk of material damage due to deformation

The tube may not be clamped in a vise while working on it.

INFORMATION

The cut must be perpendicular to the longitudinal axis of the tube (see Fig. 2).

The tube adapter must be shortened according to the patient's measurements. Special care must be taken to ensure that the tube is not deformed. A deformed tube cannot be fit easily into the tube clamp adapter.

We recommend using the 719R2 tube cutter or a similar tool.

- 1) Use the tube cutter to shorten the tube adapter to fit the patient's measurements.
- 2) Deburr the inside and outside of the cut edge using the 718R1 tube deburring tool (see Fig. 3) or a similar tool.
- 3) Loosen the tube clamp adapter cap screw by approx. 2 turns.
- 4) Slide the tube 5 mm into the tube clamp adapter while turning it slightly.
- 5) Slide the tube adapter into the tube clamp adapter until the tube rests against the bottom of the collar.

5.3 Adjustment

NOTICE

Risk of corrosion

The set screws of stainless steel can optionally be replaced by 506G5* titanium set screws, which increase the protection against corrosion.

The four set screws can be used to make static adjustments at any time during alignment, trial fittings, and even after the prosthesis is finished. The previously adjusted position can be maintained if a module is replaced or during disassembly, if only the two set screws which are located beside one another and turned in the furthest are removed.

Use the 710D4 torque wrench to tighten the set screws.

- For trial fittings, tighten the adapter set screws clockwise with a torque of 10 Nm.

5.4 Final assembly

CAUTION

Risk of falls due to changes in the prosthesis geometry

The adapter may turn while the prosthesis is in use if there is any contamination on the contact surfaces.

Prior to final assembly (i.e. prior to inserting the tube adapter into the tube clamp adapter), the inside of the tube clamp adapter and the outside part of the tube adapter which slides into the tube clamp adapter must be cleaned with a degreasing cleaning agent (e.g. 634A3 acetone)!

CAUTION

Risk of falling due to loosening of the set screws while the prosthesis is in use

The applied strains can cause the set screws to come loose.

- Secure the set screws with 636K13 Loctite® and observe the prescribed tightening torques.

- 1) Use 636K13 Loctite® to install the set screws.
- 2) Pre-tighten the set screws in the clockwise direction using the 710D4 torque wrench.
Tightening torque: **10 Nm**.
- 3) During final assembly, tighten the set screws in the clockwise direction using the 710D4 torque wrench.
Tightening torque **15 Nm** (see Fig. 1 ●)

Replace set screws which seem too long or too short for the adjusted position by set screws of suitable length.

Use the 710D4 torque wrench to tighten the cap screw!

Use the 710D4 torque wrench to tighten the cap screw, both for trial fitting and for final assembly.

Tightening torque **10 Nm** (see Fig. 1 ●)

6 Maintenance/repair

INFORMATION

The fundamental design of all Ottobock adapters allows for a possible period of use of 3 to 5 years, depending on the patient's activity level.

In general, we recommend regular annual safety inspections.

The safety inspections must be performed by a prosthetist.

6.1 Spare parts

Single components

- | | | |
|-------|--|-----|
| 1 pc. | 501T24=M5x25 (titanium cap screw) | (2) |
| 1 pc. | 4X28=3 (plastic ring) | (3) |
| 1 pc. | 4Y423 (clamping bolt) | (4) |
| 1 pc. | 4Y424 (clamping bolt with interior thread) | (5) |

Minimum order quantity

- | | | |
|-------|---|--------------------------|
| 4 pcs | 506G3=M8X12-"Niro" (stainless steel set screws) | |
| 4 pcs | 506G3=M8X14-"Niro" (stainless steel set screws) | (1) included in delivery |
| 4 pcs | 506G3=M8X16-"Niro" (stainless steel set screws) | |

Single-component pack

none

7 Disposal

7.1 Disposal Information

INFORMATION

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

8 Legal information

8.1 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

8.2 CE conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH