

















7E9

DE Gebrauchsanweisung	5
EN Instructions for use	21
FR Instructions d'utilisation	37
IT Istruzioni per l'uso	55
ES Instrucciones de uso	73
PT Manual de utilização	91
NL Gebruiksaanwijzing	107
SV Bruksanvisning	123
DA Brugsanvisning	139
NO Bruksanvisning	155
HU Használati utasítás	171
CS Návod k použití	187
RU Руководство по применению	203
JA 取扱説明書	221

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1	Product description	22
1.1	Construction and Function	22
1.2	Combination Possibilities	22
2	Application	23
2.1	Intended Use	23
2.2	Area of Application	23
2.3	Environmental conditions.....	23
2.4	Service Life.....	23
3	Safety	23
3.1	Explanation of warning symbols.....	23
3.2	General Safety Instructions.....	24
4	Scope of Delivery	25
5	Preparation for Use.....	26
5.1	Taking the Plaster Negative	27
5.2	Fabricating the Test Socket	27
5.3	Determining the Alignment Reference	27
5.4	Bench Alignment	28
5.5	Static Alignment	29
5.6	Dynamic Trial Fitting	30
5.6.1	Adjustments.....	30
5.6.2	Adjustment Possibilities.....	30
5.6.3	Gait Parameter Adjustment	31
5.6.4	Correcting Possible Problems.....	31
5.7	Finishing the Trial Prosthesis	33
5.8	Reinforcing the Prosthetic Socket.....	33
5.9	Finishing the Prosthesis	34
6	Gait Training	34
7	Cleaning	35
8	Maintenance	35
9	Disposal.....	35
10	Legal information	35
10.1	Liability	35
10.2	CE conformity	35
11	Technical Data	35

1 Product description

INFORMATION

Date of last update: 2022-03-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The prosthetic hip joint is a monocentric hip joint with hydraulic stance and swing phase control. Its design features enable the prosthesis wearer to achieve a harmonious and smooth gait pattern. The prosthetic hip joint offers the following features:

- **Extension and flexion with hydraulic damping:** The hydraulics control joint damping in all gait phases.
This permits damping of hip joint extension after stepping down, with a reduction of backward pelvic tilt and therefore hyperlordosis.
The step length (flexion of the hip joint in the swing phase) is divided into 2 ranges. In the start range there is no swing phase damping. In the end range, swing phase damping is adjustable.
- **Sitting characteristics:** Since the hip joint has a low structural height when installed, the "pelvic obliquity" in the sitting position can be reduced to a minimum. The very large flexion angle of the hip joint makes comfortable sitting possible.
- **Adjustment options for alignment:** Interior and exterior rotation, abduction and adduction, pyramid adapter on the distal end of the hip joint.

1.2 Combination Possibilities

Designation	Reference number
Tube clamp adapters (for thigh area)	Up to 125 kg: 4R82, 4R156, 4R156=1, 4R156=2 Up to 100 kg: 4R52, 4R56, 4R56=1, 4R56=2
Tubes (for thigh area)	Up to 125 kg: 2R36 Up to 100 kg: 2R30
Rotation adapter	4R57
Prosthetic knee joints	Up to 125 kg: 3R60=HD, Kenevo 3C60, C-Leg 3C98-3, Genium 3B1-3, Genium X3 3B5-3 Up to 100 kg: 3R106=HD
Torsion adapters	3R60=HD + 3R106=HD: 4R39 3C98-3 C-Leg: 2R67 3B1 Genium: 2R21
Tube adapters (for lower leg area)	Up to 125 kg: 2R57, 2R58 Up to 100 kg: 2R37, 2R38
Prosthetic feet	Up to 75 kg: 1D35 22 – 25 cm Up to 80 kg: 1C30 21 – 22 cm Up to 95 kg: 1C30 23 – 24 cm Up to 100 kg: 1D35 26 – 30 cm, 1A30 26 – 29 cm, 1C60 21 – 24 cm, 1C61 21 – 24 cm Up to 110 kg: 1C30 25 – 26 cm Up to 125 kg: 1S90 26 – 28 cm, 1E50, 1E51, 1E56,

Designation	Reference number
	1E58, 1C30 27 – 30 cm, 1C60 25 – 30 cm, 1C61 25 – 30 cm, 1D10, 1D11
Cosmetics	3S27

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.


2 Application

2.1 Intended Use

The prosthetic hip joint system is intended exclusively for the exoprosthetic fitting of patients with hip amputations such as intertrochanteric amputation, hip disarticulation or hemipelvectomy.

2.2 Area of Application

The prosthetic hip joint can be used by active prosthesis wearers as well as prosthesis wearers who require a high degree of safety. It offers dynamic comfort through hydraulic stance phase and swing phase control.

Recommended for mobility grades **2 and 3** (restricted outdoor walker and unrestricted outdoor walker). Approved for a body weight of **125 kg max.**

The prosthetic hip joint was designed for activities of daily life; it must not be used for sports.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range:

Use: -10 °C – +45 °C

Storage, Transportation: -20 °C – +60 °C

Relative humidity: 20 % – 90 % (Non-condensing)

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

2.4 Service Life

This component was tested for 2 million load cycles in accordance with ISO 15032. Depending on the user's level of activity, this corresponds to a service life of 2 to 3 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the MOBIS classification (see section "Area of Application").
- ▶ **Inform the patient.**

CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

CAUTION

Exceeding the service life and reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service Life").
- ▶ Only use the product for a single patient.
- ▶ **Inform the patient.**

CAUTION

Manipulation or modification of prosthesis components

Falling due to malfunction or loss of functionality

- ▶ Do not manipulate or modify the prosthesis components.
- ▶ Only complete the tasks described in these instructions for use.
- ▶ Have repairs completed by the manufacturer's customer service only.

► **Inform the patient.**

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- Do not reach into the joint mechanism during daily use.
- Close attention is required during assembly and adjustment tasks.
- **Inform the patient.**

⚠ CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthesis components

- Do not touch overheated prosthesis components.
- In case of changes in functionality, stop all activities and allow the overheated prosthesis components to cool down.
- In case of changes in functionality, reduce strain on the prosthesis component in order to allow cooling.
- In case of overheating or changes in functionality, have the prosthesis component inspected for damage by a service centre authorised by the manufacturer.
- **Inform the patient.**

NOTICE

Mechanical damage to the product

Damage due to change in or loss of functionality

- Use caution when working with the product.
- If the product is damaged check it for proper function and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- **Inform the patient.**

Signs of changes in or loss of functionality during use

Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

4 Scope of Delivery

The scope of delivery for the product is shown on page 2 (Figure 1).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

7E9 Hip Joint					
Fig.	No.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647G774=all_INT
1	1	–	1	Hip Joint	–
1	2	■	1	Support Plate	4G576
1	3	■	1	Lamination Plate (up to 125 kg)	7Z53=1-M10
1	4	▲	2	Oval Head Screw (up to 125 kg)	501S72=M10x25

7E9 Hip Joint					
Fig.	No.		Quantity	Designation	Reference number
1	5	▲	2	Flat Head Screw (up to 125 kg)	501S72=M10x30
1	4	▲	2	Flat Head Screw (up to 100 kg)	501F9=M8x25
1	5	▲	2	Flat Head Screw (up to 100 kg)	501F9=M8x30
1	6	■	1	Hexagon Screwdriver	710H10=2x3
1	7	■	1	Lamination Aid	7Z63
-	-	■	1	T40 Bit	709Z11
2	1	■	1	Lamination Plate (up to 100 kg)	7Z53
2	2	■	1	Bottom Plate (up to 100 kg)	4G420=S
2	3	■	1	Intermediate Piece (up to 100 kg)	4G438=G
-	-	■	1	Bottom Plate (up to 125 kg)	4G420-1=GS
-	-	■	1	Intermediate Piece (up to 125 kg)	4G438-1=G
3	-	■	0	Reference Determination Tool	743A29
	-	■	1	7E9 Quick Start Guide	646D1532=all_INT
1	8	■	1	Pyramid Dome Protector Kit	2Z11=KIT

5 Preparation for Use

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthesis components

- ▶ Observe the alignment, assembly and adjustment instructions.

INFORMATION

Correct alignment is crucial for the optimal function of the prosthesis!
Therefore, a test socket is first used for all of the described tasks.

INFORMATION

Optimum prosthesis alignment requires the use of a functional prosthesis component with torsion unit; this enables a more harmonious and comfortable gait pattern for the patient. It also reduces wear on the prosthetic joint. Either a torsion adapter (e.g. 4R39) or a prosthetic foot with torsion unit (e.g. 1C61) is recommended for this purpose (see the chapter "Combination Possibilities" – see page 22).

Fabrication

The fabrication of the prosthesis is described in the chapters that follow.

- 1) Taking the Plaster Negative (see page 27)
- 2) Fabricating the Test Socket (see page 27)
 - **The tasks that follow are completed with the test socket.**
- 3) Determining the Alignment Reference (see page 27)
- 4) Bench Alignment (see page 28)
- 5) Static Alignment (see page 29)
- 6) Dynamic Trial Fitting (see page 30)
- 7) Finishing the Check Prosthesis (see page 33)
 - **The tasks that follow are completed with the definitive socket.**
- 8) Reinforcing the Prosthetic Socket (see page 33)
- 9) Determining the Alignment Reference (see page 27)
- 10) Bench Alignment (see page 28)
- 11) Static Alignment (see page 29)

12) Dynamic Trial Fitting (see page 30)

13) Finishing the Prosthesis (see page 34)

The 646D1532=all_INT Quick Start Guide (included with the product!) contains useful information on the fabrication process in the form of illustrations.

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

5.1 Taking the Plaster Negative

> Recommended equipment:

- 743A11 Casting Apparatus, 743G5 Hip-Cast Casting Apparatus (includes 743Y26 Pressure Block)
- ▶ Take the plaster negative using the Hip-Cast plaster casting technique. Ensure that the anterior pressure block is applied with 5° exterior rotation.

5.2 Fabricating the Test Socket

- 1) With the help of the plaster negative, cast and model a plaster positive.
- 2) Vacuum form the test socket and fit it to the patient.

5.3 Determining the Alignment Reference

Determining an alignment reference in a neutral position of the hip prosthesis socket determined under load is required for the sagittal placement of the hip prosthesis socket in the prosthesis system. In this neutral position (see fig. 4), the patient perceives virtually no anterior or posterior tilting moments acting on the hip prosthesis socket in the sagittal plane, and the pelvic tilt position is as physiological as possible.

> Recommended tools:

- 743A29 Reference Determination Tool (see fig. 3)
- 1) Adjust the height of the casting table to the ischial tuberosity to ground distance of the patient.
- 2) Place the patient, with optimally applied hip prosthesis socket and without prosthetic components, on the casting table with the prosthesis side, observing the following criteria:
 - Straight stance, balanced at pelvis height
 - Foot and centre of the body are in one line in the a - p direction
- 3) Lower the plate of the casting table by the height of the reference determination tool (approx. 4.8 cm).
- 4) Position the centre of the reference determination tool between the bottom of the hip prosthesis socket and the plate of the casting table, subject to the following criteria:

The reference determination tool must be positioned approximately 5 – 6 cm lateral from the medial socket edge. Align the base plate of the reference determination tool parallel to the edge of the casting table.
- 5) **Determining the neutral position (TMS: partial mass centre of gravity):**

First, place the reference determination tool in the maximal anterior position. The patient will perceive a backward tilt of the hip prosthesis socket (see fig. 5).
Then place the reference determination tool in the maximum posterior position. The patient will perceive a forward tilt of the hip prosthesis socket (see fig. 6).
In small increments, reposition the reference determination tool in the sagittal plane between the maximum positions until the user perceives neither a forward nor a backward tilt (neutral position) (see fig. 4).

6) **Marking the neutral position that was determined:**

Sagittal on the hip prosthesis socket, project and mark a vertical line with a crosshair laser through the centre of the reference determination tool (see fig. 7). For later rotation control of the hip prosthesis socket during bench alignment, mark the line on the distal side of the hip prosthesis socket as far medial as possible and later extend it to the medial socket edge (pelvis rotation line).

Transfer the height of the trochanter major onto this TMS line from the contralateral side and mark it.

Mark the alignment reference line parallel to the TMS line at a distance of 20 mm to posterior.

7) **Mark the reference line for positioning the lamination plate in the frontal plane:**

Transfer half of the measurement y between the centre of the body and the contralateral outer side to the hip prosthesis socket side (y/2) and make a vertical mark (see fig. 8, Pos. 1).

Draw a plumb line from this mark at a distance of 10 – 20 mm (depending on patient circumference) (see fig. 8, Pos. 2).

5.4 Bench Alignment

Alignment Recommendations

(Schematic of bench alignment for the 3R60=HD see fig. 18 - The 646D1532=all_INT Quick Start Guide contains all schematics for bench alignment with the 3R60, 3R106=HD and 3C98-3 C-Leg.)

If the prosthesis is being fabricated with the 3B1 Genium prosthetic knee joint, the instructions for bench alignment in combination with the 7E10* prosthetic hip joint from the seminar documentation and the instructions for use of the prosthetic knee joint as well as X-Soft apply. Notwithstanding this, the distance from the partial gravity line to the alignment reference line from the alignment recommendations that follow must be observed.

⚠ CAUTION

Modifying and exceeding the service life of the lamination plate

Risk of falling due to breakage or deformation of the lamination plate

- ▶ Do not modify the lamination plate (e.g. by shortening or bending).
- ▶ Do not exceed the service life of the lamination plate (see the chapter "Service Life").

Both prosthetic joints must be fully extended during bench alignment.

Bench Alignment Process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
No.	Reference measurements and positions for bench alignment, assembly tasks
	Required tools and equipment: Alignment apparatus (e.g. L.A.S.A.R. Assembly)
①	a-p positioning, middle of the prosthetic foot to the alignment reference line: +30 mm
②	Heel height: effective heel height (x) + 5 mm
	External foot rotation: approx. 5°
③	Structural height of the prosthetic knee joint at the alignment reference point (anterior, upper axis of 3R60=HD + 3R106=HD; rotation axis of 3C98-3): MTP-floor measurement (z) + 20 mm
	a-p positioning of the alignment reference point to the alignment reference line: 3R60=HD + 3R106=HD (anterior, upper axis): 0 mm 3C98-3 (rotation axis): 0 to +5 mm

Bench Alignment Process	
+ = anterior placement / - = posterior placement (in relation to the alignment reference line)	
No.	Reference measurements and positions for bench alignment, assembly tasks
3	External rotation of the prosthetic knee joint: approx. 5°
4	Join the prosthetic foot and the prosthetic knee joint with the chosen adapters (e.g. torsion adapter, tube adapter, tube clamp adapter). Observe the adapter instructions for use during adaptation and installation.
5	Orientation of the two prosthetic joints: joint axes parallel to each other
6	Orientation of hip prosthesis socket with lamination plate and prosthetic hip joint: lamination plate positioned horizontally hip prosthesis socket as far in anterior direction as possible on the lamination plate TMS line on the hip prosthesis socket parallel to the alignment reference line at a distance of 20 mm anterior orientation of the hip prosthesis socket according to the dorsal pelvis inclination marking orientation of the centre front of the lamination plate according to the frontal reference line pelvis rotation line 90° to the walking direction
7	INFORMATION The 646D1532=all_INT Quick Start Guide includes selection tables for easier selection of the tube clamp adapters. Join the prosthetic hip joint and the prosthetic knee joint using the chosen adapters (e.g. rotation adapter, tube clamp adapter, tube). Observe the adapter instructions for use during adaptation and installation.

Changing the length of the tube above the prosthetic knee joint changes the distance between the partial mass centre of gravity line and the alignment reference line. This distance subsequently has to be reset to the nominal value of 20 mm.

5.5 Static Alignment

Alignment Recommendations

(Schematic of static alignment for 3R60=HD see fig. 19 – the 646D1532=all_INT Quick Start Guide contains all schematics for static alignment with the 3R60, 3R106=HD and 3C98-3 C-Leg.)
If the prosthesis is fabricated with the 3B1 Genium prosthetic knee joint, the static alignment instructions in conjunction with the 7E10* prosthetic hip joint from the seminar documentation and the instructions for use for the prosthetic knee joint as well as X-Soft apply. Compliance with the specifications in the following alignment recommendations is also required.

Static Alignment Process	
+ = anterior placement / - = posterior placement (in relation to the load line)	
No.	Reference measurements and positions for static alignment, adjustment tasks
	Required materials and tools: 743L100 L.A.S.A.R. Posture Set the prosthetic hip joint to the factory defaults (see the chapter "Dynamic Trial Fitting" – see page 30). Tighten the screws to the specified torque values (see the chapter "Finishing the Prosthesis" – see page 34). Verify the length of the prosthesis on the patient.

Static Alignment Process

+ = anterior placement / - = posterior placement (in relation to the load line)

No.	Reference measurements and positions for static alignment, adjustment tasks
①	<p>To determine the load line, position the patient on the L.A.S.A.R. Posture as follows:</p> <ul style="list-style-type: none"> • Prosthetic foot (with shoe) on force plate (sufficient load: > 35% body weight) • Other foot (with shoe) on height compensation plate • The tips of the shoes are in one line • The load line should run approx. 50 mm (±10 mm) anterior to the ankle adapter screw (depending on the type and size of the prosthetic foot) as well as approx. 0 to ±10 mm through the TMS line.
②	<p>Optimise the static alignment solely by changing the plantar flexion on the foot adapter.</p> <p>a-p positioning of the alignment reference point (rotation axis of the prosthetic knee joint) to the load line:</p> <p>3R60=HD (anterior, lower axis): -10 mm</p> <p>3R106=HD (anterior, upper axis): -35 mm</p> <p>3C98-3: -30 mm</p>

5.6 Dynamic Trial Fitting

INFORMATION

- Please make sure you are thoroughly familiar with the adjustment possibilities and their effects! Only then you will be able to optimally adjust the prosthesis to the requirements of the prosthesis wearer.

5.6.1 Adjustments

INFORMATION

Secure the prosthetic hip joint on the lamination plate solely with the accessories which are included in the scope of delivery!

- 1) **Abduction/adduction** (see fig. 14): To adjust, loosen the lower and upper screws and turn the joint around the lower screw. Then tighten the screws again. For the required torque values and tools, see the table in the chapter "Finishing the Prosthesis" see page 34.
- 2) **Rotation** (see fig. 15): Mark abduction/adduction with a pen. To adjust internal (medial) and external (lateral) rotation, loosen the lower, middle and upper screws. Then tighten the screws again. For the required torque values and tools, see the table in the chapter "Finishing the Prosthesis" see page 34.
- 3) **Distal adjustment of the prosthetic hip joint** (see fig. 16, see fig. 17): The lower pyramid adapter permits tilting of the prosthetic hip joint in ML and AP directions. This is done for correct alignment of the prosthetic hip joint in the prosthesis.

5.6.2 Adjustment Possibilities

	Swing phase damping (SW) (see fig. 20)	Stance phase damping (ST) (see fig. 21)
Parameter	The free swing phase is followed by a significant increase in damping to limit the stride.	The extension movement can be damped in the stance phase, allowing comfortable extension of the prosthesis.

	Swing phase damping (SW) (see fig. 20)	Stance phase damping (ST) (see fig. 21)
Factory setting	High The adjustment screw (SW) is turned all the way to the right.	Low The adjustment screw (ST) is turned all the way to the left.

5.6.3 Gait Parameter Adjustment

⚠ CAUTION

Trial walking without safety precautions

Falling due to lack of patient experience or incorrect product adjustment

- ▶ For safety reasons, always have the patient take the first steps between parallel bars.

NOTICE

Use of incorrect tools

Product damage due to use of incorrect tools

- ▶ Only use the tools recommended in the instructions.

INFORMATION

Before adjusting the gait parameters of the prosthetic hip joint, also observe the instructions for use of the chosen prosthetic knee joint.

Proceed as follows to make the adjustments (see fig. 22):

> **Tools required:**

710H10=2X3 Allen key

- 1) Check the factory settings (chapter "Dynamic Trial Fitting" see page 30).
- 2) Reduce the high factory setting for swing phase damping so that flexion damping is clearly noticeable but still comfortable. This will make safe heel strike possible.
Now make adjustments to the prosthetic knee joint.
- 3) **Information:** Increasing stance phase damping of the prosthetic hip joint can result in stance phase flexion of the prosthetic knee joint. In case of excessive stance phase damping of the prosthetic hip joint, this can cause uncontrollable flexion of the prosthetic knee joint.
Gradually increase stance phase damping of the prosthetic hip joint and allow the patient to get accustomed to the damping value. Adjusting stance phase damping should be done over a longer period of time. With optimally adjusted damping, the prosthesis wearer will slowly sink into the prosthetic hip joint in the stance phase and reach maximum extension of the prosthetic hip joint just before initiating the swing phase.

5.6.4 Correcting Possible Problems

Troubleshooting with the 3R60=HD and 3R106 HD prosthetic knee joints

Observe the instructions for use of the chosen prosthetic knee joint for troubleshooting.

Troubleshooting for the 3C98-3 C-Leg and 3B1 Genium prosthetic knee joints

Basic troubleshooting information is provided below. Detailed information can be found in the instructions for use for the respective knee joint.

If the chosen prosthetic knee joint does not go into extension during the stance phase after stance phase flexion when the prosthetic hip joint is used, this can have several causes. Check the possible causes and, if applicable, make the specified adjustments to achieve extension of the prosthetic knee joint in the stance phase:

No extension in the stance phase	
Cause	Possible solution
The prosthetic alignment is not optimal.	Check and optimise the sagittal prosthetic alignment under load (see the chapter "Static Alignment" see page 29).
The stride is too long.	Increase swing phase damping.
Stance phase damping of the prosthetic knee joint is too low.	Stance phase damping of the prosthetic knee joint should be optimised for sitting down as well as for alternating negotiation of stairs and ramps.
Stance phase damping of the hip joint is too high.	Reduce stance phase damping.
The chosen functional prosthesis components are not optimal for the patient.	Check the selection and use different functional prosthesis components if needed.

If the chosen prosthetic knee joint also does not switch to swing phase mode at the end of the stance phase when used with the prosthetic hip joint, this can have the following causes:

No changeover to swing phase (3C98-3 C-Leg)	
Cause	Possible solution
The forefoot load value is not reached during rollover.	Check the forefoot load with C-Soft while walking and correct if required (C-Soft Instructions for Use: 647G268).

If the prosthetic hip joint does not go into extension at heel strike, this can have several causes. Check for the following possible causes and, if needed, make the specified adjustments to achieve extension of the prosthetic hip joint at heel strike:

Hip joint: no extension at heel strike	
Cause	Possible solution
Stance phase damping of the prosthetic hip joint is too high.	Reduce stance phase damping.
The stride is too long.	Increase swing phase damping.
The prosthetic alignment is not optimal.	Check and optimise the sagittal prosthetic alignment under load (see the chapter "Bench Alignment" see page 28).
The chosen functional prosthesis components are not optimal for the patient.	Check the selection and use different functional prosthesis components if needed.

INFORMATION

When adjusting stance phase damping, change the settings gradually and keep checking whether damping can be further increased. It is important for the prosthesis wearer that weight can be put on the prosthesis side as long as possible during the stance phase.

The stance phase movement on the moving prosthetic hip joint is possible because of high damping and should be considered the goal of trial walking, even if the prosthesis wearer might feel uncomfortable at first! It takes time to understand this important procedure and to explain and practise it with the prosthesis wearer. Only then will the prosthesis wearer be able to fully benefit from the advantages of this product.

5.7 Finishing the Trial Prosthesis

Maximum Abduction/Adduction

- 1) For screwing the lamination plate to the hip prosthesis socket in the frontal plane, use the bore holes in the lamination plate marked in the illustration (see fig. 9).
This will allow maximum abduction/adduction adjustment of the hip joint, if necessary (see fig. 10).
- 2) On the bottom side of the lamination plate, all bore holes can be used for the screw connection to the pelvic socket.

Maximum interior/exterior rotation – use of the support plate

INFORMATION

The 4G576 Support Plate is intended solely for use with the trial fitting. When fabricating the definitive hip prosthesis socket, ensure that the lamination plate is attached with sufficient exterior rotation and that the pelvic socket has a laminate thickness of at least 2.5 mm. If the lower section of the prosthetic hip joint nevertheless touches the hip prosthesis socket, the laminate can be ground at this location.

INFORMATION

Ensure that the anterior pressure block is positioned with approx. 5° of exterior rotation when plastering the pelvic socket.

- For the trial fitting period, use the supplied support plate (see fig. 13).
- In this way, the maximum internal and external rotation of the prosthetic hip joint can be set and used, without the posterior linkage bar contacting the lamination plate.

INFORMATION

The trial prosthesis is only intended for trial fitting to determine correct alignment. It may only be worn by the patient during this time and is not intended for extended use.

5.8 Reinforcing the Prosthetic Socket

INFORMATION

The following reinforcing instructions only refer to the connection and stability of the lamination plate in the hip prosthesis socket.

For patients with a body weight near the allowable top limit and with expected high mobility, it is recommended to apply at least 2 extra layers of carbon fibre mesh and 2 extra layers of nylon stockinette.

> Recommended materials:

- Nylon Stockinette, 616G12 Carbon Fibre Mesh, 616G4 Fibreglass Mat, 699B1 Fibreglass Roving, 636K8 Plasta Tape
- 1) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 1st and 2nd layers of nylon stockinette are applied.
 - 2) Apply carbon fibre mesh in the area of the sitting surface and lamination plate.
 - 3) Apply fibreglass mat in the area of the sitting surface and lamination plate.
 - 4) Apply carbon fibre mesh offset at 45° to the first layer of carbon fibre mesh in the area of the sitting surface and lamination plate.
 - 5) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 3rd and 4th layers of nylon stockinette are applied.
 - 6) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 5th and 6th layers of nylon stockinette are applied.
 - 7) Pull fibreglass roving through the holes in the cleaned lamination plate.

- 8) Close the threaded holes in the lamination plate with plasta tape.
- 9) Position the lamination plate on the plaster positive.
- 10) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 7th and 8th layers of nylon stockinette are applied.
- 11) Apply carbon fibre mesh in the area of the sitting surface and lamination plate.
- 12) Apply carbon fibre mesh offset at 45° to the first layer of carbon fibre mesh in the area of the sitting surface and lamination plate.
- 13) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 9th and 10th layers of nylon stockinette are applied.
- 14) Apply carbon fibre mesh in the area of the sitting surface and lamination plate.
- 15) Apply carbon fibre mesh offset at 45° to the first layer of carbon fibre mesh in the area of the sitting surface and lamination plate.
- 16) Pull 2 layers of nylon stockinette over the plaster positive.
→ The 11th and 12th layers of nylon stockinette are applied.

5.9 Finishing the Prosthesis

The 3S107* (for 3R60 + 3R106) or 3S27* Foam Cover (for 3C98-3 + 3B1) can be used to finish the prosthesis. When cutting to length, include a 60 mm compression allowance (lower leg length +30 mm, thigh length +30 mm).

Use the 710D4 Torque Wrench and secure the grub screws with 636K13 Loctite® when finishing the prosthesis. Observe the following torque values when finishing:

Screws / grub screws	Torque
501S72=M10* Flat Head Screws on the lamination plate (with socket heads)	35 Nm (hexagon screwdriver, 6 mm)
501F9=*M8 Flat Head Screws on the lamination plate (with star heads)	25 Nm (709Z11 T40 Bit)
501T39-M4* Cap Screw (area bottom plate – intermediate piece – lamination plate)	2 Nm (709S10=2.5 Hex Key)
Grub screws on the rotation adapter	10 Nm / torque above the rotation adapter
All other grub screws	15 Nm

Finally verify the gait pattern one more time while the patient walks with the prosthesis. If necessary, readjust the prosthesis to compensate for influences of the foam cover.

6 Gait Training

Only with adequate instruction and gait training can the patient fully benefit from the features of the prosthetic hip joint system. Gait training is essential for learning how to use the prosthesis safely. The following aspects are especially important:

- ▶ The prosthesis wearer steps down with the prosthetic hip joint in flexion, which is extended with damping during the stance phase.
→ Unlike with conventional prosthetic hip joints, the prosthesis wearer does not have to perform a quick backward tilt of the pelvis. Extension of the prosthetic hip joint can be dampened during the entire stance phase, thereby reducing an abrupt backward tilt of the pelvis upon weight bearing. This is initially unfamiliar for experienced prosthesis wearers in particular and has to be practiced consciously.

7 Cleaning

NOTICE

Use of unsuitable cleaning agents

Damage to the product (e.g. bearings, seals and plastic components) due to substances contained in the cleaning agent

- ▶ Do not use aggressive cleaning agents.

8 Maintenance

INFORMATION

Due to the use of control functions which are performed hydraulically, movement noises may result when exoprosthesis hip joints are used. Such noises are normal and inevitable. They generally do not cause any problems.

If movement noises increase noticeably over the lifecycle of the hip joint, the prosthetic joint should be inspected immediately by a prosthetist.

- After the prosthetic hip joint has been worn for several weeks, it should be inspected and any necessary adjustments should be made. After an adaptation phase, prosthesis wearers usually appreciate a slightly higher extension damping in the stance phase.
- Conduct annual safety inspections.
- Check the prosthetic hip joint for wear and proper function once a year and make adjustments if necessary. Special attention should be paid to movement resistance and abnormal noises. Full flexion and extension must be ensured.
- Do not disassemble the prosthetic hip joint. Please send the complete prosthetic joint to the manufacturer for repair in case of problems.

9 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return and collection processes.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical Data

Reference number	7E9
Proximal connection	Lamination Plate
Distal connection	Pyramid Adapter
Max. flexion angle	130°
Weight of the joint [g]	695

Reference number	7E9
Weight of connecting components [g]	470
System height [mm]	81
System height, horizontal [mm]	19
Material	Aluminium
Mobility grade	2 + 3
Max. body weight [kg]	125
Operating and storage temperature	-10°C to +60°C (14°F to 140°F)

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