
3S80, 3S80=1 Sport





DE Gebrauchsanweisung	7
EN Instructions for use	20
FR Instructions d'utilisation	33
IT Istruzioni per l'uso	48
ES Instrucciones de uso	62
PT Manual de utilização	76
NL Gebruiksaanwijzing	90
SV Bruksanvisning	104
DA Brugsanvisning	117
NO Bruksanvisning	130
FI Käyttöohje	143
HR Upute za uporabu	156
JA 取扱説明書	169

Grundaufbau Bench Alignment 1E90 Sprinter

Anfänger
Beginners

Fortgeschrittene
Advanced users




5


Grundaufbau Bench Alignment

1E91 Runner


1E93 Runner junior

743A200
743A220

743Y580=

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INFORMATION

Date of last update: 2020-05-14

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

INFORMATION

These instructions for use are only valid starting with the following serial number: 201429001

1.1 Construction and Function

The 3S80* Sport is a sport knee joint designed for sports with no double support phase while running. These include in particular jogging, sprinting and related running disciplines in athletics that encompass corresponding

movement sequences. For this purpose, the prosthetic knee joint has rotation hydraulics to control the swing phase.

The prosthetic knee joint is not suitable for an everyday prosthesis since it has no stance phase control and is designed for the considerably greater flexion angles that occur during sports.

It has adjustable extension and flexion damping for adaptation to the chosen sport. The prosthetic knee joint has an integrated lock.

The hydraulic oil used in the 3S80=1 has a considerably lower viscosity than the oil used in the 3S80. The less viscous oil facilitates smooth-running movements, what can be advantageous e. g. in sport prostheses for children and especially in sprint prostheses.

1.2 Combination Possibilities

Designation	Reference number
Lamination Anchor	4R111, 4R116
Socket Adapter	4R54, 4R77, 4R55, 4R51
Double Adapter	4R72*
Foot Adapter	1E90 Sprinter: 4R204, 4R206 1E91 Runner: 4R218, 4R216 1E93 Runner junior: 4R224
Prosthetic Feet	1E90 Sprinter, 1E91 Runner, 1E93 Runner junior
Foot Soles	1E90 Sprinter: 2Z500, 2Z501 1E91 Runner: 2Z540, 2Z541 1E93 Runner junior: 2Z543

INFORMATION

In the sport prosthesis, only use prosthesis components approved for the respective sport or a body weight of 150 kg.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Approved for a body weight up to **max. 100 kg**.

Permitted sports

For use in sport prostheses (recreational sports, competitive sports) for jogging, sprinting and related running disciplines in athletics

Support for other sports

- ▶ Contact Ottobock if a prosthesis is needed for sports other than the permitted sports. Ottobock is available for enquiries regarding application in other sports.

2.3 Environmental conditions

Allowable environmental conditions
Use: Temperature range: -10 °C – +45 °C
Transportation, Storage: Temperature range: -20 °C bis +60 °C, Relative humidity 20 % – 90 %
Cleaning required after contact: Dust, Sand, Fresh water, Soap solution, Chlorine water
Immersion depth: ≤2 m, Immersion time: ≤1 h

Unallowable environmental conditions
Mechanical vibrations, Impacts
Hygroscopic particles (e.g.Talcum powder), Salt water, Sweat, Urine, Acids, Sauna
Transportation: Mechanical vibrations, Impacts

2.4 Service life

Based on ISO 10328, the prosthetic components were tested by the manufacturer for 3 million load cycles.

The service life is established for the sports listed in the section “Field of application” (see Page 21) under the heading “**Permitted sports**”, as follows:

- Competitive sports: 1 year
- Recreational sports: 2 years

3 Safety

3.1 Explanation of Warning Symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

3.2 General Safety Instructions

WARNING

Using the product

Severe falls due to improper use of the product

- ▶ Only use the product after receiving full instruction by qualified personnel.
- ▶ Do not use the product as an everyday prosthesis.
- ▶ **Inform the patient.**

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Use the prosthesis components according to the classification (see section "Area of Application").
- ▶ **Inform the patient.**

CAUTION

Unallowable combination of prosthesis components

Fall due to breakage or deformation of the product

- ▶ Only combine the product with prosthesis components that are approved according to the section "Combination Possibilities".
- ▶ Based on the instructions for use of the prosthesis components, verify that they may be combined with each other.

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if applicable (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).
- ▶ **Inform the patient.**

⚠ CAUTION

Reuse on other patients and improper maintenance

Falling due to loss of functionality or damage to the product

- ▶ Only use the product for a single patient.
- ▶ Maintain the product regularly to extend the service life of the product (see the section "Maintenance").
- ▶ **Inform the patient.**

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.
- ▶ **Inform the patient.**

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check a damaged product for proper functionality and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take suitable measures if required (e. g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

Signs of changes in or loss of functionality during use

Noticeable functional changes may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

4 Scope of Delivery

The scope of delivery for the product is shown on page 1 (cover picture).

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as

single components (■), single components with minimum order quantities (▲) or single-component packs (●):

Fig.		Item No.	Quantity	Designation	Reference number
–	■	–	1	Instructions for Use	647G813=all_INT
–	■	–	1	Quick Reference Guide	646D1563=all_INT
–	–	–	1	Sport Knee Joint	–
–	■	–	1	Adjustment Wrench	710H10=2x3
–	■	–	2	Pyramid Dome Protector Kit	2Z11=KIT

5 Preparation for Use

CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthesis components

- ▶ Observe the alignment, assembly and adjustment instructions.

CAUTION

Use of talcum

Falling, damage to the prosthetic knee joint due to lack of lubrication

- ▶ Do not use talcum on the prosthetic knee joint or other prosthesis components.
- ▶ **Inform the patient.**

The three-dimensional alignment of the prosthesis socket and the modular components affects the function of the prosthesis. The patient can only benefit from the special characteristics of the prosthetic knee joint when correctly aligned.

The residual limb position, e.g., abduction, adduction and flexion, must be considered when positioning the socket connector. Proceeding according to the manufacturer's alignment guidelines is recommended.

Alignment is performed in 2 steps:

- Bench Alignment
- Static Alignment

The prosthesis and prosthetic knee joint are individually adjusted during dynamic trial fitting.

To select and adapt the 1E90 Sprinter prosthetic foot, a trial prosthesis is first fabricated followed by the definitive prosthesis by replacing the test foot adapter (TF: 4R206) with the sport foot adapter (TF: 4R204).

INFORMATION

For more information on alignment, installation and optimising the prosthesis, please see the 647H537 Quick Reference Guide as well as the 647G839 (4R204, 4R206 TF Foot Adapter), 647G849 (1E90 Sprinter Prosthetic Foot) and 647G848 Instructions for Use (2Z500, 2Z501 Outer Sole).

Use 2Z11=KIT

INFORMATION: The connection area of the prosthetic knee joint can be protected against scratches with the protective film of the 2Z11=KIT during alignment in the workshop and testing in the trial fitting area.

- ▶ Use the protective film as illustrated in the accompanying document for the 2Z11=KIT.
- ▶ Remove the protective film before the patient leaves the fitting area.

5.1 Bench Alignment

Bench alignment of the prosthesis: 1E90 Sprinter

Optimising the Prosthetic Foot Height

The sport prosthesis is fabricated as a trial prosthesis with the test foot adapter for initial bench alignment. This verifies whether the stiffness of the prosthetic foot was correctly selected. Then the prosthetic foot height of the 1E90 is optimised to meet the needs of the patient. Movements in the hip and especially the shoulders offer important information on the optimum prosthesis length. When the prosthesis length is optimised, sagging of the shoulder on the prosthesis side while running is avoided. Further information and exercises are found in the 647H537 Quick Reference Guide.

Bench alignment is performed with the installed 2Z500 or 2Z501 Outer Sole (spikes not installed!).

Bench Alignment Process

+ = anterior placement / - = posterior placement (to the alignment reference line)

Item	see fig. 4
	Required tools and materials: 743A80 50:50 Gauge, alignment apparatus (e.g., 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly) Prosthetic knee joint, locked (note the section "Use" – see Page 30)
①	a-p positioning of the prosthetic foot alignment reference point to the alignment reference line: +40 mm (beginner) to 0 mm (advanced)

Bench Alignment Process

+ = anterior placement / - = posterior placement (to the alignment reference line)

Item	see fig. 4
②	Exterior foot rotation: approx. 5°
③	Structural height of the prosthetic knee joint at the alignment reference point (rotation axis of the prosthetic knee joint): MTP-floor measurement + 20 mm
	a-p positioning of the alignment reference point to the alignment reference line: 0 mm
	Exterior rotation of the prosthetic knee joint: approx. 5°
④	Join the prosthetic foot and the prosthetic knee joint with the help of the foot adapter (for TF trial prosthesis: 4R206; for TF sport prosthesis: 4R204) and the socket adapter. The socket adapter can be positioned on the foot adapter in three positions in the a-p direction. The exterior anterior position is recommended for advanced users, the exterior posterior position for beginners. Observe the adapter instructions for use during adaptation and installation. Additional adapters (double adapter, socket adapter) can also be used for height adjustment (see the section "Combination Possibilities" – see Page 21).
⑤	Mark 2 centre points laterally on the prosthetic socket: 1. At the level of the trochanter major 2. In the distal region
⑥	Join the two points with a line.
⑦	Join the prosthetic socket and the prosthetic knee joint with the help of the lamination anchor. Observe the lamination anchor instructions for use during adaptation and installation. Positioning of the prosthetic socket to the alignment reference line: <ul style="list-style-type: none"> • a-p position: 0 mm (proximal point of the prosthetic socket) • Socket flexion: hip flexion contracture (determined with Thomas test) + at least 4° Take the individual situation of the patient into account.

Bench Alignment Process	
+ = anterior placement / - = posterior placement (to the alignment reference line)	
Item	see fig. 4
	Additional adapters (double adapter, socket adapter) can also be used for height adjustment (see the section "Combination Possibilities" – see Page 21).

Bench alignment of the prosthesis: 1E91 Runner Runner, 1E93 Runner junior

- see fig. 5
- ▶ Please read this document carefully: 647G1145

5.2 Static Alignment

Static Alignment Process	
+ = anterior placement / - = posterior placement (to the load line)	
Item	1E90 Sprinter: see fig. 6; 1E91 Runner/1E93 Runner junior: see fig. 7
	Required materials and tools: 743L100 L.A.S.A.R. Posture Prosthetic knee joint, locked (note the section "Use" – see Page 30)
①	To determine the load line, position the patient on the L.A.S.A.R. Posture as follows: <ul style="list-style-type: none"> • Prosthetic foot on force plate (sufficient load: > 35% body weight) • Other foot (with shoe) on height compensation plate • The tips of the shoe and prosthetic foot are in one line
②	Optimise the static alignment exclusively by changing the plantar flexion on the socket adapter of the foot adapter. a-p positioning of the alignment reference point (rotation axis of the prosthetic knee joint) to the load line: <ul style="list-style-type: none"> • 1E90 Sprinter: Beginner: > -80 mm Advanced: -80 mm • 1E91 Runner/1E93 Runner junior: ~ 40 mm

5.3 Dynamic Trial Fitting

In the course of dynamic trial fitting, the prosthesis alignment is checked, plantar flexion is adjusted and flexion damping as well as extension damping

are configured. The correct use of the sport prosthesis must be practised to ensure safe running and walking. This section contains fundamental information on the exercises. See the 647H537 Quick Reference Guide for more details.

Walking Exercises

Use parallel bars for safety when trial walking with the prosthesis. Begin exercises between parallel bars after each change to the prosthesis. Always start with slow, short steps.

In order to help the patient develop a feel for the prosthetic foot, start the walking exercises with the prosthetic knee joint locked.

Then unlock the prosthetic joint in order to develop a feel for its functionality. Increase the stride length and walking speed with experience.

Conduct the walking exercises until the patient feels safe while walking.

Running Exercises

Use suitable protective equipment (e.g., for inline skating) for the initial running exercises. Choose a surface that is not too hard (e.g., tartan track, hall or short, firm lawn), since a hard surface (e.g., asphalt or concrete) can put greater strain on the patient and may cause pain.

Adapt the prosthesis to the sport and needs of the patient with corresponding adjustments and settings.

Finishing the Prosthesis

After successful optimisation, the trial prosthesis is converted to the definitive prosthesis by replacing the test foot adapter with the sport foot adapter. Further information is found in the instructions for use for the prosthetic foot and foot adapter.

5.3.1 Adjusting Flexion and Extension Damping

Valve	Symbol	Factory setting	
Flexion valve	(F)	Stop +	High
Extension valve	(E)	Stop –	Low

> Required tools and materials:

710H10=2x3 Adjustment Wrench

- 1) Use the adjustment wrench to adapt the flexion resistance to the needs of the patient.
- 2) Use the adjustment wrench to adapt the extension resistance to the needs of the patient.
- 3) Check the adjustments in dynamic trial fitting and correct as required according to the table below.

Situation	Measure	Adjustment
Lower leg swings through too far	Increase damping	Turn flexion valve to the right (+)
Lower leg does not swing through sufficiently	Reduce damping	Turn flexion valve to the left (-)
Lower leg impact too hard on the extension stop	Increase damping	Turn extension valve to the right (+)
Lower leg does not reach full extension before heel strike	Reduce damping	Turn extension valve to the left (-)

6 Use

CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthesis components

- ▶ Do not touch overheated prosthesis components.
- ▶ In case of changes in functionality, stop all activities and allow the overheated prosthesis components to cool down.
- ▶ In case of changes in functionality, reduce strain on the prosthesis component in order to allow cooling.
- ▶ In case of overheating or changes in functionality, have the prosthesis component inspected for damage by a service centre authorised by the manufacturer.
- ▶ **Inform the patient.**

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take suitable measures if required (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

INFORMATION

Secure Standing

The prosthetic joint is equipped with a lever to secure it when standing (see fig. 3).

- ▶ While standing under load (full extension), push the lever down to activate the standing lock.
- ▶ While standing under load (full extension), pull the lever up to deactivate the standing lock.

- After each use, check the prosthetic knee joint for wear and proper functionality.

6.1 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

CAUTION

Failure to follow the maintenance instructions

Injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the maintenance instructions.
- ▶ **Inform the patient.**

- The prosthesis components should be inspected after the first 30 days of use.
- Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic knee joint and adapt them to the patient's requirements again as needed.
- During the consultation, listen for unusual noises and inspect the entire prosthesis for damage, wear and functional changes.
- Arrange regular maintenance intervals with the patient depending on the level of use.
- Conduct annual safety inspections.

- As part of the safety inspections, inspect the prosthetic knee joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

Direct the patient to observe the following maintenance instructions:

- Check the entire prosthesis for damage after each use.
- In case of changes in or loss of functionality and damage, do not continue using the prosthesis and have it checked by qualified personnel.
- Comply with the maintenance intervals.

Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Only clean the product with the approved cleaning agents.
- ▶ Only disinfect the product with the approved disinfectants.
- ▶ Observe the instructions for cleaning and care.
- ▶ **Inform the patient.**

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical Data

Reference number	3S80, 3S80=1 Sport
Weight (g)	682
System height (mm)	48
Proximal system height up to the alignment reference point (mm)	28
Distal system height up to the alignment reference point (mm)	20
Proximal connection	Pyramid adapter
Distal connection	Pyramid adapter
Max. flexion angle (°)	135
Max. body weight (kg)	100

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