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

### INFORMATION

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ISO 22675 – “P” – “m” kg\*) 

- \*) Body mass limit not be exceeded!  
 For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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## INFORMATION

Last update: 2017-02-21

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

## 1 Preface

The main features of the 1C64 Triton Heavy Duty carbon fibre foot are the lightweight design and the unique interconnected functional elements which allow amputees with mobility grade 3 or 4 to walk easily in everyday life situations including participation in high energy activities.

The 1C64 is designed for use in modular prostheses for patients with different amputation levels.

The combination of spring elements offers the following functional advantages:

- Comfortable heel strike with noticeable plantar flexion.
- Smooth movement of the ankle during mid-stance phase for the natural rollover of the foot.

- Dynamic transition from the stance to swing phase.
- Controlled forefoot response and energy return.
- Compensation for uneven terrain.
- Particularly robust and corrosion-resistant design.

The foot is designed to be used with the footshell. The foam connection cap provides effective connection between the foot and the cosmetic foam cover.


### 1.1 Intended use

The 1C64 Triton Heavy Duty prosthetic foot is exclusively used for lower limb prosthetic fittings.

### 1.2 Field of application


Area of application according to Ottobock's MOBIS classification:

**1C64 = 21 – 24 cm**

Recommended for amputees with **Mobility Grade 3 and 4** (unrestricted outdoor walkers and unrestricted outdoor walkers with particularly high demands).  
**Approved for patient weight of up to 100 kg (220 lbs).**

**1C64 = 25 – 30 cm**

Recommended for amputees with **Mobility Grade 3 and 4** (unrestricted outdoor walkers and unrestricted outdoor walkers with particularly high demands).  
**Approved for patient weight of up to 150 kg (330 lbs).**

**1C64 = 25 – 30 cm**

Recommended for amputees with **Mobility Grade 3** (unrestricted outdoor walkers).  
**Approved for patient weight of up to 150 kg (330 lbs).**

Please refer to the selection chart below to determine the appropriate stiffness required according to Ottobock recommendations.

Body Weight		Foot Size																			
		21 cm		22 cm		23 cm		24 cm		25 cm		26 cm		27 cm		28 cm		29 cm		30 cm	
kg	lbs	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal	slim	normal
up to 55	up to 121	1	-	1	-	1	-	1	1	1	1	1	1	-	-	-	-	-	-	-	-
56 – 75	122 – 165	2	-	2	-	2	-	2	2	2	2	2	2	2	2	-	2	-	2	-	2
76 – 100	166 – 220	3	-	3	-	3	-	3	3	3	3	3	3	3	3	-	3	-	3	-	3
101 – 125	221 – 275	-	-	-	-	-	-	-	4	4	4	4	4	4	4	-	4	-	4**	-	4**
126 – 150	276 – 330	-	-	-	-	-	-	-	5	5	5	5	5	5*	5*	-	5*	-	5*	-	5*

\* When combining this configuration with the C-Leg, please contact Ottobock Customer Service.

\*\* When combining this configuration with the C-Leg or C-Leg compact, please contact Ottobock Customer Service.

## 2 Safety instructions

Please inform your patients about these safety issues!

### 2.1 Explanation of warning symbols

 <b>CAUTION</b>	Warnings regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warnings regarding possible technical damage.
 <b>INFORMATION</b>	Additional information on the fitting/ use.

### 2.2 Safety instructions for patients/users

 <b>CAUTION</b>
<p><b>Risk of injury as a result of excessive strain.</b></p> <p>Excessive impact loading typical with competitive sports beyond running or jumping activities, may weaken the spring elements.</p> <p>To avoid the risk of accident and possible injury to the patient, the foot must no longer be used after a noticeable loss of function has occurred.</p> <p>This noticeable function loss may be the result of reduced spring effect or forefoot resistance, in some cases accompanied with abnormal noises.</p> <p><b>Measure:</b> Have the prosthetic foot examined and, if necessary, replaced at a specialist workshop.</p>

 <b>NOTICE</b>
<p><b>Damage due to incorrect environmental conditions.</b></p> <p>Prosthetic components must not be exposed to the following elements:</p> <ul style="list-style-type: none"><li>• Corrosive elements (e. g. acids)</li><li>• abrasive substances (e. g. sand)</li></ul> <p>The footshell material is not UV-resistant.</p> <p>Do not grind the surface of this prosthetic foot, as this can lead to premature wear.</p> <p>Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.</p>

## 3 Product description

### 3.1 Scope of delivery

The 1C64 Triton Heavy Duty Prosthetic Foot is delivered with the following components:

- 1C64 Triton Heavy Duty Foot Module
- 2F60 Heel Wedge Kit
- SL = Spectra-Sock Protective Sock
- 647G824 Instructions for Use

### 3.2 Components (see Fig. 2)

- Titanium adapter (1)
- Carbon Fibre Attachment Spring (2)
- Optional Heel Wedge (3)
- Carbon Fibre Heel Spring (4)

- High Performance Polymer Base Spring (5)
- Carbon Fibre Forefoot Spring (6)

### 3.3 Single components ■ (see Fig. 1)

#### 2C6 footshell including connection cap

##### 2C6=\*S Slim footshell

2C6=(Side)(Size) / (Colour)S

Order example: 2C6=L27/4S

<b>Sides</b>	left (L), right (R)									
<b>Sizes (cm)</b>	21	22	23	24	25	26	27	–	–	–
<b>Colours</b>	beige (4), light brown (15)									

##### 2C6=\*N Normal footshell

2C6=(Side)(Size) / (Colour)N

Order example: 2C6=L27/4N

<b>Sides</b>	left (L), right (R)									
<b>Sizes (cm)</b>	–	–	–	24	25	26	27	28	29	30
<b>Colours</b>	beige (4), light brown (15)									

### 3.4 Single component pack (see Fig. 1)

#### 2F60 Triton Heel Wedge Kit

Order example: 2F60=24–25

<b>Sizes (cm)</b>	21 – 22	23 – 24	25 – 26	27 – 28	29 – 30
-------------------	---------	---------	---------	---------	---------

## 4 Assembly

### 4.1 Setup

#### 4.1.1 Footshell application / removal

#### NOTICE

##### Risk of damage as a result of improper use.

The foot is designed to be used with the footshell. Please avoid frequent removal from and insertion into the footshell as this can damage the carbon fibre elements.

#### NOTICE

##### Risk of damage as a result of improper footshell assembly.

Do not use a screwdriver for assembly and disassembly as this can cause damage to the carbon elements and/or footshell. Ottobock recommends use of 2C100 Footshell Insertion/Removal tool (accessories).

**NOTICE****Risk of damage as a result of improper use.**

The footshell must not be dusted with talcum powder as this would increase abrasion. Use silicone spray for lubrication if needed. Replace the footshell in case it shows any damage or wear.

**Inserting the foot into the footshell:**

- Attach the tube adapter to the prosthetic foot. This helps, makes it easier to insert the foot into the footshell.
- Insert the foot into the footshell.
- Make sure the heel of the polymer base spring securely snaps into the footshell heel tab.

**Removing the footshell:**

- Insert the 2C100 Footshell Insertion/Removal tool under the heel of the polymer base spring which is located under the heel tabs in the posterior part of the footshell
- Pry the heel out of the shell and pull the foot out of the footshell.

**5 Handling****5.1 Alignment****5.1.1 Bench alignment**

The prosthetic foot with cover is clamped in the alignment apparatus.

**TT-fitting**

Foot size [cm]	Mid-foot anterior to the alignment reference line [mm]	Heel height in alignment apparatus
21–25	25	Effective heel height (shoe heel height – sole thickness of forefoot) + 5mm
26–28	30	
29–30	35	

**TF-fitting**

Observe the alignment instructions for the chosen prosthetic knee joint.

Effective shoe heel height [mm]	Suitable footshell
5 – 15	Normal
10 – 20	Slim

**5.1.2 Static alignment optimization**

If necessary, the alignment recommendations (TF modular leg prostheses: 646F219\*, TT modular leg prostheses: 646F336\*) may be requested from Ottobock.

Use the L.A.S.A.R. Posture for static alignment optimisation of the prosthesis, if available.

**5.1.3 Dynamic alignment optimization**

Adapt the prosthetic alignment in the frontal plane (ML) and in the sagittal plane (AP) by angle adjustments and/or sliding adjustments, in order to ensure correct heel contact, easy roll-over, and optimal shifting of weight onto the contralateral side. For transtibial amputees please make sure there is physiological knee flexion during stance phase.

## 5.2 Adjustment and final assembly

Due to the high forces in the adapter section, it is indispensable to tighten the adjustment screws of the tube adapter with a torque of **15 Nm (135 lb-in)**.

Use 710D1 Torque Wrench. Once you have attained the final setting for the prosthesis, use 636K13 Loctite® to secure the set screws.

## 5.3 Heel wedges

### CAUTION

#### **Risk of falling due to change of the knee dynamics when using heel wedges.**

The use of heel wedges can facilitate stance phase flexion when bearing weight on the prosthesis. Therefore, particular caution is advised during the dynamic alignment optimization – especially when a knee joint without stance phase control is used.

In case the heel strike is too soft or the knee is in hyperextension during the mid-stance phase, the heel can be stiffened by means of one of the two heel wedges delivered with each 1C64 prosthetic foot.

- “Soft“ Heel Wedge (translucent)
- “Stiff“ Heel Wedge (colour: black)

### 5.3.1 Heel Wedge application

#### INFORMATION

Patients with a high activity level in MG 3 should initially test the foot without a heel wedge. For patients with a high activity level in MG 4, Ottobock recommends the application of a heel wedge.

#### INFORMATION

The heel wedge design does not require any adhesive to apply it to the prosthetic foot.

#### INFORMATION

The heel wedge is properly installed when the locking recess in the heel wedge is fully seated around the mating stop (Fig. 3, Pos. 7) on the attachment spring.

1. Orient the heel wedge (Fig. 3, Pos. 3) for insertion so the locking recess in the wedge is positioned in the proximal and posterior direction.
2. Slide the heel wedge into position from the posterior side of the prosthetic foot between the attachment spring (Fig. 3, Pos. 2) and heel spring (Fig. 3, Pos. 4).

### 5.3.2 Heel Wedge removal

#### INFORMATION

The heel wedge has a tendency to rotate during removal. Pulling it in the direction of rotation makes it easier to remove the wedge.

1. Push on the medial or lateral side of the heel wedge with your finger to free it from the mating stop (Fig. 3, Pos. 7) on the attachment spring (Fig. 3, Pos. 2).

2. After the heel wedge is partially free grab it by the respective lateral or medial side and pull it free from the prosthetic foot.

#### **5.4 Foam cover**

Both Footshell Connection Caps 2C20 (slim) and 2C19 (normal) serve as a detachable connection between the foam cover and the footshell. They snap into the proximal brim of the footshell.

1. Cut the foam cover to length, allowing for distal compression about 0.4 in (10 mm) for transtibial and about 1.6 in (40 mm) for transfemoral covers.
2. With the connection cap attached to the footshell, pull the foam cover over the prosthesis, reassemble the foot and mark the outer contours of the connection cap.
3. Remove the connection cap and clean the foam cover with Isopropyl Alcohol 634A58.
4. Apply 636N9 Ottobock Contact Adhesive or 636W17 Plastic Adhesive to the distal end of the foam cover and to the connection cap.
5. Allow the adhesive to dry (approx. 10 minutes) and cut out a space in the distal foam cover with a grinding and milling machine to allow for precise fit for the upper part of the adapter. Sufficient space should be available for the adapter to ensure proper function and eliminate noise during walking.
6. Assemble the foam cover and finish off the outer shape of the foam cover in usual way with carbon foot installed. Any compression of the foam cover by cosmetic stockings or SuperSkin needs to be taken into account.

## **6 Cleaning and care**

1. Clean with a damp, soft cloth and soap when needed.
2. Rinse the prosthesis component with clear fresh water.
3. Dry the prosthesis component with a soft cloth.

## **7 Maintenance instructions**

- Foot assembly should be inspected after 30 days (4 weeks) of use.
- Inspection of the entire prosthesis for wear during normal consultations.
- Footshell may require replacement if wear is excessive.

## **8 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **8.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **8.2 CE conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

## INFORMATION

Explanation of the warning symbol :

**ISO 22675 – “P” – “m” kg\*) **

-  \*) Body mass limit not be exceeded!  
 For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned on the left, test levels (P) are assigned to certain maximal body masses (m in kg). In some cases, which are marked with , no test level is assigned to the product-related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified test level.

## INFORMATION

This component has been tested according to ISO 22675 standard for two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks.

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**ISO 22675 – “P” – “m” kg\*) △**

- \*) Body mass limit not be exceeded!
- △ For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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**The 1C64 Triton Heavy Duty Prosthetic Foot is covered by the following patents:**

European Patent: EP 2 509 540 in DE, FR, GB, IT, SE, TR  
European Patent: EP 2 522 311 in DE, FR, GB, IT, NL, PL, SE, TR  
China: ZL201080053507.3  
Japan: Patent No. 5 781 531  
Russia: RU 2 550 003  
USA: Patent US 9,486,331

Patents pending in Brazil, China, and the USA.

**The 1C64 Triton Heavy Duty Prosthetic Foot is covered by following registered designs and design patent:**

European Design: No. 001810425  
USA: Patent US D653,759 S  
IR-Design: DM/076079 in TR

Otto Bock HealthCare GmbH  
Max-Näder-Straße 15 · 37115 Duderstadt / Germany  
T +49 5527 848-0 · F +49 5527 72330  
healthcare@ottobock.de · www.ottobock.com

Otto Bock HealthCare LP  
3820 West Great Lakes Drive  
Salt Lake City, UT 84120 · USA  
T +1 800 328 4058 · F +1 800 655 4963