Pheon 3R62, 3R62=KD, 3R62=ST

Pheon 3R62=1, 3R62=1-KD, 3R62=1-ST

DE Gebrauchsanweisung	15
EN Instructions for use	24
FR Instructions d'utilisation	32
IT Istruzioni per l'uso	41
ES Instrucciones de uso	51
PT Manual de utilização	60
NL Gebruiksaanwijzing	69
SV Bruksanvisning	78
DA Brugsanvisning	86
NO Bruksanvisning	95
FI Käyttöohje	103
PL Instrukcja użytkowania	112
HU Használati utasítás	121
CS Návod k použití	129
HR Upute za uporabu	138
TR Kullanma talimatı	146
RU Руководство по применению	155
JA 取扱説明書	164
ZH 使用说明书	173





1 **i** 5 Lieferumfang/Scope of delivery

3R62, 3R62=1

3R62=ST, 3R62=1-ST

3R62=KD, 3R62=1-KD

2Z11=KIT


2R49 (299 mm)

647G876=all_INT

710H10=2X3

646D1571=all_INT

21A37=3

21A17=
3.8X1000
(50 mm)

507S15

21A6

21A5

503F3

4G114=1

21A45=0.83

507S15


516S3

21A12


4F17=N

4G12






4G70


506G3=M8X12-V









2 i 6.2 Grundaufbau der Prothese/Bench alignment of the prosthetic (1/2)

3 i 6.2 Grundaufbau der Prothese/Bench alignment of the prosthesis (2/2)

4 **i** 6.2 Montieren des Sperren-Entriegelungszugseils/Installing the lock release cable (1/2)

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INFORMATION

Date of last update: 2021-10-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.



This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings. Ottobock seminars, information materials and service are available for continuing education and questions (for contacts, see the manufacturer's address at the end of this document).

1.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for a defined sequence		Numbering for the parts of an illustration
	Right		Wrong

2 Product description

The product (3R62*, 3R62=1*) has the following key features:

- Polycentric prosthetic knee joint
- Product components for stance phase stability:
 - Locking mechanism for the extension position
 - Adjustable EBS unit for limited knee flexion at heel strike
- Product components for swing phase control:
 - Extension assist (adjustable spring force)
 - Adjustable friction brake as movement resistance for extension and flexion

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

Unallowable amputation types

Hip disarticulation, Hemipelvectomy

3.2 Area of application

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

Allowable field of application (3R62, 3R62=KD, 3R62=ST)

Recommended mobility grade: **1 + 2**

Everyday prosthesis

Allowable body weight: **75 – 125 kg**

Allowable field of application (3R62=1, 3R62=1-KD, 3R62=1-ST)

Recommended mobility grade: **1 + 2**

Everyday prosthesis

Allowable body weight: **45 – 75 kg**

Allowable field of application (2R49)

Allowable body weight: **≤ 125 kg**

3.3 Combination possibilities

CAUTION

Improper combination of prosthetic components

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's area of application.

Unallowable combinations

With use of the lock release cable: 4R57* Rotation adapter

4G702 Spring travel limiter

3.4 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range: -10 °C – +60 °C

Use - Relative humidity: 0 % – 90 % (Non-condensing)

Storage/Transportation - Relative humidity: 100 %

Tyre pressure: 250 – 1100 mbar

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

3.5 Reuse and lifetime

⚠ CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

⚠ CAUTION

Exceeding the lifetime

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the maximum lifetime defined in this section is not exceeded.

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

4 General safety instructions

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Noticeable changes in functionality may include poor response, incomplete extension, diminished swing phase control or stance phase stability, abnormal noises, etc.

5 Scope of delivery

For the scope of delivery, see Figure **1** on page 2.

A reference number shown in bold in the illustration stands for a single component pack, in which the product components shown below it can be reordered.

A product component without a reference number in the illustration cannot be reordered separately.

6 Preparing the product for use

6.1 Information on fabrication of a prosthesis

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

CAUTION

Initial use of the prosthesis by the patient

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

6.2 Bench alignment of the prosthesis

INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images for the sagittal plane (Fig. [2](#) on page 3) and for the frontal plane (Fig. [3](#) on page 4) as well as the instructions for use of all prosthesis components that were used.

Tube adapter

CAUTION

Incorrect processing of tube

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

CAUTION

Incorrect mounting of the tube

Risk of injury due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

Installing the lock release cable

CAUTION


Incorrect installation of the lock release cable

Risk of falling due to malfunction of the lock

- ▶ Observe the following instructions for installing the lock release cable.
- ▶ After installation, check whether the lock release cable has been installed in accordance with all instructions.
- ▶ Verify the safe functioning of the lock release cable and the lock (unlocking/locking).


Text in the illustration

 4]: Quetschgefahr! \triangleq Risk of pinching!; Safe pulling direction \triangleq Safe pulling direction

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- ▶ When a prosthetic knee joint with lock is used, install the lock release cable according to the illustrations  4 on page 5 and  5 on page 6.
- ▶ **INFORMATION:** Keep the screws you have removed (see Fig.  4, step  3) in case the lock needs to be permanently deactivated later on.
- ▶  4 In order for the lock pin to function reliably, install the lock release cable in the area marked for a safe pulling direction.
- ▶ To avoid malfunctions, set the length of the lock release cable so that it is not under tension when the lock slide is in the rest position. Reliable unlocking when operating the lock slide is essential.
 - Operating the lock slide: see Fig.  5, item  1
 - Rest position of the lock slide: see Fig.  5, item  2
- ▶  9 Do not break out the nose in the lock slide because the locking position of the lock slide is not allowed to be used.
- ▶  11 Check the reference dimension to ensure reliable unlocking (marking on the locking pin) according to the illustration.
INFORMATION: The reference dimension must be maintained so that the locking pin can be pulled out far enough to reliably lock the prosthetic knee joint.
- ▶ If the dimension is too small, shorten the lock release cable until the reference dimension is met when operating the lock slide.

6.3 Optimising the static alignment

- ▶ To optimise the static alignment, position the patient on the measuring device (see Fig.  6 on page 7).
- ▶ In the sagittal plane, check the course of the load line relative to the marked alignment reference point and correct this if needed (see Fig.  6 on page 7).
- ▶ In the frontal plane, check the course of the load line relative to the marked points (prosthetic foot, prosthetic knee joint, Spina iliaca anterior superior) and correct if needed (see Fig.  7 on page 8).

6.4 Optimising during dynamic trial fitting

CAUTION

Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

Checking the factory settings

Text in the illustration

[8]: Vorbringer \triangle Extension assist; Friktionsbremse \triangle Friction brake; EBS-Einheit \triangle EBS unit

- ▶ Prior to the dynamic trial fitting, restore the factory settings (see Fig. [8] on page 9).
- ▶ If a suitable torque wrench is not available, first set the adjustment screw of the friction brake so that high axle friction is felt and then slowly turn it back to minimal axle friction.
- ▶ **CAUTION!** For use as a locking knee joint, leave the friction brake set to minimal axle friction so the patient can safely use the locking function of the prosthetic knee joint.

The following settings on the prosthetic knee joint are required for use as a free-running knee joint.

Adjusting the EBS unit

INFORMATION

The EBS function must be explained to the patient in detail. Through intensive training (active loading and unloading of the prosthesis side while standing between parallel bars), the patient must become familiar with the functionality of the EBS and find the optimum setting by trial and error.

Text in the illustration

[9]: Erhöhen des Standphasenflexionswiderstands \triangle Increasing the stance phase flexion resistance; Standphasenflexion verringert \triangle Reduced stance phase flexion; Verringern des Standphasenflexionswiderstands \triangle Reducing the stance phase flexion resistance; Standphasenflexion erhöht \triangle Increased stance phase flexion

- ▶ For adjusting the EBS unit, note Fig. [9] on page 10.
- ▶ **CAUTION!** Explain the EBS function to the patient and familiarise the patient with the EBS between parallel bars.
- ▶ With the adjustment wrench, adjust the EBS so that the prosthetic knee joint goes into stance phase flexion at heel strike, not too fast and not too slow in accordance with the needs of the patient.

Adjusting the extension assist

Text in the illustration

[10]: Erhöhen der Federvorspannung \triangle Increasing the spring tension; Schwungphasenflexion verringert \triangle Reduced swing phase flexion; Schwungphasenextension verstärkt \triangle Stronger swing phase extension; Verringern der Federvorspannung \triangle Reducing the spring tension; Schwungphasenflexion erhöht \triangle Increased swing phase flexion; Schwungphasenextension schwächer \triangle Weaker swing phase extension

- ▶ For adjusting the extension assist, note Figure [10] on page 11.
- ▶ **CAUTION!** The lower leg of the prosthesis must not swing through too far in the dorsal direction at the beginning of the swing phase.

Adjust the extension assist so that the lower leg of the prosthesis swings through in the dorsal direction according to the patient's requirements, but reaches full extension in a timely manner for the next heel strike.

Adjusting the swing phase friction

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- ▶ For adjusting the swing phase friction, note Figure 11 on page 12.
- ▶ **CAUTION!** The prosthetic knee joint has to reach full extension, even at slower walking speeds. Observe the information under the heading "Fine-tuning the settings".
- ▶ Adjust the swing phase friction so that the prosthetic knee joint does not swing too hard against the extension stop and so that the prosthetic foot swings through sufficiently according to the needs of the patient.

Fine-tuning the settings

- ▶ **CAUTION!** The various settings of the prosthetic knee joint cannot be adjusted entirely independently of each other. When the settings cannot be fully adjusted to the comfort needs of the patient, the settings should be established primarily based on safety aspects.
- ▶ Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.
- ▶ If the extension position is not reached with every step and at any walking speed, reduce the braking effect and/or increase the initial tension of the extension assist.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and the annual safety inspection.
Advise the patient to check the prosthesis if changes in function occur.

6.5 Attaching the cosmetic cover

CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

- ▶ Note the following points for attaching the cosmetic cover:
 - The cosmetic cover must not interfere with the movement of the lock release cable.
 - After completion, verify functioning of the lock.
 - Check the settings from the dynamic trial fitting.
- ▶ To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

6.6 Finishing the prosthesis

- ▶ **CAUTION!** To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long, and by tightening all screw connections. In doing so, note the instructions for use for all prosthesis components regarding torque values and thread lock.

7 Use

7.1 Information for use

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

7.2 Using the lock slide

CAUTION

Malfunctions and improper use of the lock

Falling due to unexpected activation or deactivation of the prosthetic knee joint lock

- ▶ Fully extend the prosthetic knee joint when you want to activate the lock. Then carefully check whether the lock has actually engaged.
- ▶ Only deactivate the lock when you are able to control flexion of the prosthetic knee joint.
- ▶ Ensure that the lock is not activated or deactivated unintentionally while walking.

7.3 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

Cleaning light soiling

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

Cleaning heavier soiling

- > **Required materials:** cleaning cloth, 634A58 isopropyl alcohol, soft cloth
- ▶ **NOTICE!** Verify material compatibility. Do not remove lubricant from bearing points.
- ▶ Clean the product with a cleaning cloth and isopropyl alcohol.
- ▶ Dry the product with a cloth.

Cleaning with Disinfectants

- > **Required materials:** colourless, alcohol-free disinfectant (verify material compatibility!), soft cloth
- ▶ Disinfect the product with the disinfectant.
- ▶ Dry the product with a cloth.
- ▶ Allow to air dry in order to remove residual moisture.

8 Maintenance

CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.
- ▶ **NOTICE!** Do not lubricate and grease the prosthetic joint.
- ▶ **NOTICE!** Repair work must be performed exclusively by the manufacturer's technical service.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

- ▶ Check the initial tension settings (EBS unit, extension assist) to make sure they are free of play. Increase the initial tension if needed.

Replacing the stops

- ▶ Check the stops shown in Fig. 12 on page 13 for wear and integrity every time the prosthesis is inspected.
- ▶ Check the serial number according to the steps ① – ③).
- ▶ **CAUTION!** If the serial number is smaller than shown in the illustration, only have the stops replaced by the manufacturer's technical service.
- ▶ If the serial number is equal to or greater than shown in the illustration, replace the stops according to the subsequent steps ④ – ⑥).
- ▶ **CAUTION!** If the control measurement in the illustration is exceeded, do not use the product. Contact the manufacturer's technical service.

9 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

10 Legal information

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

11 Technical data

For the technical data, see Figure 13 on page 14.

Abbreviations in the technical data					
①	h _{s-prox}	Proximal system height	②	h _{s-dist}	Distal system height

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[9]: Increasing the stance phase flexion resistance ⚠ 提高站立期屈曲阻力; Reduced stance phase flexion ⚠ 站立期屈曲已减小; Decreasing the stance phase flexion resistance ⚠ 减小站立期屈曲阻力; Increased stance phase flexion ⚠ 站立期屈曲已增大

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[10]: Increasing the spring tension ⚠ 提高弹簧预加应力; Reduced swing phase flexion ⚠ 摆动期屈曲已减小; Stronger swing phase extension ⚠ 摆动期伸展已增强; Reducing the spring tension ⚠ 减小弹簧预加应力; Increased swing phase flexion ⚠ 摆动期屈曲已提高; Weaker swing phase extension ⚠ 摆动期伸展更弱

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