
3R46

DE Gebrauchsanweisung	5
EN Instructions for use	14
FR Instructions d'utilisation	22
IT Istruzioni per l'uso	32
ES Instrucciones de uso	41
PT Manual de utilização	50
NL Gebruiksaanwijzing	58
SV Bruksanvisning	67
DA Brugsanvisning	75
PL Instrukcja użytkowania	84
HU Használati utasítás	92
CS Návod k použití	101
RO Instrucțiuni de utilizare	109
HR Upute za uporabu	118
TR Kullanma talimatı	126
EL Οδηγίες χρήσης	134
RU Руководство по применению	144
JA 取扱説明書	153
ZH 使用说明书	161

- Einzelteile-Pack 4D18
- Single-component pack 4D18
- Pack de pièces détachées 4D18
- Kit componenti 4D18
- Kit de componentes 4D18
- Conjunto de componentes 4D18
- Onderdelenpakket 4D18
- Komponentpakke 4D18
- Zestaw naprawczy 4D18
- 4D18 alkatrészosomag
- Σετ μεμονωμένων εξαρτημάτων 4D18
- Комплект деталей 4D18
- 交換パーツセット 4D18
- 维修组件 4D18




- ▲ Mindestmenge
- Minimum order quantity
- Quantité minimum
- Quantità minima per l'ordine
- Cantidad mínima
- Quantidade mínima
- Minimikvantitet
- Mindste mængde
- Ilość minimalna
- Minimális mennyiség
- Ελάχιστη ποσότητα
- Минимальное количество
- 最低注文個数あり
- 最低起订量

- Einzelteile
- Single components
- Pièces détachées
- Singoli componenti
- Componentes
- Componentes
- Onderdelen
- Komponenter
- Podzespoły
- alkatrészek
- Μεμονωμένα εξαρτήματα
- Отдельные детали
- 個別販売
- 单个部件

①

INFORMATION

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INFORMATION

Last update: 2020-04-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document in a safe place.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

1 Single components (Fig. 1)

1.1 4D18 Single-component pack ●

Consisting of: 2 stops, 4 truss head screws, 4 two-hole nuts, 4 slotted bushings, 4 Belleville spring washers, 3 set screws, 2 lock rings, 2 hexagon nuts, 1 shock guard, 2 attachment nipples (short), 1 attachment nipple (long).

1.2 Minimum order quantity ▲

502R1=M5x16 Two-hole Nut (19), 501T1=M5x16 Truss Head Screw (18), 506G3=M8x12-V Set Screw (12), 506G3=M8x10 Set Screw (12a)

502S19=M6 Lock Nut (14), 4B82 Slotted Bushing (9), 4B86 Axle Bushing (8) 513T4=18.8x10.2 Belleville Spring Washer (7), 507S16=10x1 Lock Ring (13), 4Z43=H Stop (15), 507U16=6.4 Spherical Washer (10), 501T15=M6x40 Socket Head Cap Screw (11)

1.3 Single component ■

4G70 Lamination Anchor (1), 4X16 Setting Aid (23)

2 Description

2.1 Indications for use

The polycentric 3R46 Modular Knee Joint is to be used **exclusively** for the prosthetic fitting of the lower extremities.


2.2 Field of application

⚠ CAUTION**Excessive strain on the product**

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

Field of application according **to the Ottobock MOBIS Mobility System:**

Recommended for amputees with **mobility grade 3 and 4** (unrestricted outdoor walkers and unrestricted outdoor walkers with particularly high demands)

Approved for a patient weight of up to 125 kg / 275 lbs.

The hydraulic swing phase control allows to optimize the gait, in particular for dynamic patients.

The Ottobock MOBIS is not intended to be used as a guide to obtaining reimbursement for prosthetic components in the USA.

2.3 Function

The upper joint section (3) with lamination anchor and lower joint section (2) are connected to one another by axle brackets (4+5) and linkage bars (6). The cosmetic foam cover is only extended by a small amount due to the polycentric design. The lamination anchor (1) is attached to the residual limb bedding and constitutes the adjustable connection piece to the coupling core. The pivot point (instantaneous center of rotation) therewith changes its position in relation to the flexion. In an extended position, the pivot point lies at the height of the femoral condyles and moves downward and forward with increasing flexion.

The integrated compact hydraulic regulator (17) serves to generate movement resistances. These prevent a swing when flexing that is too wide or a stop when extending that is too hard. Contrary to mechanical gait regulators, the hydraulic unit adapts automatically to different walking speeds, thereby making for a harmonious gait. The movement resistance (damping resistance) can be regulated by adjusting the flow cross section, whereby flexion and extension can be set independently of each other (see Section 4).

2.4 Reuse and service life

CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

CAUTION

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the service life defined in this section is not exceeded.

As a basic principle, all Ottobock knee joints are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

2.5 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions
Temperature range for use: -10 °C – +60 °C
Use - Relative humidity: 0 % – 90 % (Non-condensing)
Storage/Transportation - Relative humidity: 100 %
Tyre pressure: 250 – 1100 mbar
Commercially available, solvent-free cleaning agents

Unallowable environmental conditions
Storage/Transportation: Mechanical vibrations, Impacts
Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine
Cleaning agents containing solvents

2.6 General safety instructions

CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

CAUTION

Overheating of the hydraulics due to excessive strain

Burns, injuries caused by falling due to changes in functionality and damage to prosthetic components

- ▶ Do not touch overheated prosthetic components.
- ▶ Reduce all activities if there are changes in functionality, in order to allow the overheated prosthetic components to cool down.
- ▶ If overheating or changes in functionality occur, have the prosthesis checked by authorised qualified personnel.

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

3 Technical data 3R46

Article number	3R46
Proximal connection	Lamination anchor
Distal connection	Adjustment pyramid
Knee flexion angle	110°
Weight	740 g
System height	99 mm
System height until alignment reference point	82 mm
System height from alignment reference point	17 mm
System height 90°	30 mm
Max. user weight	125 kg / 275 lbs
Mobility Grade	3 + 4

4 Handling

CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

CAUTION

Initial use of the prosthesis by the patient

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position of the axes.

The advantages of the 3R46 Knee Joint can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor.

To align the prosthesis please proceed in two steps:

- 1) First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly.
- 2) Then use 743L100 L.A.S.A.R. Posture for static alignment optimisation.

4.1.1 Bench alignment with alignment tool (the following steps refer to Fig. 7)

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

- 1 Position the middle of the foot **30 mm anterior** to the alignment reference line.
- 2 Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- 3 Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run through the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- 4 Connect the foot to the modular knee joint using a tube adapter.
- 5 Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- 6 Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.

If the residual limb flexion is not taken into account, the joint will be positioned too far to the front. This will lead to malfunction and premature wear.

4.1.2 Static alignment optimization using L.A.S.A.R. posture 743L100

(the following steps refer to Fig. 8)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- 1 To make the load line visible, the amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- 2 Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **35 mm anterior to the upper front knee axis** (see Fig. 8).
- 3 After step 2, dynamic optimisation can take place between parallel bars.

⚠ CAUTION

Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

When making extension damping adjustments, make sure that the leg still completely extends with each step. When adjusting the swing phase control, the influence of a cosmetic foam cover must be compensated for.

4.2 Combination possibilities

⚠ CAUTION

Failure to observe the manufacturer requirements regarding combination possibilities

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

INFORMATION

When fabricating a prosthesis for the lower limbs, all of the prosthetic components used have to meet the requirements in regard to the patient's body weight and activity level.

4.3 Adjustments and final assembly

4.3.1 Socket fabrication

4.3.1.1 Lamination prior to trial fitting

Take a piece of 623T3 Perlon Stockinette twice the length of the plaster model and pull half of it over the isolated socket. Twist the other half of the stockinette on the distal end and reflect it over the socket. In order to absorb strong forces in the area of the modular knee joint, the socket must be gradually reinforced with 616G13 Fibreglass Stockinette. To do this, pull over the first layer to 2/3 of the socket length, tie it off and then reflect it up to half of the socket length. Put 2 layers of 616G12 Carbon-Fibre Cloth on the distal area so as to provide an underlying base with a projecting brim of 3 cm around the 4G70 Adapter, which will be mounted later in the procedure. Pull on two layers of 623T3 Perlon Stockinette. The socket lamination is performed in two stages. In the first stage the socket is laminated up to 2/3 of its length with 617H19 Orthocryl Lamination Resin. In the second stage the proximal part of the socket is laminated with 617H17 Orthocryl Flexible Resin. Once the first casting has hardened, pull on 2 more layers of 623T3 Perlon Stockinette before laminating the proximal part of the socket with 617H17 Orthocryl Flexible Resin.

The lamination anchor serves to restrict flexion by stopping at the shock guard (20) (Fig. 2). This stop function must also be taken into consideration after the final lamination. To ensure a proper flexion stop and to prevent damage, the set screw (12), which may be too long for certain adjustment positions, should be replaced by a shorter one (12a).

4.3.1.2 Applying the lamination anchor

The correctly fitted lamination anchor must be prepared before it is glued to the socket. To do so, fill the coupling opening with 636K8 Plastaband (Fig. 3). Clamp the residual limb bedding and the joint into the alignment apparatus.

Remove any soft, flexible or porous materials from underneath the supporting surface of the lamination anchor. **Mix the putty exclusively from 617H21 Orthocryl Sealing Resin and 639A1 Talcum Powder. Glue on the lamination anchor.**

Secure with 627B2 Adhesive Tape for trial fitting. Test the stop function. If necessary, form the required contact surface with putty. If need be, glue a stop made of Pedilin onto the outer laminate.

4.3.1.3 Finalizing the socket

Once the trial fitting has been completed, fasten the lamination anchor to the socket with truss head screws (17) and two-hole nuts (18). Then laminate.. Then proceed with reinforcement as follows: Pull one layer of 623T3 Perlon Stockinette over the entire socket and tie circularly on the top. Apply two layers of 616G12 Carbon Fiber Cloth over the arms of the 4G70 Adapter. Then pull the remaining half of the 623T3 Perlon Stockinette back over the socket. The carbon fiber layers are thus sandwiched between the stockinette layers. Again reinforce the socket in increments with 616G13 Fiberglass Stockinette (see Section 4.3.1.1). Apply two layers of 623T3 Perlon Stockinette. Proceed with lamination as described in Section 4.3.1.1.

Deviating from the processing instructions or the recommended materials for the lamination of the lamination anchor can result in a loosening or breaking of the adapter.

Upon completion of assembly, tighten the set screws (12, 12a) with the 710D1 Torque Wrench. For trial fitting: Tighten the 506G3 Set Screws to 10 Nm. For completion: Secure the 506G3 Set Screws with 636K13 Loctite and tighten to 10 Nm.

4.3.2 Adjustment of the bearing

The axis pins of the linkage bars, encased in bushings, are connected to the upper and lower joint section. The type of bearing and screw connection allows to increase or decrease the friction and to remove play resulting from ordinary wear and tear. For this, attach lock nut (14) using enclosed 709S18 Socket Wrench (16). Then tighten socket head cap screw (11) with 5-mm Allen wrench (Fig. 4). It is of utmost importance that an adjustment to remove play take place after a breaking-in period of 2 to 4 weeks. At every maintenance check, the adjustment must be verified and the 4Z43=H Stops (15) replaced. To lubricate the axle bushings (8) and slotted bushings (9) use 633G6 Ottobock Special Lubricant.

4.3.2.1 Regulating the movement resistance

The flexion and extension resistance can be set independently of each other by sliding the pins in the transverse slots at the back and front side of the piston rod of the hydraulic regulator (17). To do so, use the enclosed setting aid (23).

In direction + = more resistance

In direction – = less resistance

4.3.2.2 Regulating the flexion resistance (Fig. 5)

Sliding the back pin (at the flexion side) to the **right** = decreases resistance = facilitates flexion. Sliding it to the **left** = increases resistance = impedes flexion.

4.3.2.3 Regulating the extension resistance (Fig. 6)

Flex joint so that the front pin is accessible.

Sliding pin to the **right** = **decreases** resistance = **facilitates** extension.

Sliding pin to the **left** = **increases** resistance = **impedes** extension.

4.4 Cosmetic foam cover

CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

Use Foam Cover 3R6, 6R6 and 3S107TF for Joint 3R46.

After affixing the cosmetic elements, verify proper functioning of the prosthesis.

To reduce friction and to eliminate noise, apply silicone spray directly onto the contact surfaces of the cosmetic foam cover.

To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

The elastic sleeve of the hydraulic regulator must be protected from mechanical damage.

The hydraulic regulator should not be disassembled as its position within the joint is fine-tuned. In case of failure, please send the entire joint to Ottobock Service.

4.5 Information for use

CAUTION

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

4.6 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
 - ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.
-
- ▶ Clean the product with a damp, soft cloth.
 - ▶ Dry the product with a soft cloth.
 - ▶ Allow to air dry in order to remove residual moisture.

4.7 Maintenance instructions

CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

Ottobock recommends readjusting the knee joint settings once the patient has gotten used to the prosthesis. The duration of this breaking-in period may differ from patient to patient. Check the knee joint for wear and proper function at least once a year and make adjustments if necessary. When doing so, pay special attention to the resistance to movement and to the development of unusual noise. A complete flexion and extension must be ensured.

We recommend carrying out regular safety checks once a year.

Proceed **with caution** and in **small increments** when adjusting the resistances, since there is the risk of falling!

5 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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