
1E95 Challenger

DE Gebrauchsanweisung	4
EN Instructions for use	10
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	23
ES Instrucciones de uso	29
PT Manual de utilização	35
NL Gebruiksaanwijzing	41
SV Bruksanvisning	48
DA Brugsanvisning	53
NO Bruksanvisning	59
FI Käyttöohje	65
PL Instrukcja użytkowania	72
HU Használati utasítás	78
CS Návod k použití	84
RO Instrucțiuni de utilizare	91
HR Upute za uporabu	97
SL Navodila za uporabo	103
SK Návod na použitie	109
BG Инструкция за употреба	115
TR Kullanma talimatı	121
EL Οδηγίες χρήσης	127
RU Руководство по применению	134
JA 取扱説明書	141
ZH 使用说明书	147
KO 사용 설명서	152

1

INFORMATION

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INFORMATION

Date of last update: 2020-07-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 1E95 Challenger is a prosthetic foot for fitness activities involving rapid changes of direction and speed (such as tennis, basketball or jogging) as well as for normal daily use. The contour of the spring is modelled on a sprinting foot. It also has a base spring with heel wedge and permits normal walking.

The prosthetic foot is worn directly in the shoe without a footshell.

1.2 Combination possibilities

INFORMATION

The prosthetic socket must be able to withstand the higher demands that arise when practising fitness activities.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

- Only use prosthetic components approved for the desired activity or a body weight of 150 kg.
- **Compatible Ottobock prosthetic knee joints:** Due to its specific kinematics, the prosthetic foot is only compatible with Genium and Genium X3.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthesis fittings.

2.2 Area of application

- The maximum approved body weight is specified in the technical data (see Page 16).

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60 °C
Storage and transport conditions: -10 °C to +60 °C, relative humidity: no restrictions
Moisture: fresh water, salt water, chlorinated water, sweat, urine
Dust, sand, highly hygroscopic particles (e.g. talcum)
Clean the product after it has come into contact with moisture/sand/dirt.

Unallowable environmental conditions
Acids

2.4 Service life

Prosthetic foot

An exact service life cannot be determined when the product is being used for fitness activities since a very wide range of loads occur depending on the activity and intensity.

Heel wedge, forefoot pad

The product is a wear part, which means it is subject to normal wear and tear.

3 Safety

3.1 Explanation of warning symbols

 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

3.2 General safety instructions

 CAUTION
Usage under extreme loads Risk of injury due to excessive loads during extreme activity
<ul style="list-style-type: none">▶ Although the product was developed for fitness activities, it may be damaged if exposed to extreme loads.▶ Inspect the product for damage before and after each activity.▶ If damage is apparent or in case of doubt, do not continue using the product.

 CAUTION
Use under unallowable environmental conditions Risk of injury due to damage to the product
<ul style="list-style-type: none">▶ Do not expose the product to unallowable environmental conditions.▶ If the product has been exposed to unallowable environmental conditions, check it for damage.▶ If damage is apparent or in case of doubt, do not continue using the product.

- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 11).

⚠ CAUTION

Reuse on other patients and improper maintenance

Fall due to loss of functionality or damage to the product

- ▶ Only use the product for a single patient.
- ▶ Maintain the product regularly to ensure a long service life (see Page 15).

⚠ CAUTION

Mechanical overload

Falling due to breakage of load-bearing components, impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the carbon spring are indications of loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	647G973
1	Prosthetic foot	–

Quantity	Designation	Reference number
1	Heel wedge kit	2F95=*
1	Forefoot pad	2Z95=*
1	Shoehorn	2C101

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

5.1 Information on use

CAUTION

Using the prosthetic foot without a shoe

Falling due to lack of surface grip

- ▶ Always use the prosthetic foot with a shoe.

INFORMATION

The prosthetic foot has no footshell

The chosen shoe is pulled directly over the prosthetic foot. The forefoot pad on the prosthetic foot prevents sliding and cushions the shoe.

5.1.1 Shoe selection

INFORMATION

The spring of the prosthetic foot has a posterior projection. This means the prosthetic foot is not suitable for all shoes. If the upper of the shoe is too high, the projection of the prosthetic foot spring pushes on the upper. This can damage the shoe.

- ▶ Select shoes where the posterior projection of the prosthetic foot spring does not push on the upper of the shoe.
- ▶ **When walking down stairs with the prosthetic foot, note the projecting spring.**

5.1.2 Use in water

The prosthetic foot may come into contact with water, for example, when jogging on the beach or when swimming.

- 1) **NOTICE! Always use the prosthetic foot with a shoe when in water. Otherwise, the heel wedge may be lost.**
Before water contact: Check the firm fit of the shoe.
- 2) **After water contact:** Clean the prosthetic foot (see Page 15).

5.2 Alignment

INFORMATION

There is a plastic adapter cover on the proximal connection of the prosthetic foot. This protects the connecting section from scratches during the alignment and trial fitting.

- ▶ Remove the adapter cover before the patient leaves the workshop/fitting area.

5.2.1 Selection of stiffness

The stiffness of the prosthetic foot is selected depending on the body weight.

- The maximum approved body weight is specified in the technical data (see Page 16).

5.2.2 Bench Alignment

INFORMATION

Only use the prosthetic foot with the heel wedge.

INFORMATION

- For determining the overall system height of all prosthetic components for bench alignment, add 5-10 mm to the socket-floor measurement to compensate for the spring effect of the prosthetic foot. Use the larger value when the patient's weight is at the upper limit of the maximum permissible patient weight.

TT bench alignment

Bench alignment process

Recommended materials: 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)

Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:

Sagittal plane

①	Heel height: effective heel height (shoe heel height - sole thickness of forefoot) +10 mm
②	Exterior foot rotation: approx. 5°
③	Sagittal alignment: alignment reference line on a-p markings of the prosthetic foot (see fig. 1)
④	Connect the prosthetic foot and prosthetic socket using the chosen adapters. Observe the adapter instructions for use during adaptation and installation.
⑤	Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: individual residual limb flexion +5°

Frontal plane

⑥	Prosthetic foot alignment reference line: set exterior rotation according to the patient's requirements Alignment reference line of prosthetic socket: along the lateral patella edge
⑦	Observe the individual abduction or adduction position of the residual limb.

TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

5.2.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.2.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.

- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

5.2.4.1 Optimising the heel characteristics

INFORMATION

Only use the prosthetic foot with the heel wedge.

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted by replacing the heel wedge. The prosthetic foot is delivered with 2 heel wedges.

Heel wedge degrees of hardness: The heel wedges are marked with numbers. The smaller the number, the softer the heel wedge.

- 1) Pull the prosthetic foot apart slightly and pull out the existing heel wedge.
- 2) Align the other heel wedge so the Ottobock lettering is legible and the tip points to anterior.
- 3) Insert the heel wedge in the prosthetic foot (see fig. 2, see fig. 3, see fig. 4).

6 Cleaning

> **Recommended cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)

- 1) Clean the product with clear fresh water and a pH neutral soap.
- 2) Rinse the soap away with clear fresh water.
- 3) Dry the product with a soft cloth.
- 4) Allow to air dry in order to remove residual moisture.

7 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

► Observe the following maintenance instructions.

- The prosthetic components should be inspected after the first 30 days of use.
- Inspect the entire prosthesis for wear during normal consultations.
- Conduct annual safety inspections.

7.1 Mounting/changing the forefoot pad

The prosthetic foot comes with a forefoot pad and a hook-and-loop strip. Use the hook-and-loop strip to fasten the forefoot pad (see title image).

- 1) Remove the forefoot pad from the prosthetic foot.
- 2) Place a new forefoot pad on the prosthetic foot.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical data

Size [cm]	23/24	25/26	27/28	29/30
Heel height [mm]	10			
System height [mm]	156			
Weight [g]	770	785	800	820
Max. body weight [kg] for prosthetic foot stiffness 2	70			-
Max. body weight [kg] for prosthetic foot stiffness 3	-	90		
Max. body weight [kg] for prosthetic foot stiffness 4	-	110		

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