
8E51

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	12
FR Instructions d'utilisation (Personnel spécialisé)	21
IT Istruzioni per l'uso (Personale tecnico specializzato)	32
ES Instrucciones de uso (Personal técnico especializado)	41
PT Manual de utilização (Pessoal técnico)	51
NL Gebruiksaanwijzing (Vakmensen)	61
SV Bruksanvisning (Fackpersonal)	70
FI Käyttöohje (Ammattihenkilöstö)	79
SK Návod na používanie (Odborný personál)	89
EL Οδηγίες χρήσης (Τεχνικό προσωπικό)	98
RU Руководство по применению (Квалифицированный персонал)	109
JA 取扱説明書 (有資格担当者)	119

INFORMATION

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INFORMATION

Date of last update: 2017-01-18

- ▶ Please read this document carefully before using the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

The Electric Hand 2000 is referred to simply as the product/prosthesis/terminal device below. These instructions for use provide you with important information on the use, adaptation and handling of the product. Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Product description

2.1 Function

The product is a myoelectrically controlled terminal device of the MyoBock prosthesis system.

On the terminal device, the thumb and fingers rotate around the same axis and carry out a spatial movement in the process. Therefore, objects can be gripped without unnatural compensating movements of the forearm and upper arm.

Size 5 is driven by one motor.

Sizes 5 ½, 6 and 6 ½ are driven by two motors with a multi-stage overriding gear mechanism. The first motor opens and closes the hand at high speed and low force. The second motor is activated in order to develop the gripping force required to pick up objects.

2.2 Combination possibilities

This product can be combined with the following Ottobock components:

Use of the 9E369 or 9E370 4 in 1 Controller

- 9E369 or 9E370 4 in 1 Controller
- 13E182 coding plug set or 13E184=* individual coding plug
- 757B13 replaceable battery
- 757L13 battery charger

Use of the 9E420* 7 in 1 Controller

- 9E420=L or 9E420=R 7 in 1 Controller
- 757B35=0 or 757B35=1 or 757B35=3 MyoEnergy Integral
- 757L35 MyoCharge Integral

Wrist joints

- 10V51=1 or 10V51=2 MyolinoWrist 2000

Elbow

- 12K19=40 elbow component
- 12K12 MovolinoArm Friction

Lamination rings

- 10S16=34 or 10S16=38 lamination ring
- 10S18=40 lamination ring
- 11D12=34 lamination ring

3 Application

3.1 Indications for use

The Electric Hand 2000 is intended **exclusively** for exoprosthetic fittings of the upper limbs.

3.2 Use/Field of Application

The Electric Hand 2000 can be used for unilateral or bilateral amputations, and in case of dysmelia.

3.3 Environmental conditions

The prosthesis system was developed for everyday activities and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the prosthesis system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see Page 20).

3.4 Qualification

Fitting a patient with the product may only be carried out by orthopaedic technicians who have been authorised by Ottobock after completion of a corresponding training course.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

 CAUTION
<p>The heading describes the source and/or the type of hazard</p> <p>The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:</p> <ul style="list-style-type: none">> E.g.: Consequence 1 in case of failure to observe the hazard> E.g.: Consequence 2 in case of failure to observe the hazard▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

 WARNING
<p>Non-observance of safety instructions</p> <p>Injury due to damage or unexpected product behaviour.</p> <ul style="list-style-type: none">▶ Follow the safety instructions and the precautions specified in this accompanying document.
 WARNING
<p>Operating a vehicle with the prosthesis</p> <p>Accidents due to unexpected actions of the prosthesis.</p> <ul style="list-style-type: none">▶ The prosthesis should not be used for the operation of motor vehicles and heavy equipment (e.g. construction machines).
 WARNING
<p>Operating machines with the prosthesis</p> <p>Injury due to unexpected actions of the prosthesis.</p> <ul style="list-style-type: none">▶ The prosthesis should not be used for the operation of industrial machines or motor-driven equipment.
 WARNING
<p>Operating the prosthetic system near active implanted systems</p> <p>Interference with active implantable systems (e.g., pacemaker, defibrillator, etc.) due to electromagnetic interference of the prosthetic system.</p> <ul style="list-style-type: none">▶ When operating the prosthetic system in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.

- ▶ Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.

⚠ WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

⚠ WARNING

Charging the prosthesis without taking it off

Risk of electric shock due to defects in the power supply unit or in the battery charger.

- ▶ For safety reasons, remove the prosthesis prior to charging it.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Signs of wear to system components

Injury due to faulty control or malfunction of the product.

- ▶ In the interest of the patient's safety and in order to maintain operating reliability, the product should be serviced at regular intervals.

4.4 Information on Alignment/Adjustment

⚠ CAUTION

Operator errors during the adjustment process with the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During this product training course you will receive a password giving you access to the adjustment software. Additional product training courses may become necessary to qualify for software updates.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

⚠ CAUTION

Insufficient skin contact of the electrodes

Injury due to faulty control or malfunction of the product.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible.

- ▶ In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary.
- ▶ If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.

⚠ CAUTION

Incorrect electrode settings

Injury due to unexpected product behaviour.

- ▶ Set the electrode sensitivity as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

⚠ CAUTION

Incorrect electrode settings due to muscle fatigue

Injury due to faulty control or malfunction of the product.

- ▶ Allow the patient to rest during the adjustment of the electrodes.

⚠ CAUTION

Establishing a connection between the skin and metal components of the prosthesis with the use of carbon fibre

Injury due to unexpected product behaviour as a result of interference conducted by carbon fibres.

- ▶ During fabrication, make sure that no connection is formed between the skin and metal prosthesis components by carbon fibres.

INFORMATION

Do not use silicone spray when putting on the prosthetic glove. This impairs the secure fit of the glove.

Observe the instructions for use provided for the prosthetic glove and the process for putting it on and taking it off described in the same.

4.5 Information on Proximity to Certain Areas

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the prosthetic system.

⚠ CAUTION

Distance to HF communication devices is too small (e.g. mobile phones, Bluetooth devices, WiFi devices)

Injury due to unexpected behaviour of the prosthesis system caused by interference with internal data communication.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
 - Mobile phone GSM 850 / GSM 900: 0.99 m
 - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
 - DECT cordless phones incl. base station: 0.35 m
 - WiFi (routers, access points ...): 0.22 m
 - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

⚠ CAUTION

Wearing in extreme ambient temperatures

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Avoid wearing the prosthesis in areas with temperatures outside of the permissible range (see Page 20).

4.6 Information on Use

⚠ CAUTION

Using the prosthesis with pointed or sharp objects (e.g. knives in the kitchen)

Injury due to unintentional movements.

- ▶ Use extra caution when using the prosthesis for handling pointed or sharp objects.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Mechanical stress on the product

Injury due to faulty control or malfunction of the prosthetic system.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

⚠ CAUTION

Penetration of dirt and humidity into the prosthetic components

Injury due to faulty control or malfunction of the prosthesis system.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the prosthesis system or prosthetic components (e.g. the terminal device).

NOTICE

Improper care of the prosthesis system

Damage to the prosthesis system due to the use of incorrect cleaning agents.

- ▶ Clean the prosthesis system with a damp cloth and mild soap only (e.g. 453H10=1 Ottobock DermaClean).
- ▶ Use only the following products for cleaning/disinfecting the inner socket:
Cleaning: 453H10=1 Ottobock DermaClean
Disinfection: 453H16 Ottobock DermaDesinfect

4.7 Information on the Power Supply/Battery Charging

⚠ CAUTION

Charging the prosthesis system with soiled or damaged contacts

Injury due to unexpected behaviour of the prosthesis system caused by insufficient charging.

- ▶ Make sure that the contacts are always clean and free of grease.
- ▶ Clean the electrical contacts of the charging plug and charging receptacle regularly using cotton swabs and a mild soap solution.
- ▶ Take care to avoid damaging the contacts with pointed or sharp objects.

NOTICE

Charging the prosthesis system with an incorrect power supply/battery charger

Damage to the prosthesis system due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for the prosthesis system by Ottobock (see instructions for use and catalogues).

NOTICE

Contact of the charging plug with magnetic data carriers

Wiping of the data carrier.

- ▶ Do not place the charging plug on credit cards, diskettes, audio or video cassettes.

INFORMATION

The prosthesis system cannot be used as long as the charging plug is connected to the charging receptacle. The prosthesis system is deactivated for the duration of charging.

4.8 Notes on using an elbow joint (e.g. 12K19=40, 12K12)

⚠ CAUTION

Risk of pinching in elbow joint flexion region

Injury due to pinching of body parts.

- ▶ Ensure that fingers and other body parts are not in this area when bending the elbow joint.

⚠ CAUTION

Manual unlocking of elbow lock under load

Injury by release of elbow lock under load.

- ▶ Particular caution should be exercised when unlocking the elbow lock while lifting heavy loads.
- ▶ Be careful when unlocking the lock under such conditions due to the possibility of injury.

⚠ CAUTION

Penetration of dirt and moisture in the prosthetic arm

Injury due to unexpected behaviour of the prosthesis system or malfunction.

- ▶ Do not let foreign particles or liquids get into the prosthetic arm.
- ▶ Do not expose the prosthetic arm, especially the elbow joint, to dripping or splashing water.
- ▶ Wear the prosthetic arm and especially the elbow joint under suitably resistant clothing in the rain.

NOTICE

Coating, gluing or painting the prosthesis

Damage or fracture due to chemical processes.

- ▶ The prosthesis must not be coated, glued or painted.

5 Scope of delivery

- 1 pc. Electric Hand 2000
- 1 pc. 647H58 instructions for use (qualified personnel)

6 Preparation for use

6.1 Settings

The control settings depend on the controller being used. For more information, see the instructions for use of the controller in question.

6.2 Assembly

6.2.1 Installing the O-rings

- 1) Be sure to install the outer 623F3=26x1.0 O-ring (1).
- 2) Insert the number of inner O-rings (2) according to the desired friction level (see fig. 1).

6.2.2 Installing the controller

See the instructions for use of the controller being used.

6.2.3 Installing the hand

- 1) Put on the hand (see fig. 2).
- 2) Firmly tighten the wrist nut against the O-ring.

7 Cleaning

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1 Ottobock DermaClean).
- 1) Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.
Ensure that no liquid penetrates into the system component(s).
 - 2) Dry the product with a lint-free cloth and allow it to air dry fully.

8 Maintenance

Since all moving mechanical parts are subject to wear, regular service inspections are recommended. These service inspections include an inspection of the product and, when necessary, the replacement of wear parts.

9 Legal information

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9.3 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

Ambient conditions	
Storage and transport in original packaging	$\pm 0\text{ }^{\circ}\text{C}/+32\text{ }^{\circ}\text{F}$ to $+60\text{ }^{\circ}\text{C}/+140\text{ }^{\circ}\text{F}$ max. 93 % relative humidity, non-condensing
Storage and transport without packaging	$\pm 0\text{ }^{\circ}\text{C}/+32\text{ }^{\circ}\text{F}$ to $+60\text{ }^{\circ}\text{C}/+140\text{ }^{\circ}\text{F}$
Operation	$\pm 0\text{ }^{\circ}\text{C}/+32\text{ }^{\circ}\text{F}$ to $+60\text{ }^{\circ}\text{C}/+140\text{ }^{\circ}\text{F}$

Reference number	8E51			
Size	5	5½	6	6½
Operating voltage [V]	4.8/7.4			
Average current consumption [mA]	Approx. 200			
Electric shut-off "closed" at [mA]	Approx. 400			
Electric shut-off "open" at [mA]	Approx. 400	Approx. 250		
Opening width [mm]	28	37	52	58
Minimum possible gripping force [N]	8	25	35	
Weight without prosthetic glove [g]	86	115	125	130
Service life [years]	5			

11 Appendices

11.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

SN.YYYY.WW.NNN Serial number

11.2 Troubleshooting

Event	Cause/required action
Lack of gripping force	Insufficient contact between hand and switch unit. <ul style="list-style-type: none">Firmly screw the wrist nut onto the outer O-ring.
Excessive friction in pronation and supination direction	Excessive frictional forces on the O-rings. <ul style="list-style-type: none">Take out O-ring(s).
Insufficient friction in pronation and supination direction	Insufficient frictional forces on the O-rings, or no O-rings installed. <ul style="list-style-type: none">Insert O-ring(s).
Hand operates too slowly in open direction	Second motor in operation because the prosthetic glove is too stiff. <ul style="list-style-type: none">Warm up the prosthetic glove (max. 60°C/140°F), then open the hand wide. Allow to cool in open position. Battery charge level low <ul style="list-style-type: none">Charge battery
Hand operates too slowly in close direction	Second motor in operation because the prosthetic glove is too stiff. <ul style="list-style-type: none">Close hand and warm prosthetic glove from the outside (max. 60°C/140°F). Battery charge level low <ul style="list-style-type: none">Charge battery
Finger cover damaged	Cut open and pull off finger cover. Install new finger cover.

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Otto Bock Healthcare Products GmbH