
9X37

DE Gebrauchsanweisung	3
EN Instructions for Use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	16
ES Instrucciones de uso	20
PT Manual de utilização	24
NL Gebruiksaanwijzing	28
SV Bruksanvisning	32
DA Brugsanvisning	36
NO Bruksanvisning	40
PL Instrukcja użytkowania	44
CS Návod k použití	48
TR Kullanma talimatı	52
EL Οδηγίες χρήσης	56

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English

Date of the last update: 2014-01-24

- Please read this document carefully.
- Follow the safety instructions and the precautions specified in this accompanying document.
- If assistance is required during the start-up, use or maintenance of the product, as well as in the event of unexpected operating behaviour or circumstances, contact the manufacturer (see manufacturer's address on the back of this accompanying document).

1 Intended Use

The Pressure Switch 9X37 is to be used **solely** in exoprosthetic fittings of the upper extremities for myoelectric control of arm prostheses.

The Pressure Switch 9X37 can be used for both right or left arm amputees as well as bilateral amputees. It enables control of a System Electric Hand 8E38=*, a System Electric Greifer 8E33=* or a Wrist Rotator 10S17 in **combination** with a MyoRotronic 13E205. The fitting of a patient with the Pressure Switch 9X37 may only be carried out by a prosthetist.

2 Safety Instructions

Failure to follow the safety instructions given below can lead to damage to or malfunction of the product. Follow the safety instructions and the precautions specified in this accompanying document.

- When mounted into a provisional socket, the pressure switch may only be worn by an amputee within an orthopedic workshop or as an inpatient in a rehabilitation facility.

- When installing the connection cable, make sure there are no tight angles or kinks.
- The pressure switch should not be incorporated into bathing prostheses.
- Use the Patient Information section of this pamphlet (Section 3) to inform your patient about the function and use of the pressure switch.

3 Patient Information

- Handle the pressure switch with care and check the casing and pivoting mount for damage (e.g. cracks or breakage) before every use. Should the pressure switch be damaged or not function properly, please contact your prosthetist immediately.
- Driving a vehicle or operating machinery is strictly prohibited while wearing a provisional socket fitted with a pressure switch.
- Do not allow debris or liquids to get into the pressure switch.
- The pressure switch should not be subjected to intense smoke, dust, mechanical vibrations, shocks or high temperatures.
- The pressure switch as well as damaged components may only be opened or repaired by the certified Ottobock Myo-Service.
- Regularly clean the pressure switch with a cloth **slightly** dampened with soap solution. Do not allow moisture to enter into the pressure switch or into the contacts!

4 Scope of Delivery

1 pc.	Pressure Switch	9X37
2 pcs.	Oval Countersunk Head Screw	501S75=M2×8
2 pcs.	Oval Countersunk Head Screw	501S75=M2×5
1 pc.	Instructions for Use	647H71

5 Description

The pressure switch is used to control the prosthetic components with digital control or to switch between the prosthetic components for which the patient cannot provide his/her own myo-signal. It can be used in combination with MyoBock electrodes, linear transducer and 4-step transducer. Should the pressure switch be mounted in the provisional socket, the patient can practice the control or switching as early as the first stage of the rehabilitation. This will ensure a smooth transition to the use of the definitive prosthesis.

6 Function

The Pressure Switch 9X37 features 2 control points that are arranged and controlled one after another and independently. Each control point is assigned to a function for control of the prosthetic components. The control point 1 should always be assigned to the 'close hand' function.

Independent setting of pressure required for activating both control points can be performed by turning the adjustment screws on the reverse side (Fig. 1). The right-hand adjustment screw (Fig. 1, Item 1) is used for adjustment of the control point 1. The left-hand adjustment screw (Fig.1, Item 2) is used for adjustment of the control point 2.

Note: The adjustment screws should be screwed in and out with a fine touch. When turning the adjustment screw in counter-clockwise direction, please make sure that the switch is not locked and can be operated.

By pressing the switch, it will achieve control point 1 (S1) and interconnect the control signal A.

By pressing the switch more deeply, it will achieve control point 2 (S2) and interconnect the control signal B. Now the control signal B is solely active.

7 Assembly

Shorten the enclosed attachment screws for mounting the pressure switch to the socket so as not to exceed the screw depth 2.3 mm in the switch compartment.

8 Connecting the Connection Cable

The pressure switch is delivered without Connecting Cable 13E99=1200. The latter must be ordered separately. Apply abundantly Silicone Grease 633F11 before connecting the contacts of the connecting cable. The grease will protect the connector from humidity.

9 Technical Data

Service life of pressure switch 5 years

10 Disposal

	<p>These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.</p>
--	---

11 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Ottobock technicians.

12 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

13 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

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