

Evidence Essentials

Microprocessor-Controlled Hydraulic Ankle-Foot

	Mobility need or deficit of the patient	Evidence for benefits of the Meridium or other hydraulic or MP ankles compared to standard feet
Safety	Patient trips and falls repeatedly	<ul style="list-style-type: none"> - Significant improvement in minimal toe clearance (De Asha et al., 2014; Rosenblatt et al., 2014) - Significant increase in minimal toe clearance during slope ascent (Lamers et al., 2019) - Reduction of the risk of tripping over an unseen obstacle of 5 mm (1/5 in) height from 1 in 166 steps to 1 in 3,169 steps (Rosenblatt et al., 2014) - Patient reported improved perceived safety and balance as well as 35% fewer stumbles and 23% fewer falls with Meridium (Hahn et al., 2018)
Mobility	Patient feels limited or restricted in his/her overall mobility	<ul style="list-style-type: none"> - Cross-sectional study found that users of MP controlled ankles had significantly higher overall mobility compared to users of other types of prosthetic feet (Wurdeman et al., 2019) - Significant increase in self-selected walking speed and walking distance in K2 patients (Barnett et al., 2018) - Significant improvement in overall patient-reported prosthetic function in K2 and K3 patients as measured with the PEQ (Kaluf et al., 2020; Moore et al., 2017)
Mobility	Patient has difficulty negotiating uneven terrain	<ul style="list-style-type: none"> - Significantly improved gait symmetry and reduction in kinematic and kinetic gait deviations and, thus, unloading of the locomotor system (Bai et al., 2017; Schmalz et al., 2019) - Patients perceived walking on uneven terrain as safer and easier (Bai et al., 2017; Hahn et al., 2018)

Evidence Essentials

Microprocessor-Controlled Hydraulic Ankle-Foot

Mobility	Patient feels limited or restricted in ambulating on level ground	<ul style="list-style-type: none"> - Significant reduction in braking forces (“dead spot phenomenon”) in level walking with consecutive significant increase in self-selected walking speed (De Asha et al., 2013A and 2014) - Significant reduction in kinetic compensations of the sound limb (De Asha et al., 2013B) - Significantly improved gait symmetry and reduction in kinematic and kinetic gait deviations in K2 and K3 as well as bilateral patients (Moore et al., 2016) - Significantly reduced metabolic energy consumption (Askew et al., 2019; Delussu et al., 2013)
Mobility	Patient has difficulty ascending and descending stairs	<ul style="list-style-type: none"> - Improved gait symmetry and reduction in kinematic and kinetic gait deviations during stair ascent and descent and, thus, unloading of the locomotor system (Alimusaj et al., 2009)
Residual limb health	Patient suffers from residual limb pressure sores or ulcers due to frequent ambulation on non-level surfaces	<ul style="list-style-type: none"> - Significant reduction in interface pressures between the residual limb and the socket when walking on level and uneven terrain, slopes and stairs (Moore et al., 2018; Portnoy et al., 2012; Wolf et al., 2009) - Significantly improved socket comfort when walking on non-level terrains (Kaluf et al., 2020)

Evidence Essentials

Microprocessor-Controlled Hydraulic Ankle-Foot

Mobility	Patient has difficulty standing on and/or negotiating slopes and hills	<ul style="list-style-type: none"> - Significant increase in minimal toe clearance during slope ascent (Lamers et al., 2019) - Significant increase in self-selected walking speed during slope ascent (Askew et al., 2019) - Significantly slowed shank rotational velocity and, thus, improved control over walking speed during slope descent (Struchkov et al., 2016) - Significantly improved gait symmetry reduction in kinematic and kinetic gait deviations and, thus, unloading of the locomotor system during slope descent (Darter et al., 2014; Fradet et al., 2010; McGrath et al., 2018) - Patients perceived slope ambulation and standing on slopes as safer and easier (Darter et al., 2014; Fradet et al., 2010; Hahn et al., 2018) - Significantly improved symmetry and reduction in kinematic and kinetic gait deviations and, thus, unloading of the locomotor system during standing on slopes (Ernst et al., 2017; McGrath et al., 2019; Thomas-Pohl et al., 2019) - Significantly reduced metabolic energy consumption for slope ascent (Askew et al., 2019)
-----------------	--	---

References

Alimusaj M, Fradet L, Braatz F, Gerner HJ, Wolf SJ. Kinematics and kinetics with an adaptive ankle foot system during stair ambulation of transtibial amputees. *Gait Posture* 2009; 30: 356-363. [Download](#)

Askew G, McFarlane LA, Minetti AE, Buckley JG. Energy cost of ambulation in trans-tibial amputees using a dynamic-response foot with hydraulic versus rigid `ankle`: insights from body centre-of-mass dynamics. *J Neuroeng Rehabil* 2019 Mar 14;16(1):39. [Download](#)

Bai X, Ewins D, Crocombe AD, Wei X. Kinematic and biomimetic assessment of a hydraulic ankle/foot in level ground and camber walking. *PLoS One* 2017; DOI:10.1371/journal.pone.0180836 July 13, 2017. [Download](#)

Evidence Essentials

Microprocessor-Controlled Hydraulic Ankle-Foot

Barnett C, Brown O, Bisele M, Brown M, De Asha A, Strutzenberger G. Individuals with unilateral trans-tibial amputation and lower activity levels walk more quickly when using a hydraulically articulating versus rigidly-attached prosthetic ankle-foot device: A case series. *J Prosthet Orthot* 2018;30(3):158-164. [Download](#)

Darter BJ, Wilken JM. Energetic consequences of using a prosthesis with adaptive ankle motion during slope walking in persons with a transtibial amputation. *Prosthet Orthot Int* 2014; 38: 5 published online 22 March 2013. [Download](#)

De Asha AR, Johnson L, Munjal R, Kulkarni J, Buckley JG. Attenuation of centre-of-pressure trajectory fluctuations under the prosthetic foot when using an articulating hydraulic ankle attachment compared to fixed attachment. *Clin Biomech* 2013A;28: 218-224. [Download](#)

De Asha AR, Munjal R, Kulkarni J, Buckley JG. Walking speed related joint kinetic alterations in trans-tibial amputees: impact of hydraulic ‘ankle’ damping. *J Neuroeng Rehabil* 2013B; 10: 107-121. [Download](#)

De Asha AR, Munjal R, Kulkarni J, Buckley JG. Impact on the biomechanics of overground gait using an ‘Echelon’ hydraulic ankle-foot device in unilateral trans-tibial and trans-femoral amputees. *Clin Biomech* 2014; 29: 728-734. [Download](#)

Delussu AS, Brunelli S, Paradisi F, Iosa M, Pellegrini R, Zenardi D, Traballes M. Assessment of the effects of carbon fiber and bionic foot during overground and treadmill walking in transtibial amputees. *Gait Posture* 2013; 38: 876-882. [Download](#)

Ernst M, Altenburg B, Bellmann M, Schmalz T. Standing on slopes – how current microprocessor-controlled prosthetic feet support transtibial and transfemoral amputees in everyday tasks. *J NeuroEngin Rehabil* 2017 Nov 16;14(1):117. [Download](#)

Fradet L, Alimusaj M, Braatz F, Wolf SI. Biomechanical analysis of ramp ambulation of transtibial amputees with an adaptive ankle system. *Gait Posture* 2010; 32: 191-198. [Download](#)

Hahn A, Sreckovic I, Reiter S, Mileusnic M. First results concerning the safety, walking and satisfaction with an innovative, microprocessor-controlled four axes prosthetic foot. *Prosthet Orthot Int* 2018 Jun;42(3):350-356. [Download](#)

Johnson L, De Asha AR, Munjal R, Kulkarni J, Buckley JG. Toe clearance when walking in people with unilateral amputation: Effects of passive hydraulic ankle. *J Rehabil Res Dev* 2014;51(3):429-438. [Download](#)

Kaluf B, Duncan A, Bridges W. Comparative effectiveness of microprocessor-controlled and carbon-fiber energy-storing-and-returning prosthetic feet in persons with unilateral transtibial amputations: Patient-reported outcome measures. *J Prosthet Orthot* 2020; in print. [Download](#)

Lamers EP, Eveld ME, Zelik KE. Subject-specific responses to an adaptive ankle prosthesis during incline walking. *J Biomech* 2019 Oct 11;95:109273. doi: 10.1016/j.jbiomech.2019.07.017. Epub 2019 Jul 26. [Download](#)

McGrath M, Laszczak P, Zahedi S, Moser D. The influence of a microprocessor controlled hydraulic ankle on the kinetic symmetry of trans-tibial amputees during ramp walking: a case series. *J Rehabil Assist Technol Eng* 2018 Oct 8;5:2055668318790650. doi: 10.1177/2055668318790650. eCollection 2018 Jan-Dec. [Download](#)

McGrath M, Davies KC, Laszczak P, Rek B, McCarthy J, Zahedi S, Moser D. The influence of hydraulic ankles and microprocessor control on the biomechanics of trans-tibial amputees during quiet standing on a 5° ramp. *Can Prosthet Orthot J* 2019;2(2):2. [Download](#)

Evidence Essentials

Microprocessor-Controlled Hydraulic Ankle-Foot

Moore R. Effect on stance phase timing asymmetry in individuals with amputation using hydraulic ankle units. *J Prosthet Orthot* 2016;28(1):44-48. [Download](#)

Moore R. Patient evaluation of a novel prosthetic foot with hydraulic ankle aimed at persons with amputation with lower activity levels. *J Prosthet Orthot* 2017;29(1):44-47. [Download](#)

Moore R. Effect of a prosthetic foot with a hydraulic ankle unit on the contralateral foot peak plantar pressures in individuals with unilateral amputation. *J Prosthet Orthot* 2018;30(3):165-170. [Download](#)

Portnoy S, Kristal A, Gefen A, Siev-Ner I. Outdoor dynamic subject-specific evaluation of internal stresses in the residual limb: Hydraulic energy-stored prosthetic foot compared to conventional energy-stored feet. *Gait Posture* 2012; 35:121-125. [Download](#)

Rosenblatt NJ, Bauer A, Rotter D, Grabiner MD. Active dorsiflexing prostheses may reduce trip-related fall risk in people with transtibial amputation *J Rehabil. Res Dev* 2014; 51(8): 1229-1242. [Download](#)

Schmalz T, Altenburg B, Ernst M, Bellmann M, Rosenbaum D. Lower limb amputee gait characteristics on a specifically designed ramp: Preliminary results of a biomechanical comparison of two prosthetic foot concepts. *Gait Posture* 2019 Feb;68:161-167. doi: 10.1016/j.gaitpost.2018.11.017. Epub 2018 Nov 16. [Download](#)

Struchkov V, Buckley JG. Biomechanics of ramp descent in unilateral trans-tibial amputees: Comparison of a microprocessor controlled foot with conventional ankle-foot mechanisms. *Clin Biomech (Bristol, Avon)* 2016 Feb;32:164-70. [Download](#)

Thomas-Pohl M, Villa C, Davot J, Bonnet X, Facione J, Lapeyre E, Bascou J, Pillet H. Microprocessor prosthetic ankles: comparative evaluation of people with transtibial traumatic amputation during standing on level ground and slope. *Disabil Rehabil Assist Technol* 2019 Sep 19:1-10. doi: 10.1080/17483107.2019.1629112. [Epub ahead of print] [Download](#)

Wolf SI, Alimusaj M, Fradet L, Siegel J, Braatz F. Pressure characteristics at the stump/socket interface in transtibial amputees using an adaptive prosthetic foot. *Gait Posture* 2009; 24: 860-865. [Download](#)

Wurdeman SR, Stevens PM, Campbell JH. Mobility analysis of Amputees (MAAT 5): Impact of five common prosthetic ankle-foot categories for individuals with diabetic/dysvascular amputation. *J Rehabil Assist Technol Eng* 2019 Feb 13;6:2055668318820784. doi: 10.1177/2055668318820784. eCollection 2019 Jan-Dec. [Download](#)