## 1C30-1 Trias

|           |                                  |     |
|-----------|----------------------------------|-----|
| <b>DE</b> | Gebrauchsanweisung .....         | 3   |
| <b>EN</b> | Instructions for use .....       | 10  |
| <b>FR</b> | Instructions d'utilisation ..... | 17  |
| <b>IT</b> | Istruzioni per l'uso .....       | 25  |
| <b>ES</b> | Instrucciones de uso .....       | 33  |
| <b>PT</b> | Manual de utilização .....       | 41  |
| <b>NL</b> | Gebruiksaanwijzing .....         | 48  |
| <b>SV</b> | Bruksanvisning .....             | 56  |
| <b>DA</b> | Brugsanvisning .....             | 63  |
| <b>NO</b> | Bruksanvisning .....             | 70  |
| <b>FI</b> | Käyttöohje .....                 | 77  |
| <b>PL</b> | Instrukcja użytkowania .....     | 85  |
| <b>HU</b> | Használati utasítás .....        | 93  |
| <b>CS</b> | Návod k použití .....            | 100 |
| <b>RO</b> | Instrucțiuni de utilizare .....  | 107 |
| <b>HR</b> | Upute za uporabu .....           | 115 |
| <b>SL</b> | Navodila za uporabo .....        | 122 |
| <b>SK</b> | Návod na používanie .....        | 129 |
| <b>BG</b> | Инструкция за употреба .....     | 137 |
| <b>TR</b> | Kullanma talimatı .....          | 145 |
| <b>EL</b> | Οδηγίες χρήσης .....             | 152 |
| <b>RU</b> | Руководство по применению .....  | 160 |
| <b>JA</b> | 取扱説明書 .....                      | 168 |
| <b>ZH</b> | 使用说明书 .....                      | 175 |
| <b>KO</b> | 사용 설명서 .....                     | 181 |

1

## INFORMATION

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### INFORMATION

Date of last update: 2021-01-29

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 1C30-1 Trias prosthetic foot is suitable for walking on various surfaces and at variable speeds.

The functional properties of the prosthetic foot are achieved through the carbon spring elements. The prosthetic foot permits perceptible plantar flexion at heel strike and a natural rollover movement. Stored energy is released by the spring elements.

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 2 (restricted outdoor walker) and mobility grade 3 (unrestricted outdoor walker).

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and the size of the prosthetic foot.

| Spring stiffness                |                           |        |        |        |        |
|---------------------------------|---------------------------|--------|--------|--------|--------|
| Body weight [kg]                | Prosthetic foot size [cm] |        |        |        |        |
|                                 | 21, 22                    | 23, 24 | 25, 26 | 27, 28 | 29, 30 |
| 45 to 60 (99 lbs to 132 lbs)    | 1                         | 1      | –      | –      | –      |
| 61 to 80 (134 lbs to 176 lbs)   | 2                         | 2      | 1      | 1      | –      |
| 81 to 95 (178 lbs to 209 lbs)   | –                         | 3      | 2      | 2      | 1      |
| 96 to 110 (211 lbs to 242 lbs)  | –                         | –      | 3      | 3      | 2      |
| 111 to 125 (244 lbs to 275 lbs) | –                         | –      | –      | 4      | 3      |

### 2.3 Environmental conditions

| Storage and transport  |
|--|
| Temperature range $-20^{\circ}\text{C}$ to $+60^{\circ}\text{C}$ ( $-4^{\circ}\text{F}$ to $+140^{\circ}\text{F}$ ), relative humidity 20 % to 90 %, no mechanical vibrations or impacts |
| Allowable environmental conditions   |
| <b>Temperature range:</b> $-10^{\circ}\text{C}$ to $+45^{\circ}\text{C}$ ( $14^{\circ}\text{F}$ to $113^{\circ}\text{F}$ )   |
| <b>Chemicals/liquids:</b> fresh water, salt water, perspiration, urine, soapsuds, chlorine water   |
| <b>Moisture:</b> Submersion: max. 1 h in 3 m depth, relative humidity: no restrictions   |
| <b>Solids:</b> Dust, occasional contact with sand  |

#### Allowable environmental conditions

**Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (Cleaning).**

#### Unallowable environmental conditions

**Solids:** Highly hygroscopic particles (e.g. talcum), continuous contact with sand

**Chemicals/liquids:** Acids, continuous use in liquid media

## 2.4 Lifetime

The product was tested by the manufacturer with 2 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 3 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

**CAUTION!**

#### Risk of injury and risk of product damage

- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 11).
- ▶ Note the combination possibilities/comboination exclusions in the instructions for use of the products.
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.

**NOTICE!**

## Risk of product damage and limited functionality

- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

## Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

## 4 Scope of delivery

| Quantity | Designation          | Reference number  |
|----------|----------------------|-------------------|
| 1        | Instructions for use | -                 |
| 1        | Prosthetic foot      | -                 |
| 1        | Protective sock      | SL=Spectra-Sock-7 |

### Spare parts/accessories (not included in the scope of delivery)

| Designation    | Reference number |
|----------------|------------------|
| Footshell      | 2C3-1            |
| Connection cap | 2C10, 2C20       |

## 5 Preparing the product for use

### CAUTION

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### NOTICE

#### Grinding the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not grind the prosthetic foot or footshell.

## 5.1 Alignment

### INFORMATION

There is a plastic protector on the product's pyramid. It protects the connecting section from scratches during the alignment and fitting of the prosthesis.

- ▶ Remove the protector before the patient leaves the workshop/fitting area.

### 5.1.1 Applying/removing the footshell

- ▶ Apply or remove the footshell as described in the footshell instructions for use.

### 5.1.2 Bench Alignment

#### TT bench alignment

| Bench alignment process   |   |                 |   |          |    |          |    |          |    |
|---|---|-----------------|---|----------|----|----------|----|----------|----|
| <b>Required materials:</b> 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly) |   |                 |   |          |    |          |    |          |    |
| Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:   |   |                 |   |          |    |          |    |          |    |
| <b>Sagittal plane</b>   |   |                 |   |          |    |          |    |          |    |
| ①   | Heel height: <b>Effective heel height</b> (shoe heel height – sole thickness of fore-foot) + <b>5 mm</b>  |                 |   |          |    |          |    |          |    |
| ②   | Exterior foot rotation: <b>approx. 5°</b>   |                 |   |          |    |          |    |          |    |
| ③   | a-p positioning, middle of the prosthetic foot to the alignment reference line:   |                 |   |          |    |          |    |          |    |
|   | <table border="1"> <thead> <tr> <th>Foot size [cm]:</th> <th>Anterior placement, middle of foot to alignment reference line[mm]:</th> </tr> </thead> <tbody> <tr> <td>21 to 25</td> <td>25</td> </tr> <tr> <td>26 to 28</td> <td>30</td> </tr> <tr> <td>29 to 30</td> <td>35</td> </tr> </tbody> </table> | Foot size [cm]: | Anterior placement, middle of foot to alignment reference line[mm]: | 21 to 25 | 25 | 26 to 28 | 30 | 29 to 30 | 35 |
| Foot size [cm]:   | Anterior placement, middle of foot to alignment reference line[mm]:   |                 |   |          |    |          |    |          |    |
| 21 to 25  | 25  |                 |   |          |    |          |    |          |    |
| 26 to 28  | 30  |                 |   |          |    |          |    |          |    |
| 29 to 30  | 35  |                 |   |          |    |          |    |          |    |
| ④   | Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed.   |                 |   |          |    |          |    |          |    |
| ⑤   | Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line.<br>Socket flexion: <b>Individual residual limb flexion + 5°</b>  |                 |   |          |    |          |    |          |    |
| <b>Frontal plane</b>  |   |                 |   |          |    |          |    |          |    |
| ⑥   | Alignment reference line of prosthetic foot: <b>between big toe and second toe</b>  |                 |   |          |    |          |    |          |    |

| Bench alignment process |  |
|-------------------------|--|
|                         | Alignment reference line of prosthetic socket: <b>along the lateral patella edge</b> |
| <b>7</b>                | Observe the abduction or adduction position.   |

## TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

### 5.1.3 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219\***, TT modular leg prostheses: **646F336\***) may be requested from Ottobock.

### 5.1.4 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the plastic protector from the pyramid after completing the dynamic fitting and the walking exercises.

## 5.2 Optional: Installing the foam cover

The foam cover sits between the prosthetic socket and prosthetic foot. It is cut longer in order to compensate for the movements of the prosthetic foot and prosthetic knee joint. During flexion of the prosthetic knee joint, the foam cover undergoes posterior compression and anterior elongation. The foam cover should be stretched as little as possible in order to increase its service life. There is a connecting element (such as a connection plate, connection cap or connection cover) on the prosthetic foot.

- > **Required materials:** degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive
- 1) Measure the length of the foam cover on the prosthesis and add the length allowance.  
**TT prostheses:** Distal allowance for movement of the prosthetic foot.  
**TF prostheses:** Allowance proximal of the knee rotation point for flexion of the prosthetic knee joint and distal allowance for movement of the prosthetic foot.
  - 2) Cut the pre-shaped foam cover to length and fit it in the proximal area on the prosthetic socket.
  - 3) Pull the foam cover over the prosthesis.
  - 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or rests on the foot adapter.
  - 5) Install the prosthetic foot on the prosthesis.
  - 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
  - 7) Disassemble the prosthetic foot and remove the connecting element.
  - 8) Clean the connecting element using a degreasing cleaner.
  - 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
  - 10) Let the glue dry (approx. **10 minutes**).
  - 11) Install the prosthetic foot and adapt the exterior cosmetic shape. Take into account compression caused by cosmetic stockings or SuperSkin.

## 6 Maintenance and repair

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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