

Body Jacket by Measurement

Body shapes – information is required

Finished trim measurements

| ANATOMICAL MEASUREMENTS | | Body shape description |
|---|-----------|------------------------|
| | | |
| | | Special instructions |
| FINISHED TRIMLINES (TLSO or LSO) | | |
| ANTERIOR | POSTERIOR | |
| | | |

In order to provide the best possible fit for your patient, we have several measurements that differ slightly from normal orthometry. Instead of estimating the TLSO or LSO trim lines from the sternal notch and the pubis, we request your exact measurements from the waist trim line.

Body Jacket by Measurement

Customer information

Order Date _____ Ship Date _____

OB Customer No. _____ PO No. _____

Practitioner _____

Facility _____

Address _____

City/Province/Postal code _____

Patient information

Patient Name (optional) _____

Male Female, Age _____ Height _____ Weight _____

Diagnosis _____ Level of involvement _____

Shipping instructions

Base finished product is shell trimmed, vented and strapped

Opening

Bivalve

- Anterior under Posterior.
- Posterior under Anterior.

Plastic

Jacket is constructed of $\frac{5}{32}$ " low density polyethylene and a $\frac{3}{16}$ " foam liner.

Strapping

- Anterior Chafe & Loop.
 - Extended Chafe & Loop.
- Posterior Chafe & Loop.
 - Extended Chafe & Loop.
- Dacron© reinforced straps.
- Abdominal binder/support.
- Over the shoulder straps.
- Other _____

Trim

- LSO
- TLSO
- CTLSO
 - Cervical components
- Distal flaring (Additional modification and heat flaring, outside of the standard order, will be subject to an additional labour charge of \$120.00 per hour)

Options

- Sternal Pad with Adjustable Aluminium Bar
 - Ventilation holes
 - Yes
 - No
- Cervical Attachment

Special instructions

When your order form is faxed to Ottobock Fabrication Services, please call for conformation of receipt.

P 800 665 3327 (ext 4864 or 4878)

F 800 891 4514