

CUSTOM SILICONE – LINERS

If you have a client with skin breakdown problems, a uniquely shaped residual limb, or comfort and durability issues, Otto Bock has the answer.

FEATURES AND BENEFITS

Otto Bock Custom Silicone Liners use a proprietary silicone significantly more durable than the standard silicones used in the off-the-shelf liners, allowing us to vary both thickness and durometer anywhere you choose. It also lets us place softer gel pads right into the liner so you can target bony prominences or areas susceptible to breakdowns before they become a problem. To have us fabricate a Custom Silicone Liner please refer to the checklist below.

Additional features available include customized thickness, durometer, distal connection, side pins, anti-elongation material, gel pads, valve and custom colour. A lycra cover can be added to aid in donning if requested. The result is a custom liner designed to meet individual needs.

Otto Bock Custom Silicone liners offer:

- A completely customized fit
- High durability
- Cushioning wherever you need it
- Superior suspension
- Protection against skin irritation

One year warranty on manufacturing defects

CUSTOM SILICONE CHECKLIST – LINERS

— For Clinical Use Only —

PLEASE MAKE SURE TO INCLUDE THE FOLLOWING:

Cast of affected side

Please mark the following on the cast:

- Client name/ID
- Left or right
- Circumferential measurements
- Trim lines of liner
- Alignment lines
- Electrode sites (for upper extremity liners)
- Gel pad borders and apex points
- Trim lines of fabric

Completed order form on page 22



Can I make special requests? Yes! If you have any special requests or design features that you want incorporated, please call us at 800.665.3327 x4842 and we will do everything we can to help make your ideas a reality.