

**ottobock.**



## 2015 Release

# 3B1-2 Genium microprocessor-controlled knee

Reclaim your strength

From the powerful new walk-to-run function to enhanced Optimized Physiological Gait to the new Cockpit app, the Genium strengthens your options for improved outcomes.

Part of the only family of microprocessor knees, Genium capitalizes on more than 15 years of real-life experience, offering the choice of what knee is exactly right for you and your patient.

Contact your Sales Rep if you are interested in a demo or trial opportunity. Call (US) 800 328 4058 or (Canada) 800 665 3327.

## New

- ▶ Enhanced Optimized Physiological Gait
- ▶ Choice of two stance functions
- ▶ More control with new Cockpit app
- ▶ Manage rehab with an activity report
- ▶ Quick battery level indicator
- ▶ Free updated X-Soft 1.4
- ▶ New colors

# 3B1-2 Genium microprocessor-controlled knee

## Features and benefits

- Revolutionary Optimized Physiological Gait with improved software
- Confidently descend stairs and step over obstacles
- Chose from Intuitive Stance or Deliberate Stance functions
- Easily pick up the pace with the Walk2run function
- Your patient can easily access functions with the new Cockpit App for Android devices
- Maximize your patient's mobility with up to 5 days of smart power management

## Indications

- For bilateral or unilateral patients with knee disarticulation, transfemoral, and hip disarticulation
- For K Level 2 to 4
- Maximum weight is 150 kg/330lbs (125kg/275lbs with torsion unit 2R21)
- Contraindicated if the patient's physical and/or mental condition or living situation do not allow the correct handling of Genium

## Technical Data

<b>Suitable for</b>	K2 – K4
<b>Max. body weight</b>	330 lbs / 150 kg

## Ordering Information

3B1-2 Volcano Shadow, Genium with pyramid  
3B1-2=9.2 Desert Pearl, Genium with pyramid  
3B1-2=ST Volcano Shadow, Genium with screw top  
3B1-2=ST-9.2 Desert Pearl, Genium with screw top

## Your Genium purchase includes

- Genium knee joint
- 2R20 or 2R21 AXON Tube Adapter
- Inductive charger
- AC adapter
- Remote control
- Instruction manuals for clinicians (Qualified Personnel) and patients (User)
- Case for charger & AC adapter
- Wallet for cards
- Bluetooth® PIN Card
- Cockpit App
- Prosthesis Passport

**Download free X-Soft:** <http://ottobock.com/en/datastation>



With a splashproof IP 54 rating, you can eliminate the worry of rainy days and other incidental water contact.



Designed for extra stability, Deliberate Stance takes 1/8 second to lock and unlocks when unloaded or with 1° of knee extension. Bilateral patients, people who move more slowly, and those who kneel frequently preferred Deliberate Stance over Intuitive stance.



Choose between Volcano Shadow or Desert Pearl colors



Opt for a protective cover to keep the Genium at its best