

SiOCX TR Socket 8T350=1

Order form

Contact	<input type="text"/>	Customer number	<input type="text"/>	Date	<input type="text"/>
Customer			Shipping address (if different from customer address)		
Company	<input type="text"/>	Company	<input type="text"/>		
Street	<input type="text"/>	Street	<input type="text"/>		
Postal code/city	<input type="text"/>	Postal code/city	<input type="text"/>		
Email	<input type="text"/>	Phone	<input type="text"/>		
Patient ID	<input type="text"/>	PO#	<input type="text"/>	US Tax ID#	<input type="text"/>

Affected side:

- Left Right

Has the patient had a SiOCX socket before?

- yes no

Note: A well fitting trial prosthesis must be sent along with this order form, this includes: a well fitting check socket, fitting frame with wrist unit dummy attached in proper alignment and length.

Silicone inner socket

Colour of the inner socket

- OB swatch #.....
 Other solid colour

Attachment points

- Standard
 Own specification

Electrode

- Without 13E202 13E200

Myoelectric contact surfaces - Conductive Silicone

- Yes* (only used with 13E200 electrodes)
 No

Prepreg outer socket

Surface design

- Finished carbon design
 Water transfer finish*
 Lycra fabric colour

Componentry


- Attach lamination ring (wrist unit)
 Enclose
 Complete assembly of all components
 (please send all components with order)

Pull tube or valve options

- 12V10 Pull tube and valve for suction socket*
 21Y139=1 Exclusion valve* (if not using pull tube)

Rigid foam shaping

Olecranon-wrist measurement:.....mm



*Surcharge

Comments: