

ottobock.

FITNESS REVOLUTION

Foot and knee options
for active individuals



Image	K-Level	Material Number	Product Name	Description	Max Pt Weight	Inches - Clearance	L-Codes
	K4	1E2/1E2=1	ProCarve™ Foot	Sports	220 lbs	5 3/8" (138 mm) (w/ parts for ski) 4 3/4" (121 mm) (w/ footshell)	L5987
	K4	1E90	Springlite™ Sprinter	Sports	275 lbs	13 3/4" – 15 3/4" (348 – 394 mm) (depending on adapter)	L5980
	K4	1E91	Runner	Sports	275 lbs	11 3/4" (Approx. 300 mm)	L5980
	K4	1E93	Runner junior	Sports	99 lbs	7 1/8" (180 mm) 7 7/8" (201 mm) (w/ 4R224)	L5980
	K3, K4	1E95	Challenger™	Sports	240 lbs	6 7/8" (174 mm)	L5987
	K4	3S80/3S80=1	Modular Sports Knee Joint	Sports	220 lbs (100 kg)	3 3/4" (95 mm)	L5824 + L5925 + L5930
	K4	3R2	ProCarve™ Knee Joint	Sports	220 lbs (100 kg)	10 5/8" (270 mm)	L5830 + L5845 + L5848 + L5925 + L5930