

# Experience the difference of Ottobock Liners

**ottobock.**

Ottobock Liner	Similar competitor	Drop-in Replacement*
<b>SIL</b> 6Y42 Skeo	Össur Original and Synergy with cover	Yes
<b>SIL</b> 6Y43 Skeo Pure	Össur Original and Synergy uncovered	Yes
<b>SIL</b> 6Y70/75 Skeo/Skeo SkinGuard	Össur Iceross Comfort, Össur Iceross Dermo	Yes
<b>SIL</b> 6Y77 Skeo 3D TT	No competitor with 3D and flexion	Unique design requires new socket
<b>SIL</b> 6Y80/85 Skeo/Skeo Skinguard	Össur TF Standard	Yes
<b>SIL</b> 6Y87 Skeo 3D TF	Össur TF Conical	Yes but 6Y87 one size smaller
<b>TPE</b> 6Y90/92 Caleo	Willow Wood Alpha-Classic 3mm, ALPS 3mm uniform	Yes
<b>TPE</b> 6Y93 Caleo 3D	Willow Wood Alpha-Classic 6mm,	Yes
<b>PUR</b> 6Y512 Uneo 3D	No competitor with 3D and flexion	Unique design requires new socket

\* Liners with similar dimensions not requiring a new socket.

**SIL** Silicone

**TPE** Copolymer

**PUR** Polyurethane



# Unmatched liner selection

**ottobock.**



## **The right material for the right situation**

Ottobock is the **only** company to offer three different material choices  
Silicone, Urethane, TPE

## **Need easier flexion and protection over bony anatomy?**

Exclusive Anatomic Specific Geometry in our TT 3D liners (6Y512, 6Y77)

## **TF tissue control challenges?**

6Y87 Skeo 3D TF liner has unique conical distal shape with cylindrical proximal shape for a 3D fit and KISS compatible for rotation control.

## **Smelly Liners?**

Only Ottobock offers Skinguard® Antibacterial Technology (6Y512, 6Y522, 6Y523, 6Y75, 6Y85, Custom)

## **Skin irritation or loose feeling during sweating?**

Unique textured inner surface for a great connection (6Y512, 6Y43, 6Y75, 6Y77, 6Y85, 6Y87)

## **Difficulties visualizing the socket fit?**

Try our Pure “see-through” liner that’s also waterproof (6Y43)

## **Uncomfortable elongation effect on the limb?**

Integrated stretch reducing matrix on locking liners (6Y43, 6Y75, 6Y77, 6Y85, 6Y87, 6Y90)

## **When your patient’s limb is not “off-the-shelf”**

Only Ottobock has individual custom-made urethane liners