

F21 Maverick Xtreme AT

Product features

- Developed for extremely high loads
- Split keel design delivers +/- 16 degree inversion/ eversion for all-terrain control
- Smooth rollover during stance phase
- Lightweight fiberglass
- Innovative heel / keel design
- Waterproof



Technical data

K-Level	3, 4
Max. body weight	365 lbs / 166 kg
Sizes	23-30 cm
Heel height	3/8" / 10 mm
Weight with footshell*	28.15 oz / 798 g
Build height*	6 1/16" / 170 mm

* Reference size is 26 cm.

F21 Maverick Xtreme AT

Order form

Customer		Shipping address (if different from customer)	
Customer no.		Customer no.	
Company		Company	
Street		Street	
Postal code/city		Postal code/city	
O&P professional		Customer ref.	

Stiffness Chart

Body Weight	Activity Level		
	Low	Moderate	High
100-115 lbs / 44-52 kg	1	1	2
116-130 lbs / 53-59 kg	1	2	3
131-150 lbs / 60-68 kg	2	3	4
151-170 lbs / 69-77 kg	3	4	5
171-195 lbs / 78-88 kg	4	5	6
196-220 lbs / 89-100 kg	5	6	7
221-225 lbs / 101-116 kg	6	7	8
256-285 lbs / 117-130 kg	7	8	9
286-325 lbs / 131-147 kg	8	9	-
326-365 lbs / 148-166 kg	9	-	-

Scope of delivery

- F21* Maverick Xtreme AT
- 2F11=1* Set of heel wedges
- SL=SPECTRA-SOCK-7* Spectra protective sock black

Accessories & spare parts

- SL=SPECTRA-SOCK-7 Spectra-Sock, black
- 2F11=1 Set of Heel Wedges

Available footshell variants

- FTC-2A-1*-S*: Footshell Maverick sandal toe - available in size 23-30 cm

Foot module part number

Reference number	Stiffness	Size	Toe Configuration
F21 - 00 - 0	A	-	

Stiffness	Size (cm)	Toe Configuration
1-9	23, 24, ..., 30	SL Sandal toe, left SR Sandal toe, right

Footshell part number

Reference number	Generation	Size	Color	Toe	Side
FTC -	A - 1			4 -	

Generation	Size (cm)	Color	Toe	Side
2 Sandal toe	23, 24, ..., 30	4 Beige 15 Light Brown	S Sandal	L Left R Right