


Meridium. Reclaim your way.

Advantages at a glance:

- Meridium adapts to the current situation in real-time, tailored to the user's individual needs
- Maximum adaptation in real-time, with every step:
 - on varying terrain such as grass, forest paths, cobblestone or gravel
 - on ramps and inclines
 - to suit various walking speeds and loads
- Adapts to standing, walking, and walking backwards
- Enables full-surface foot contact when walking down stairs for enhanced safety
- Intuitive stance function offers stability regardless of the incline of the surface
- Automatic adaptation to different heel heights up to 5 cm
- Relief function for noticeably more comfortable and natural sitting or standing



Technical data

Mobility grade / K-Level	K3
Amputation level	unilateral TT, KD, TF, bilateral TT
Max. body weight	275 lbs (125 kg)
Sizes	24-29 cm
Weight (with footshell)	Sizes 24-25 approx. 1250 g Sizes 26-27 approx. 1500 g Sizes 28-29 approx. 1550 g
Structural height	S: 24-25 166 mm M: 26-27 175 mm L: 28-29 178 mm
Range of motion (size 27)	36.5° (22° PF; 14.5° DF)
Recommended knee joints for TF and KD amputees*	C-Leg, C-Leg compact, Genium, Genium X3

*All components are sold separately and are available Ottobock products that are compatible with the Meridium foot. To help insure optimal performance practitioners need to select components based upon individual patient criteria.



Meridium. Order form.

Customer		Shipping Address (if different from customer address)	
Customer no.		Customer no.	
Company		Company	
Address		Address	
Zip Code/City		Zip Code/City	
Practitioner		Com.	

Prosthesis user


Gender male female

Size _____ cm

Side of amputation left right bilateral

Color footshell translucent beige light brown

Mobility grade

Mobility grade 3 

Amputation level

Transtibial
 Knee Disarticulation
 Transfemoral

Prosthesis fitting

Initial fitting of prosthesis (of a microprocessor controlled foot by Ottobock)

Follow-up fitting (of a microprocessor controlled foot by Ottobock)

Prefitting with: _____
(Article number & Serial number present prosthetic foot)

Secondary fitting¹: _____
(Article number & Serial number present prosthetic foot)

¹ Order an additional foot)

Beside the 4E50-2 battery charger and 757L16-2 power supply, the order for the 1B1-2 Meridium also includes the 2C101 donning aid and the 2C7=* functional footshell (including 4G872=* cover cap set).

Warranty

SP-1B1-2=3
 three-year warranty package

SP-1B1-2=6
 six-year warranty package

Accessories

704G30 clamping tool

757P48 Y-adapter cable



Date: _____ Place: _____ Signature: _____