

Agilium® Vantage

Less Pain. More Life.

Meet the Agilium Vantage, a soft knee osteoarthritis (OA) brace for everyday use, recreational activities, and sports. The Agilium Vantage is a low-profile, wraparound knee orthosis designed to provide lasting pain relief for patients with mild unicompartmental osteoarthritis. The dynamic Y force strap system unloads the affected compartment of the knee joint with one pull, while reducing rotation and providing stability.

Features & Benefits

- Dynamic Y force strap system provides unloading to minimize rotation and reduce pain
- Single, low-profile metal upright with hinge for fixed flexion and extension limits
- Lightweight, breathable material prevents popliteal bunching for user comfort and increased compliance
- Four adjustable, numbered closure straps for quick and easy application
- Trimmable sleeves to fit a broader range of patients with three sizes to reduce inventory

Indications

- Mild to moderate unicompartmental knee osteoarthritis
- Medial or lateral knee pain
- Rheumatoid knee arthritis
- Post-operative meniscus reconstruction or ligament injuries needing unilateral relief



Sizing

	Article Number	Side	Size	A (in)	B (in)	Type
	50K306	L/R	XS-S	14-18	11-14	MO/LO
		L/R	M-L	18-23.5	14-18	MO/LO
		L/R	XL-3XL	23.5-32	18-24	MO/LO



Example format: 50K306=L-XS-S-MO or 50K306=R-M-L-LO
Place hinge on the side of the affected compartment