

Post-Operative Knee Braces

50K326=1 Post-Op Knee Brace, Cool

50K326=2 Post-Op Knee Brace, Full Foam

Suggested HCPCS Coding: L1832 Custom Fit / L1833 OTS

Universal telescoping ROM knee brace

Our telescoping post-operative knee braces are designed for range of motion control or knee immobilization after surgery.

Features & Benefits

- Easy-to-use telescoping system with adjustable strapping to provide a quick, custom fit for patients
- ROM hinge controls desired knee flexion and extension: limits motion between 0° to 120° and locks between 0° to 90° in 10° increments
- Quick-release buckles ease donning and doffing
- Soft grip padding limits migration and rotation
- Breathable, moisture-wicking material for user comfort
- Universal sizing with uprights adjusting from 19" to 26"
- Available in cool and full foam versions

Indications

- For locked or limited range of motion control following injury or operative procedures of the knee
- For use including but not limited to: Following major ligament surgery, meniscal repairs, patella realignment, regenerative chondroplasty, and stable femoral fracture

Tech Specs

Part Numbers	50K326=1 Post-Op Knee Brace, Cool 50K326=2 Post-Op Knee Brace, Full Foam
Size	Universal Measure leg circumference 6" above mid-patella Fits thigh circumferences up to 38"

