

# Malleo Immobil Walker

Quality for life



# Malleo Immobil Walker

The Malleo Immobil Walker product line is intended for conservative and post-operative foot and ankle fittings. From ligament, soft tissue and tendon injuries to fractures of the forefoot, metatarsus, ankle and tarsus, the Malleo Immobil Walker line offers a wide range of therapy options.

Uncomplicated adaptation using hook-and-loop closures makes for easy handling. The rounded rollover sole ensures a physiological gait pattern for all versions. All models in the Malleo Immobil Walker line are universal for the right and left side.

**1 50S11 Malleo Immobil Walker, low**

Indications for the 50S11 Malleo Immobil Walker low include the post-operative and conservative treatment of stable forefoot, metatarsus and / or ankle fractures, exostosectomy in case of hallux valgus as well as soft tissue injuries of the foot.

**2 50S10 Malleo Immobil Walker, high**

The 50S10 Malleo Immobil Walker high is intended for the treatment of ligament, soft tissue and tendon injuries, conservative treatment of forefoot, metatarsus and tarsus injuries, post-operative and conservative treatment of distal fibula fractures and stable forefoot, metatarsus and ankle fractures.

**3 50S14 Malleo Immobil Air Walker, low**

Post-operative and conservative treatment of stable forefoot, metatarsus and / or ankle fractures, exostosectomy in case of hallux valgus as well as soft tissue injuries of the foot. The integrated pump (Fig. 6) is used to increase compression to the desired level. This enhances foot stability and provides a pleasant cushioning effect for a high level of wearer comfort.

**4 50S12 Malleo Immobil Air Walker, high**

Indications for the 50S12 Malleo Immobil Air Walker high include the post-operative and conservative treatment of stable forefoot, metatarsus and / or ankle fractures, post-operative and conservative treatment of distal fibula fractures, soft tissue injuries of the foot and severe ankle distortion. An integrated pump (Fig. 6) is used to increase compression on the affected area for enhanced foot stability and a pleasant cushioning effect.

**5 50S15 Malleo Immobil ROM Walker**

The 50S15 Malleo Immobil ROM (range of motion) Walker is intended for the conservative treatment of forefoot, metatarsus and tarsus injuries, post-operative and conservative treatment of Achilles tendon ruptures, ligament, soft tissue and tendon injuries, post-operative and conservative treatment of distal fibula fractures and stable forefoot, metatarsus and / or ankle fractures. The joint of the ROM Walker (Fig. 7) is adjustable in 7.5° increments from 45° dorsal flexion to 45° plantar flexion, allowing it to be customised to patient requirements. The range of motion is extremely easy to adjust by means of a pin that is moved to a different setting.

**Malleo Immobil Walker**

	<b>Walker, high 50S10-1=*</b>	<b>Walker, low 50S11-1=*</b>	<b>Air Walker, high 50S12-1=*</b>	<b>Air Walker, low 50S14-1=*</b>	<b>ROM Walker, high 50S15-1=*</b>
<b>S</b>	31 – 35	31 – 35	31 – 35	31 – 35	31 – 35
<b>M</b>	36 – 42	36 – 42	36 – 42	36 – 42	36 – 42
<b>L</b>	43 – 45	43 – 45	43 – 45	43 – 45	43 – 45
<b>XL</b>	46 – 50	46 – 50	46 – 50	46 – 50	46 – 50



6 Integrated pump of the Malleo Immobil Air Walker



7 Joint of the Malleo Immobil ROM Walker

Otto Bock HealthCare GmbH  
Max-Näder-Straße 15 · 37115 Duderstadt/Germany  
T +49 5527 848-1706 · F +49 5527 72330  
export@ottobock.de · www.ottobock.com