

Residual Limb Compression Socks

Socket Technologies



Quality for life

Residual Limb Compression Socks

Our residual limb compression socks are an important complement to post-operative residual limb care. They ensure constant pressure on the residual limb, which decreases from distal to proximal, while at the same time offering a high level of wearer comfort for the user.

Advantages at a glance

- High level of wearer comfort and easy handling
- Good fixation of soft tissue and rapid dissipation of tissue swelling
- Can be used to prevent post-traumatic swelling after the amputation
- Even compression
- Low-cost, effective early care for the residual limb

Suspension band with silicone nubs (lower leg)

- Keeps the residual limb compression sock from sliding down
- Breathable with simultaneous adhesion

Area of Application

The residual limb compression sock is a post-operative replacement for wrapped bandages. It is used after amputation in order to prevent post-traumatic swelling, improve circulation in the residual limb and prevent oedemas after the prosthesis is removed.


We offer the residual limb compression socks for the thigh and lower leg in two compression classes and in various sizes in order to meet the individual requirements of your patients.

Contraindications

- For example arterial occlusive disease, congestive heart failure – the decision is made by your doctor

Order Information

For the thigh (with waist band)



1 Compression class selection

compression class 1 451F12
compression class 2 451F11

2 Size selection

Size	Circumference g	Circumference e
XS	41 – 44 cm	29 – 31 cm
S	44 – 48 cm	31 – 34 cm
M	48 – 52 cm	34 – 37 cm
L	52 – 56 cm	37 – 40 cm
XL	56 – 60 cm	40 – 43 cm
XXL	60 – 64 cm	43 – 46 cm


3 Length selection

Length (L)=g-e: 20, 25, 30, 35 cm

4 Order example

451F12=XS-20

For the lower leg (suspension band with silicone nubs)



1 Compression class selection

compression class 1 451F13
compression class 2 451F10

2 Size selection

Size	Circumference f	Circumference e	Circumference c
XS	39 – 41 cm	29 – 31 cm	27 – 29 cm
S	41 – 44 cm	31 – 34 cm	29 – 32 cm
M	44 – 47 cm	34 – 37 cm	32 – 35 cm
L	47 – 50 cm	37 – 40 cm	35 – 38 cm
XL	50 – 53 cm	40 – 43 cm	38 – 41 cm
XXL	60 – 64 cm	43 – 46 cm	41 – 44 cm

3 Length selection

Length (L)=f-c: 30, 38, 46 cm

4 Order example

451F13=XS-30-N



Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt / Germany
T +49 (0) 5527 848-0 · F +49 (0) 5527 848-1414
healthcare@ottobock.de · www.ottobock.com